Imuno-Rápido DENGUE NS1

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ENGLISH

SUMMARY

Dengue is an arbovirus caused by a *Flavivirus*, belonging to the *Flaviviridae* virus, and transmitted by the *Aedes aegypti* mosquito, and has been shown to have four subtypes (Den-1, Den-2, Den-3 e Den-4).

Its incubation period is 3 to 15 days, with an average of 5 to 6 days.

The clinical symptoms includes high fever (39°C to 40°C), a sudden onset, followed by headache, myalgia, prostration, arthralgia, retro-orbital pain, asthenia, anorexia, nausea, vomiting, rash, and itching skin. The hemorrhagic form of the disease are the most serious and can include: gingival bleeding, petechiae and ecchymoses, gastroenterorragia, shock and death. The important laboratory finding is thrombocytopenia. The Dengue

hemorrhagic fever occurs in 2 to 4% of reinfected individuals. Other serious complications of the disease are the encephalitis or dengue shock syndrome by and hepatitis.

The Flavivirus are enveloped RNA viruses, single stranded, whose genome encodes three structural proteins: C (core protein), M (membrane protein) and E (envelope protein). The non-structural proteins (NS) are NS1, NS2a, NS2b, NS3, NS4a, NS4b and NS5.

The NS1 is a non-structural glycoprotein that is essential for viral replication. During the acute phase of infection, the NS1 is found circulating in the serum of patients with detectable concentrations by immunological methods, being currently considered as a marker of infection by dengue virus before the appearance of IgM and IgG antibodies, allowing for early detection of the virus, 24 hours after the onset of symptoms, in addition to being found in primary and secondary infections.

Comparative studies show high correlation between the NS1 antigen and the molecular biology RT-PCR test in the detection of 4 subtypes of the dengue virus after 24 hours of the onset of symptoms. Thus, the NS1 test is quite useful for the diagnostic confirmation, not having cross-reactions with other flaviviruses, allowing for rapid and differential diagnosis from other diseases that manifest with similar febrile conditions, offering opportunity for a quick start on the treatment and, consequently, decrease in the mortality due to complications of the disease.

The **Imuno-RÁPIDO Dengue NS1 from WAMA** is a solid phase immunochromatographic test for the qualitative detection of the presence of NS1 antigen in whole blood, serum or plasma, with the purpose of early identification of individuals infected with the dengue virus.

PRINCIPLE OF THE METHOD

Samples of serum, plasma or whole blood containing NS1 antigen bind to conjugated anti-NS1 labeled with colloidal gold, whose immuno-complex migrates through capilarity through the test membrane binding to the anti-NS1 monoclonal antibodies present in the test area (T), determining the appearance of a colored band in this area. If the NS1 antigen is not present in the sample, there is no appearance of color band in the test area. In the control area (C) a colored band will always appear, demonstrating that the reagents are working correctly.

KIT PRESENTATION

[REF**]** **659010-R:** - **(10 determinations)**

1. Test Device: 10 Units

2. Instructions for use

[REF**]** **659020-R:** - **(20 determinations)**

1. Test Device: 20 Units

2. Instructions for use

[REF**]** **659025-R:** - **(25 determinations)**

1. Test Device: 25 Units

2. Instructions for use

[REF**]** **659040-R:** - **(40 determinations)**

1. Test Device: 40 Units

2. Instructions for use

[REF**]** **659080-R:** - **(80 determinations)**

1. Test Device: 80 Units

2. Instructions for use

PREPARATION AND REAGENT STABILITY

TEST DEVICE (1): Ready to use. Stable until expiration date when kept at 2 - 30°C. Do not freeze. If the device is stored in a refrigerator, let the device reach room temperature before performing the tests.

Note: The kit has the same performance after first use, and is stable until the expiration date described on the label, as long as it is kept in the indicated temperature (2-30°C).

MATERIALS REQUIRED BUT NOT PROVIDED

- Automatic pipette

- Absorbent paper

- Container for material disposal

- Timer

SPECIMENS

Use freshly-collected specimens of whole-blood, serum, or plasma free of hemolysis, lipemia, and contamination.

The sample can be stored covered in refrigerator at 2-8 °C for 48 hours. For long term storage, the serum or plasma specimen must be kept in a freezer at -20 °C.

Avoid repeated freezing and thawing as this may lead to false results.

WARNING: If the serum or plasma specimen is frozen, it must be completely thawed and thoroughly homogenized; after this, place it in upright position to allow for sedimentation of any particle left in the solution. Do not vortex the sample. Avoid making foam. **Diluted Sample can cause false-negative results.**

PROCEDURE


- Leave the test device and samples reach room temperature, if they were stored in the

refrigerator, before performing the test.

2. Remove the test device from the packaging and use immediately.

3. Place the test device on a flat dry surface.

4. Dispense 100 µl of the sample or control (no air bubbles) in the cavity of the sample on the test device (►) - see figure.


- Read the results at 15-20 minutes. Do not interpret the results after 20 minutes.

INTERPRETATION OF THE RESULTS

NEGATIVE


One colored band (reddish) appears in the control area (C).

POSITIVE


Two colored bands (reddish) appear, one in the test area (T) and the other in the control area (C).

Note: Any intensity of pink in the test area should be considered as positive.

INVALID

No band appears in the test area (T) and in the control area (C) or no band appears in the control area (C). If this happens the specimen must be tested again.

NOTE 1: results read after 20 minutes should not be considered.

NOTE 2: It´s recommended running a positive and a negative control of known value for each batch of tests to ensure the good performance of the assay.

PERFORMANCE CHARACTERISTICS

This comparative study with other commercial Kit, using the same methodology, showed that the **Imuno-RÁPIDO Dengue NS1 from WAMA** has 94.6% sensitivity and 98.6% specificity.

SENSITIVITY AND SPECIFICITY

	56 positive samples for NS1		Sensitivity
	Positives	Negatives	
IMUNO-Rápido Dengue NS1	53	3	94,6%

	74 negative samples for NS1		Especificidad
	Positives	Negatives	
IMUNO-Rápido Dengue NS1	1	73	98,6%

The **Imuno-RÁPIDO Dengue NS1 from WAMA** was compared with a RT-PCR molecular test using 180 samples of individual suspect of infection with the dengue virus, with the goal of determining its sensitivity and specificity. The results of this evaluation are shown

N = 180	RT-PCR		TOTAL
	Positive	Negative	
<i>Imuno</i> -RÁPIDO Dengue NS1	102	1	103
	6	71	77
TOTAL	108	72	180

This comparative study showed that the **Imuno-RÁPIDO Dengue NS1 from WAMA** has 94.4% sensitivity and 98.61% specificity.

PRECISION

Intra-Assay

The **Imuno-Rápido Dengue NS1** kit was tested for repeatability using 3 samples (low reactive, strong reactive and not reactive). 20 repetitions of the test were done by the same operator on the same day. The tests were done following the instructions for use.

According to the tests carried out, it was verified that the kit **Imuno-Rápido Dengue NS1** repeated satisfactorily the results with the samples tested under the same conditions. The shade of bands is within the expected variation and does not compromise the effectiveness of the test.

Inter-Assay

The **Imuno-Rápido Dengue NS1** kit was tested for reproducibility using 3 samples (low reactive, strong reactive and not reactive). 10 repetitions were done by diferente technicians on different days. The tests were done following the instructions for use.

According to the tests carried out, it was verified that the kit **Imuno-Rápido Dengue NS1**repeated satisfactorily the results with the samples tested even under different conditions. The shade of bands is within the expected variation and does not compromise the effectiveness of the test.

LIMITATIONS OF USE

The **Imuno-RÁPIDO Dengue NS1** is a screening test to characterize the presence of NS1 antigen, therefore, a result repeatedly positive with this test is a presumable evidence of the presence of the dengue virus in the sample.

Positive results must be confirmed by other tests.

As it occurs in other methodologies, false positive and false negative results may be found with this test whose proportion will depend on the prevalence of the disease on the assayed population.

As in any all diagnostic procedure, the results must be interpreted considering the clinical and epidemiological data for each case, along with other serological findings.

PRECAUTIONS AND WARNINGS

- For in vitro Diagnostic use only.
- Read carefully the instructions for use prior to performing the test.
- The expiry date refers to the last day of the month indicated on the labels of the test device and of the kit box. Do not use kit components after the expiration date.
- Avoid exposing the reagents to high temperatures and direct sunlight.
- Do not freeze the test device, for it may lead to irreversible deterioration.
- Allow the reagents to reach the room temperature before starting the test.

7. Do not replace components of this kit with components from manufacturers, or use components from different lots or different catalog numbers.

8. Discard the material following local regulations.

9. Use the Good Laboratory Practices (GLPs) in storage, handling and disposal of materials.

10. Single Use

11. Professional Use

12. Don´t use in case of damaged protective packaging

WARRANTY

WAMA Diagnóstica guarantees the exchange of the diagnostic kit, provided it is within the expiration date and is proven by its technical support that there were no technical mistakes in the execution, handling and storage of this product. **WAMA** and its distributors are not liable for failures on the performance of the kits under these conditions.

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BIBLIOGRAFIA / BIBLIOGRAPHY / BIBLIOGRAFIA

- ALCON-LEPODER, S. et al. The secreted form of dengue virus nonstructural protein NS1 is endocytosed by hepatocytes and accumulates in late endosomes: implications for viral infectivity. **J. Virol.**, v. 79, n. 17, p. 11403-11411, 2005.
- WORLD HEALTH ORGANIZATION. **Dengue: Guidelines for diagnosis, treatment, prevention and control**. Geneva: WHO, 2009.
- PEELING, R. W. et al. Evaluation of diagnostic tests: dengue. **Nat. Rev. Microbiol.**, v. 8, p. S30-S38, 2010.
- WHITEHORN, J.; FARRAR, J. Dengue. **Clin. Med.**, v. 11, n. 5, p. 483-487, 2011.
- SIMMONS, C. P et al. Dengue. **N. Engl. J. Med.**, v. 366, n. 15, p. 1423-1432, 2012.
- GANGULY, A. et al. Dengue diagnostics: current scenario. **Research Biotech.**, v. 4, n. 2, p. 19-25, 2013.

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