PCR TURBILÁTEX

PCR Turbilátex | CRP Turbilátex
Ref. 20.015.00

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ENGLISH

INTENDED USE

Kit intended for determination of the C Reactive protein (CRP) in serum.

PRINCIPLE OF THE METHOD

Latex particles covered with antibodies anti-human CRP reacts with CRP present in the sample. The agglutination produces increment of the absorbance and is proportional to the CRP concentration in the sample and may be determined at 540nm by comparison with a calibrator of known concentration.

REAGENT COMPOSITION

R DIL	TRIS buffer 20 mmol/L pH 8.2, preservatives	
R LAX	Latex particles covered with goat IgG anti C Reactive Protein, pH 7.3, preservatives	 
CAL	Human serum containing CRP; preservatives	 

All the components of human origin show negative results to HBs antigen, and anti-HIV. However, it may be treated as potentially infectious material.

STORAGE AND STABILITY

- After first opening the reagents are stable until the expiration date printed on the label if stored 2-8 °C tightly closed
- Do not use reagents over the expiration date.
- Store at 2 to 8 °C.
- Reagents must remain outside the specified temperature only the time required for testing.
- Do not freeze and protect from light.

ADDITIONAL EQUIPMENT

- Spectrophotometer or photometer with a thermostated cuvet able to read at 540 (530-550) nm.
- Thermostatic water bath at 37 °C.
- Pipettes and micropipettes.
- Clock or chronometer.
- Assay tubes.

WARNINGS AND PRECAUTIONS

- Perform the sample collection according to Good Laboratory Practices
- The kit is intended only for in vitro diagnostic use.
- The samples must be considered as potentially infective.
- Wear the individual protection equipment accordingly to Good Laboratory Practice.
- Discard the reactions surplus according to Good Laboratory Practice in a proper place for potentially infective material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product. The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, also is available on the site www.biotechnica.ind.br or calling for 55-35-3214-4646.
- Do not combine reagents from different batches
- Do not exchange the caps from different reagents in order to avoid cross contamination.
- Do not use the reagent if it displays any sign in disagreement with the one specified in MSDS of the product
- Use glass pipettes and disposable tips for each sample, control, standard/calibrator and reagent.
- The level of water in the thermostatic bath must be over the level of mix reaction within the assay tubes.
- The laboratories must establish the chemical requirements, microbiological and particulate matter for water prior to and define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
- The appropriate cleaning and drying of the lab material are fundamental factors for the stability of the reagents to reach correct results.

SAMPLE – PREPARATION AND STABILITY

Fresh serum.

Stable 7 days between 2 to 8 °C or 3 months at –20 °C.

Discard lipemic and/or hemolyzed samples.

Samples containing fibrin must be centrifuged for 10 min. at 3000 rpm before assay.

TECHNICAL PROCEDURE

A) REAGENT PREPARATION

Reagents (RDIL and RLAX) ready for use.

CRP Calibrator

The calibrator is lyophilized.

- Remove the refrigerator and leave at 10 30 °C for 10 minutes.
- Reconstitute with 1.0 ml of purified water.
- Dissolve with gentle shaking and allow to stand for 10 minutes to prevent foaming.
- Gently mixed by inversion until complete dissolution. After reconstitution, the calibrator is stable for 30 days if stored at 2 to 8 °C and 3 months if kept at -20.

B) PROCEDURE

- Heat the reagent at 37 °C.
- Zero the reading equipment with purified water.
- Pipette in the cuvet:

RDIL	800 µL
RLAX	200 µL
CAL or SAMPLE	5 µL

- Homogenize and insert in the thermostated cuvet at 37 °C. Trigger the timer.
- Measure the initial absorbance (A1) and 120 seconds after (A2) at 540nm (530-550).

CALIBRATION/CALCULATION PROCEDURES

CRP (mg/L) = $\frac{(A2 - A1) \text{ Sample} \times \text{X Calibrator concentration (mg/L)}}{(A2 - A1) \text{ Calibrator}}$

Exemple:

A1 - sample = 1,825

A2 - sample = 1,877

A1 - calibrator = 1,896

A2 - calibrator = 1,990

Calibrator concentration = 55 mg/L

CRP = $\frac{(1,877 - 1,825)}{(1,990 - 1,896)} \times 55 = \frac{0,052}{0,094} \times 55$

CRP = 30,4 mg/L

SENSITIVITY AND LINEARITY	
Analytical sensitivity: 1,92 mg/L. Linearity: 150 mg/L. For values higher than 150 UI/mL, dilute the sample with NaCl solution 150 mM (0.9%), repeat the assay and multiply the obtained result by the dilution factor. Prozone: The reaction has no prozone effect < 800 mg/L.	
TECHNICAL LIMITATIONS	
<ul style="list-style-type: none"> Anticoagulants: They cause measurement errors. Hemolyzed, Jaundice and Lipemic Sera 	
Bilirrubin > 20 mg/dL / Hemoglobin > 5g/L / Lipids > 10 g/L	
<ul style="list-style-type: none"> Others Rheumatoid Factors ≥ 800 UI/mL 	
QUALITY CONTROL	
Any clinical laboratory must keep an internal quality control program, which defines the aids, procedures and criteria for the tolerance limits, corrective actions and registration of the activities. Also, it must be kept a defined system for verification of the analytical variability that occurs in any measuring system. The use of controls to evaluate the imprecision of the analysis must be a routine practice in the lab. The lab must participate of external quality control programs. For the internal quality control of the lab it's indicated the use of the control sera below: CRP CAL (included) 20.015.00 Rheumatic Control Serum N1 21.007.00 Rheumatic Control Serum N2 21.008.00	
REF	
REFERENCE RANGES	
Normal Values: up to 6 mg/L These values are intended for orientation only. It's recommended that each lab establish its own reference range.	
PERFORMANCE CHARACTERISTICS	
Intra-Assay The realization of 20 determinations of the same sample at the same day showed a Coefficient of Variation of 2,75%. Inter-Assay The realization of 10 determinations of the same sample at different days showed a Coefficient of Variation of 2,81%.	
Analytical Specificity A comparison with a reference method showed a correlation coefficient (r) of 0.95 obtained from ambulatory samples. The result equation of the linear regression is Y = 0,955 X - 0,513 .	
CLINICAL SIGNIFICANCE	
The C-reactive protein (CRP) is an acute phase protein present in the serum of healthy individuals, which can significantly increase the majority of bacterial and viral infectious processes, tissue damage, inflammation and malignant neoplasms. The elevation of CRP concentration occurs after a few hours of the onset of the inflammatory process, reaching levels of 300 mg / L in 12-24 hours.	
QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE	
Before being approved for use BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@biotecnicatda.com.br	
AUTOMATION	
This product is compatible to the most types of biochemical automatic analysers. The applications are available at www.biotecnica.tda.ind.br	
ESPAÑOL	
FINALIDAD	
Kit destinado a la determinación cuantitativa de Proteína C Reactiva en suero humano.	
PRINCIPIO DEL MÉTODO	
Las partículas de látex recubiertas con anticuerpos anti-PCR humana son aglutinadas por la Proteína C-Reactiva (PCR) presente en la muestra. La aglutinación provoca un aumento de la absorbancia proporcional a la concentración de PCR en la muestra y por comparación con un calibrador de PCR de concentración conocida, se puede determinar el contenido de PCR en la muestra ensayada en 540 nm.	
COMPOSICIÓN DE LOS REACTIVOS	

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