TIRA DE URINA

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ENGLISH

INTENDED USE

The Urine Strip Kit consists of plastic strips on which are affixed reagent areas. The test is intended for the qualitative and semi-quantitative analysis of the following analytes: glucose, bilirubin, ketone, Specific Gravity, blood, pH, protein, Urobilinogen, nitrite, and leukocytes.

PRINCIPLE OF THE METHOD

Leukocytes (LEU): This test reveals the presence of esterase of granulocytes. The esterase reacts with the amino pyrazol acid ester releasing hydroxypyrazole. The diazonium salt is reacted with pyrazol producing a violet color. The intensity of the color developed is proportional to the number of leukocytes present in the urine specimen.

Nitrite (NIT): The test is specific to nitrites, which are formed by the reduction of nitrates due to the action of reductases produced by Gram negative bacteria. In an acidic medium, nitrite in the urine reacts with p-arsanilic acid to form a diazonium compound. The same binds with 1N- (1-naphthyl) ethylenediamine to produce a pink color.

Urobilinogen (URO): This test is based on modified Ehrlich reaction between p-diethylaminobenzaldehyde and urobilinogen in a strongly acidic medium which produces a pink color.

Protein (PRO): This reaction is based on the color change from blue indicator tetrabromofenol in the presence of protein. Results vary from yellow to green.
pH: The pH determination is based on the Methyl double red indicator system / bromothymol blue. Results vary from orange to blue-green.

Blood (BLO): The test is based on hemoglobin pseudoperoxidase activity that catalyzes the reaction diisopropylbenzene dihydroperoxide and 3,3', 5,5'-tetramethylbenzidine. A uniform color ranging from yellow to blue indicates the presence of myoglobin, hemoglobin or hemolyzed erythrocytes. Green-blue spots scattered or compact spots indicate intact erythrocytes according to the results of scale.

Specific Gravity (SG): The specific gravity of the urine is related to the ionic concentration, based on the apparent pKa change. The resulting color ranges from dark bluish green to yellowish green.

Ketone (KET): This test is based on the reaction of ketones with nitroprusside and acetoacetic acid resulting a change in color, varying from light pink to purple. Ketones are normally not present in urine. The presence of ketone may occur in urine during physiological stress such as: fasting, pregnancy and excessive exercise 4-6. In diets with food restriction or abnormal carbohydrate metabolism situations, ketones appear in the urine in high concentration7.

Bilirubin (BIL): This test is based on binding reaction of bilirubin with diazotized dichloroaniline in a strongly acid medium. The variation in bilirubin produces a color intensity proportional to its concentration in urine.

Glucose (GLU): glucose, when present in the urine is oxidized by glucose oxidase to form gluconic acid and hydrogen peroxide. The hydrogen peroxide reacts with the chromogen potassium iodide in the presence of peroxidase. The amount of oxidised chromogen produced determines the color, ranging from green to brown. Small amounts of glucose can be excreted by the kidney 3. Glucose concentrations near 100 mg / dl, which are repeated in various analyses can be considered abnormal. The test is specific for glucose not reacting with lactose, galactose, fructose or other metabolic substances, not even with drug metabolites (salicylates and nalidixic acid).

PERFORMANCE CHARACTERISTICS / COMPOSITION

Reagent surface	Composition	Performance Characteristics
Leukocytes (LEU)	Derivado de éster ácido amino pirazol0,5%;saldia zônio0,4%; tampão borato 32% and stabilizer 67,1%.	Negative at 500 Leu/ μ L (+++).
Nitrite (NIT)	p-arsanilic acid 4,5%; N-(1-naphthyl) ethylenediamine 0,5%; Citrate Buffer 45% e polyvinylpyrrolidone 50%.	Negative and Positive
Urobilinogen (URO)	p-diethylaminobenzaldehyde 2,5%; Sulfosalic acid 40 and stabilizer 57,5%.	3,5 a 200 μ mol/L.
Protein (PRO)	tetrabromophenol blue 0,3%; Citrate Buffer 94,7%; and stabilizer 5%.	Negative at 2000 mg/dL (20 g/L).
pH	methyl red sodium salt 0,5 %; bromthymol blue 5% and stabilizer 94,5%.	5,0 and 9,0.
Blood (BLO)	3,3',5,5'-tetramethylbenzidine (TMB) 4%; diisopropylbenzenedi hydroperoxide6%; phosphate Buffer 45% and stabilizer 45%.	Negative at +++ for hemolysis and 5-10 a 50 Ery/ μ L.
Specific Gravity (SG)	bromthymol blue indicator 2,5%; poly (methyl vinyl ether/maleic anhydride) 55%; sodium hydroxide 25%and stabilizer 17,5%.	1,000 and 1,030.
Ketone (KET)	sodium nitroprusside 3%; Citrate Buffer 30 and stabilizer 67%.	Negative at 160 mg/dL (16 mmol/L).
Bilirubin (BIL)	2,6-dichloroaniline 0,5%;Citrate Buffer 69,5 and stabilizer 30%.	Negative at 4,0 mg/dL (70 μ mol/L).
Glucose (GLU)	glucose oxidase 0,5% ;peroxidase0,5,potassiumiodide4,0%; Citrate buffer 80%and stabilizer 15%.	Negative at 2000 mg/dL (110 mmol/L).

STORAGE AND STABILITY

- Stability: Stable until the kit expiration date which is printed on the package label.
- Do not use strips whose expiration date has expired..
- Keep from 15° to 30°C.

WARNINGS AND PRECAUTIONS

- Keep the tube closed until time of usage and the desiccant inside.
- Do not touch the reagent areas of the strip.
- Discard any strip with color change.
- The kit is intended only for in vitro diagnostic use.