FINALIDADE

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ENGLISH

INTENDED USE

Kit for the determination of anti-streptolysin O antibodies (ASLO) in serum. In vitro diagnostic use.

STORAGE AND HANDLING

- Store at 2 to 8 ° C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE

ASLO Latex is a plate agglutination method for qualitative and semi-quantitative detection, where the latex particles coated with streptolysin O are agglutinated by anti-streptolysin O antibodies present in the sample.

SAMPLES – TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION

Sample Type: Serum.

Collection, handling and preparation: Perform the sample collection according to Good Clinical Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.

Preservation: ASLO in serum, is stable for 7 days stored at 2 to 8 ° C or 3 months at a temperature of -20 ° C.

PRODUCT DESCRIPTION

R LAX	Suspension of Latex particles coated with Streptolysin O, pH 8.2. Preservative. The sensitivity of RLAX is traceable to NIBSC 97/662 (National Institute for Biological Standards and Control).	
CONTROL +	Human serum with a concentration of ASO> 200 IU / ml. Preservative. (INCLUDED IN THE ASLO LATEX KIT)	
CONTROL -	Human serum in buffer Phosphate, stabilizer and preservative. (INCLUDED IN THE ASLO LATEX KIT)	

STABILITY IN USE

- After opening, the product (RLAX, CONTROL (+) and CONTROL (-)) in use is stable until the expiration date printed on the label, since recommended storage conditions are followed (2 to 8 ° C).
- Reagents should remain outside the specified temperature only for the time required to perform the tests.

TREATMENT AND HANDLING OF THE PRODUCT

A) REAGENT PREPARATION

RLAX: Ready for use.
CONTROL (+) and CONTROL (-): Ready for use.

B) OPERATING INTERVAL FOR ASLO LÁTEX RLAX

Up to 800 IU / mL.

QUALITY CONTROL

The use of controls should be routine practice in the laboratory. It is recommended to use Positive and Negative Controls to monitor the performance of the procedure and to serve as a comparative model for a better interpretation of the results. For Internal Quality Control Laboratory it is recommended to use the controls below:

CONTROL (+) (ASLO LÁTEX KIT)	REF	30.001.00
CONTROL (-) (ASLO LÁTEX KIT)		30.001.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

Qualitative Method

- 1) Before starting the tests, allow the RLAX, controls and samples to reach room temperature. Homogenize before use.
- 2) On the reaction plate use a different circle for each sample or control.
- 3) Place 20 µL of the sample / controls.
- 4) In each circle add 20 µL of RLAX.
- 5) Mix each determination using disposable rods (one for each sample).
- 6) Homogenize the plate for 2 minutes in circular motions, tilting it back and forth. Observe the result under artificial light: Presence or absence of agglutination.

Semi-quantitative method

- 1) Dilute the sample with saline (0.9% NaCl) as follows:

Dilution	Saline (µL)	Sample(µL)	Dilution Factor
1/2	100	100 (pure)	2
1/4	100	100 (sol.1/2)	4
1/8	100	100 (sol.1/4)	8
1/16	100	100 (sol.1/8)	16

- 2) Carry out the test of each dilution according to the qualitative method.

B) CALCULATIONS

The approximate concentration in the semi-quantitative test is obtained by multiplying the dilution factor by the sensitivity of the test (200 IU/mL).

Example:

AgAgglutination up to 1/4 dilution.

DiDilution Factor = 4

ASASLO = 800 IU / mL (4 x 200)

C) INTERPRETATION

Streptolysin O is a toxic immunogenic exoenzyme produced by β-hemolytic *Streptococcus* groups A, C and G. The detection of ASLO antibodies is used for the diagnosis and treatment of diseases such as rheumatic fever, acute glomerulonephritis, and other streptococcal infections. These infections can later lead to impairment of the heart or kidneys. Early diagnosis, efficient treatment and patient monitoring can reduce these risks.

INTERFERING AND LIMITATIONS

Hemolysis, Jaundice, Lipemia and Rheumatoid Factor: Hemoglobin> 1000 mg/dL / Bilirubin> 40 mg/dL / Triglycerides> 2000 mg/dL and Rheumatoid Factor> 300 IU/mL may interfere with the determination. Medications: consult recommended literature reference (Young, 2000).

PERFORMANCE CHARACTERISTICS

Sensitivity: Detection limit: 200 IU/mL.

Pro-zone effect: Pro-zone high-dose effect is not observed at concentrations of 800 IU/mL.

Analytical Specificity: The product specifically determines Anti-Streptolysin O antibodies in the presence of interferents in the sample, to the concentrations reported above.

Accuracy: To estimate the concordance percentages, the method was compared with a similar one by the determination of 100 samples, which were obtained: Percentage of Positive Agreement: 98% Negative Agreement Percentage: 100%

Precision: In order to estimate the intra-assay precision, the study was performed in two daily runs, in duplicate, for 10 days and for the inter-run one daily run of four replicates, for 5 days. In both, four samples were used, two positive and two negative with 100% agreement.

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange reagent bottle caps to avoid cross-contamination.
- Do not use the reagent when it has a visual characteristic that is not in accordance with the product's MSDS requirements.
- Avoid leaving the reagents out of the specified storage conditions when they are not in use.
- Use specific glass pipettes and disposable tips for each sample, control, standard / calibrator and reagent.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the used material are fundamental factors for the stability of the reagents and obtaining correct results.

REFERENCE RANGES

Adults	Until 200 UI/mL
--------	-----------------

These values are for guidance only and it is recommended that each laboratory establishes its own reference range.

MATERIALS REQUIRED TO CARRY OUT THE TEST

- Reaction plate.
- Disposable rod.
- Automatic pipette.
- Clock or Stopwatch

ALERTS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Disposal, Safety and First Aid information are described in the Individual Product Safety Data Sheet (MSDS) for this product, available at www.biotechnica.ind.br or by telephone (35) -3214-4646.
- Dispose of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS).

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@biotechnicaltda.com.br

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us of its Immunology, Diagnostic Criteria, and Prophylaxis. *Quartely Journal Of Medicine*, New Series 84, n° 305, p.641-658, 1992.

- KASSEM, Ahmed Samir; ZAHER, Salah Rafik. An International Comparison of the Prevalence of Streptococcal Infections and Rheumatic Fever in Children. *Pediatric Annals*. v. 21, n° 12, p.835-842, 1992.

- KLEIN, George C.; BAKER, Carolyn N.; JONES, Wallis L. "Upper Limits of Normal" Antistreptolysin O and Antideoxyribonuclease B Titers. *American Society for Microbiology*, v. 21, n° 6, p.999-1001, 1971.

- SPAUN, J. et al. International Standard for Antistreptolysin-O. *Bull. Wild Hlth Org.* v. 24, p.271-279, 1961.

- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. *Clin. Chem.* v.27 p.493-501, 1981.

- YOUNG, D.S. Effects of drugs on clinical laboratory tests - vol. 2, 5 ed. Washington DC: AAC Press, 2000.

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