ASO TURBILÁTEX

ASO Turbilátex | ASO Turbilátex
Ref. 20.004.00

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ENGLISH

INTENDED USE

Kit for the determination of anti-streptolysin O (ASO) antibodies in serum. Use *in vitro* diagnosis.

STORAGE AND HANDLING

- Store at 2 to 8 °C.
- Protect from light.
- Kit expiration date is printed on the label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE

Immunoturbidimetric method in which latex particles coated with Streptolysin O are agglutinated by anti-streptolysin O (ASO) antibodies present in the sample. The agglutination causes a turbidity whose intensity is proportional to the concentration of ASO.

SAMPLES – TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION

Sample Type: Serum.

Collection, handling and preparation: Perform the sample collection according to Good Clinical Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.

Preservation: ASO in Serum is stable for 7 days at a temperature of 2 to 8 °C or 3 months at a temperature of -20 °C.

PRODUCT DESCRIPTION

R DIL	Tris buffer 20 mmol / L pH 8.2, preservative.	
R LAX	Suspension of latex particles coated with Streptolysin O, pH 10,0, preservative. Human serum containing anti-streptolysin O, preservative. The ASO value is traceable to the International Standard for anti-streptolysin O, NIBSC (97/662).	
CAL		

STABILITY IN USE

- Once opened, the product (RDIL, RLAX and CAL) in use is stable up to the expiration date printed on the label, provided that it is followed by recommended storage conditions (2 to 8 °C).
- The stability of the working reagent is 20 days, provided that they are followed by recommended conditions of preparation and storage (2 to 8 °C).
- Reagents should remain outside the specified temperature only for the time required to perform the tests.

TREATMENT AND HANDLING OF THE PRODUCT

A) REAGENT PREPARATION

CAL: Ready-to-use.

Working Reagent (RT): Mix in proportion: 4 parts RDIL + 1 part RLAX. Homogenize gently.

B) OPERATING INTERVAL

The operating range of the product is from 34,2 UI/mL to 736 UI/mL

For higher values, dilute the sample with 150 mM NaCl (0.9%), carry out a new dosage and multiply the result obtained by the dilution factor.

QUALITY CONTROL

The use of controls should be routine practice in the laboratory. It is suggested to use one control in the reference range or at the decision level and another control with value in another range of clinical significance. For internal control of laboratory quality it is recommended to use the control sera below:

ASO Turbilátex CAL (Kit Included)	20.004.00
Rheumatic Control Level 1	21.007.00
Rheumatic Control Level 2	21.008.00

REF

TEST PROCEDURE, CALCULATION AND INTERPRETATION

REACTION PARAMETERS	
Type of reaction	Fixed time
Wave-length	540 nm
Temperature	37 °C
Sample Volume / CAL *	3 µL
RT* Volume	300 µL
Reading 1 (A ₁)	Immediately
Reading 2 (A ₂)	2 min.

* Sample volumes / CAL and RT can be modified as long as the established ratio is maintained.

B) CALCULATIONS

The equipment automatically calculates the analyte concentration in each sample.

C) INTERPRETATION

Streptolysin O is a toxic immunogenic exoenzyme produced by β-hemolytic Streptococci from groups A, C and G. The detection of ASLO antibodies is used for the diagnosis and treatment of diseases such as rheumatic fever, acute glomerulonephritis, and other streptococcal infections. These infections can later lead to impairment of the heart or kidneys. Early diagnosis, efficient treatment, and patient monitoring can reduce these risks.

INTERFERING AND LIMITATIONS

Hemolysis, Jaundice and Lipemia: Hemoglobin > 1000 mg / dL, Bilirubin > 40 mg / dL and Triglycerides > 1000 mg / dL interfere in the determination.

Medications: consult recommended literature reference (Young, 2000).

PERFORMANCE CHARACTERISTICS

Sensitivity: Limit of detection: 7,2 UI/mL / Limit of quantification: 34,2 UI/mL
Prozone High Dose Effect: No prozone effect is observed up to 4000 IU / mL of ASO.

Analytical Specificity: The product specifically determines ASO in the presence of other interfering substances in the sample, up to the concentrations reported above.

Accuracy: The method was compared with a similar method by the determination of 40 samples in duplicate. The regression equation $Y = 0,9999X - 0,053$ and correlation coefficient $r = 1,000$ was obtained. Using this equation, the estimated total systematic error of -0.15% for a level of 120 IU / mL and -0.12% for a level of 500 IU / mL.

Precision: It was determined using samples at 03 decision levels, two runs in duplicate per day for 20 days, obtaining:

Samples (UI/mL)	Repetitions	In-run Precision		Total Precision	
		SD (UI/mL)	%CV	SD (UI/mL)	%CV
79,3	80	1,33	1,7	1,56	2,0
317	80	8,14	2,6	8,51	2,7
457	80	8,92	2,0	9,54	2,1

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange reagent bottle caps to avoid cross-contamination.
- Do not use the reagent when it exhibits a visual characteristic in disagreement with that specified in the product's MSDS.
- Avoid leaving reagents out of the specified storage conditions when they are not in use.
- Use specific glass pipettes and disposable tips for each sample, control, standard / calibrator, and reagent.
- The water level of the water bath should be higher than that of the test tubes containing the reactions.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the requirements by means of periodic checks.

- Proper cleaning and drying of the used material are fundamental factors for the stability of the reagents and obtaining correct results.

REFERENCE RANGES

Grown ups	Up to 200 UI/mL
-----------	-----------------

These values are for guidance only and it is recommended that each laboratory establishes its own reference range.

MATERIALS REQUIRED TO CARRY OUT THE TEST

- Spectrophotometer or photometer for reading at 540 (530 – 550) nm
- Glass and / or automatic pipettes.
- Clock or Stopwatch.
- Test tubes

ALERTS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Disposal, Safety and First Aid information are described in the Individual Product Safety Data Sheet (MSDS) for this product, available at www.biotechnica.ind.br or by telephone (35) -3214-4646.
- Dispose of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS).

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@biotechnicaltda.com.br

AUTOMATION

This product is compatible to the most types of biochemical automatic analyzers. The applications are available at www.biotechnicaltda.ind.br

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- GÁMBARO, M.A.S.; *et al.* Intra- and Inter-Individual Biological Variability Data Blank. *Eur J Clin Chem Clin Biochem.* v. 35, p. 845-852, 1997.
- HAFEEJEE, I.; Rheumatic Fever and Rheumatic Heart Disease: The Current Status of its Immunology, Diagnostic Criteria, and Prophylaxis. *Quarterly Journal of Medicine.* v. 84, n. 305, p. 641-658, 1992.
- KLEIN, G.C.; *et al.* "Upper Limits of Normal" Antistreptolysin O and Antideoxy-ribonuclease B Titers. *Applied Microbiology.* v. 21, n. 6, p. 999-1001, 1971.
- SPAUN, J.; *et al.* International Standard for Antistreptolysin-O. *Bull Wild Hlth Org.* v. 24, p. 271-279, 1961.
- WANNAMAKER, L.W.; *et al.* Antibody Titers in Acute Rheumatic Fever. *Circulation.* v. 21, p. 598-614, 1960.
- YOUNG, D.S. Effects of drugs on clinical laboratory tests - vol. 2, 5 ed. Washington DC: AACC Press, 2000.

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