Alfa-1 Glicoproteína Ácida Plus

Alpha-1 Glycoprotein Acid Plus |
Alfa-1 Glicoproteína Ácida Plus
Ref. 20.018.00

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ENGLISH

INTENDED USE

Kit for the determination of alpha-1 acid glycoprotein in serum and plasma. In vitro diagnostic use.

STORAGE AND HANDLING

- Store at 2 to 8 °C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE

Immunoturbidimetric method in which anti-human alpha-1-glycoprotein antibodies form insoluble complexes when mixed with samples containing alpha-1-acid glycoprotein. The turbidity intensity is proportional to the acid alpha-1-glycoprotein concentration of the sample.

SAMPLES – TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION

Sample Type: Serum and plasma (heparin or EDTA).

Collection, handling and preparation: Perform the sample collection according to Good Clinical Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.

Preservation: Alpha-1 acid glycoprotein in serum and plasma is stable for 7 days at 2 to 8 °C or 3 months at -20 °C.

PRODUCT DESCRIPTION

R 1	Goat anti-alpha-1 antibodies glycoprotein-human acid; protein; Tris buffer 20 mmol / L pH 8.2; preservative.	
------------	--	--

STABILITY IN USE

- After opening, the product (R1) in use is stable to the expiration date printed on the label, since recommended storage conditions (2 to 8 °C) are followed.
- Reagents should remain outside the specified temperature only for the time required to perform the tests.

TREATMENT AND HANDLING OF THE PRODUCT

A) REAGENT PREPARATION

R1: Reagent ready for use.

B) OPERATING INTERVAL

The operating range of the product is from 7.4 mg / dL to 209 mg / dL. For higher values, dilute the sample with 150 mM NaCl (0.9%), carry out a new dosage and multiply the result obtained by the dilution factor.

QUALITY CONTROL

The use of controls should be routine practice in the laboratory. For Calibration and Internal Quality Control Laboratory it is recommended to use the calibrator and control below:

Multicalibrador Plus	REF	21.011.00
Medium Multiparameter Control		21.003.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

REACTION PARAMETERS	
Type of reaction	Final Point / Growing
Wave length	340 nm
Temperature	37 °C
Sample volume/CAL*	2 µL
R1 volume*	285 µL
Incubation time	4 minutes
Calibration	Calibration curve**

* Sample volumes / CAL and R1 can be modified as long as the established ratio is maintained.

**CALIBRATION CURVE:

Dilute the PLUS Multicalibrador in NaCl 9g/L as follows:

Dilution	1	2	3	4	5	6
Calibrador (µL)	-	10	25	50	75	100
NaCl 9g/L (µL)	100	90	75	50	25	-
Factor	0	0,1	0,25	0,5	0,75	1

Multiply the concentration of alpha-1 acid glycoprotein indicated in the Multicalibrador PLUS use instruction by the corresponding factor to obtain the concentration of alpha-1 acid glycoprotein at each dilution. The concentration of alpha-1-glycoprotein of the calibrator is traceable to the Standard Reference Materials ERM DA470k / IFCC.

B) CALCULATIONS

The equipment automatically calculates the analyte concentration in each sample.

C) INTERPRETATION

Alpha-1-acid glycoprotein is an acute phase protein that exerts a regulatory effect influencing the cascade of inflammation. It is synthesized by hepatic parenchymal cells, also granulocytes and monocytes contribute significantly to increase plasma levels during sepsis. The plasma concentration of alpha-1-acid glycoprotein increases by 3 to 4 times in most conditions associated with inflammatory processes and tissue necrosis. Levels may also be increased, along with haptoglobin and pre-albumin, by the effect of both endogenous and exogenous glucocorticoids.

INTERFERING AND LIMITATIONS

Hemolysis, Jaundice and Lipemia: Hemoglobin > 400 mg / dL / Bilirubin > 20 mg / dL, Triglicéridos > 500 mg/dL, interfere in the dosage.

Medications: see recommended literature reference (Young, 2000).

PERFORMANCE CHARACTERISTICS

Sensitivity: Limit of detection: 0.9 mg / dL / Limit of quantification: 7.4 mg / dL.

Prozone High Dose Effect: There is no prozone effect up to 750 mg / dL of Alpha-1 acid glycoprotein.

Analytical Specificity: The product specifically determines Alpha-1 acid glycoprotein in the presence of other interfering substances in the sample up to the concentrations reported above.

Accuracy: The method was compared with a similar method by the determination of 19 samples in duplicate. The regression equation $Y = 1.017X - 5.60$ and correlation coefficient $r = 0.9767$ was obtained. Using this equation the estimated total systematic error of -9.5% for a level of 50 mg / dL and -1.8% for a level of 160 mg / dL.

Accuracy: It was determined using samples at 03 decimal levels, two runs in duplicate per day for 20 days, being obtained.

Samples (mg/dL)	Repetitions	In-run Precision		Total Precision	
		SD (mg/dL)	%CV	SD (mg/dL)	%CV
55	80	0,91	1,7	1,24	2,3
123	80	2,03	1,6	2,29	1,9
228	80	3,64	1,6	6,11	2,7

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange reagent bottle caps to avoid cross-contamination.
- Do not use the reagent when it has a visual characteristic that is not in accordance with the product's MSDS requirements.
- Avoid leaving the reagents out of the specified storage conditions when they are not in use.
- Use specific glass pipettes and disposable tips for each sample, control, standard / calibrator, and reagent.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the used material are fundamental factors for the stability of the reagents and obtaining correct results.

REFERENCE RANGES

Adults	50 - 120 mg/dL
--------	----------------

These values are for guidance only and it is recommended that each laboratory establishes its own reference range.

MATERIALS REQUIRED TO CARRY OUT THE TEST

- Spectrophotometer or photometer for reading at 340 nm.
- Glass and / or automatic pipettes.
- Multicalibrator Plus Biotécnica – CAT BT 21.011.00

ALERTS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Disposal, Safety and First Aid information are described in the Individual Product Safety Data Sheet (MSDS) for this product, available at www.biotechnica.ind.br or by telephone (35) -3214-4646.
- Dispose of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS).

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use Biotécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiry date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@biotechnicaltda.com.br

AUTOMATION

This product is compatible to the most types of biochemical automatic analysers. The applications are available at www.biotechnicaltda.ind.br

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