Bilirrubina Direta e Total
 Direct & Total Bilirubin | Responsável Técnico:
 Dr. Gilson Sérgio Pizzo
 Bilirrubina Directa y Total
 Ref. 10.003.00
 CRF MG – 5310
MS 80027310205

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Recém-nascido		Bilirrubina Total (mg/dL)	
Idade	Prematuro	A termo	
0-1 dia	1,0 a 8,0	2,0 a 6,0	
1-2 dias	6,0 a 12,0	6,0 a 10,0	
3-5 dias	10,0 a 14,0	4,0 a 8,0	

Adultos		mg/dL
Bilirrubina Direta		0,0 a 0,2
Bilirrubina Total		0,0 a 2,0

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ENGLISH
INTENDED USE
 Kit for determination of Direct and Total Bilirubin in serum. In vitro diagnostic use.
STORAGE AND HANDLING
 • Store at 2 to 8 ° C.
 • Protect from light.
 • Stable until the kit expiration date that is printed on the label.
 • Do not use reagents whose date has expired.

WORKING PRINCIPLE
 The determination is based on the diazo reaction, which colorimetrically quantifies the red staining of azobilirubin when bilirubin reacts under certain conditions with diazotized sulfanilic acid. The intensity of the formed color is proportional to the concentration of bilirubin in the sample. In the presence of dimethylsulfoxide (accelerator), the indirect form (unconjugated) and the direct form (monoconjugate, diconguate and protein bound) are reacted to determine Total Bilirubin. The direct form is determined without the accelerator.
SAMPLES – TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION
Sample Type: Serum.

Collection, handling and preparation: Perform the sample collection according to Good Clinical Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material. Immediately after collection, protect the sample from the action of light.
Preservation: Serum bilirubin is stable for 7 days at 4 to 8 ° C and 6 months at -20 ° C.

PRODUCT DESCRIPTION

R	BD	Sulfamic acid ≥ 10mmol / L, hydrochloric acid ≥ 0.5% v / v.	
R	BT	Sulfamic acid ≥ 10mmol / L, hydrochloric acid ≥ 0.1% v / v, Dimethylsulfoxide ≥ 30% v / v.	
R	NI	Sodium nitrite ≥ 10 mmol / L, sucrose ≥ 2 mol / L, preservative.	

STABILITY IN USE

- Once opened, the product (RBD, RBT and RNI) in use is stable to the expiration date printed on the label, since recommended storage conditions (2 to 8 ° C) are followed.
- Reagents should remain outside the specified temperature only for the time required to perform the tests.

TREATMENT AND HANDLING OF THE PRODUCT
A) REAGENT PREPARATION
 RBD, RBT e RNI: Ready-to-use reagents.
B) OPERATING INTERVAL
 The operating range of Direct Bilirubin is from 0.17 mg / dL to 15.0 mg / dL. The operational range of Total Bilirubin is from 0.18 mg / dL to 15.0 mg / dL.
 For higher values, dilute the sample with 150 mM NaCl (0.9%), carry out a new dosage and multiply the result obtained by the dilution factor.

QUALITY CONTROL
 The use of controls should be routine practice in the laboratory. It is suggested to use one control in the reference range or at the decision level and another control with value in another range of clinical significance. For Calibration and Internal Laboratorial Quality Control it is recommended to use the calibrator serum and the control sera below:
 Serum Calibrator - Autocal H 13.002.00
 Normal Control Serum- Quantinorm 13.003.00
 Serum Pathological Control - Quantialt 13.004.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION
A) DIRECT BILIRRUBIN
 1 - Pipette into test tube.

	Blank	Blank Sample	Sample
RBD	1,0 mL	1,0 mL	1,0 mL
RNI	---	---	30 µL
Sample/Calibrator	---	50 µL	50 µL

2. Homogenize well and keep the tubes for 5 minutes at room temperature.
 3. Measure the absorbance of the Calibrator, Calibrator Blank, Sample and Sample Blank against white at 550 nm (530-570). Color is stable for 5 minutes.

B) TOTAL BILIRRUBIN
 1. Pipette into test tubes:

	Blank	Blank Sample	Sample
RBT	1,0 mL	1,0 mL	1,0 mL
RNI	---	---	30 µL
Sample/Calibrator	---	50 µL	50 µL

2. Homogenize well and keep the tubes for 5 minutes at room temperature.
 3. Measure the absorbance of the Calibrator, Calibrator White, Sample and Sample Blank against white at 550 nm (530-570). Color is stable for 30 minutes.

C) CALBRATION
 Use Calibrator Autocal H Ref.13.002.00. The concentration of Bilirubin in the calibrator is traceable to the Standard Reference Material (SRM) 916a of the National Institute of Standards and Technology (NIST).

D) CALCULATIONS
 Bilirubin (mg / dL) = $\frac{(\text{Sample Abs} - \text{Sample White Abs})}{(\text{Calibrator Abs} - \text{White Calibrator Abs})} \times \text{Calibrator Value (mg / dL)}$

Example:
 Calibrator Concentration = 4,62 mg / dL
 Sample Absorbance = 0,194
 Sample Blank Absorbance = 0,040
 Calibrator Absorbance = 0,234
 Absorbance Blank Calibrator = 0,024
 Bilirubin (mg / dL) = $\frac{0,194 - 0,040}{0,234 - 0,024} \times 4,62 = 3,39$ mg / dL

With Calibration Factor
 Calibration Factor = $\frac{\text{Calibrator Concentration (mg / dL)}}{(\text{Calibrator Abs} - \text{White Calibrator Abs})}$

Bilirubin (mg / dL) = (Sample Abs - Abs Sample White) x Calibration Factor

Example:
 Calibration Factor = $\frac{4,62}{0,234 - 0,024} = 22$

Bilirubin (mg / dL) = (0,194-0,040) x 22 = 3,39 mg / dL

C) INTERPRETATION
 Bilirubin is a pigment derived from heme degradation. Bilirubin is transported to the liver, poorly associated with albumin, in its unconjugated form. Carried into the hepatocyte it is rapidly conjugated to produce the bilirubin glucuronides, which are then excreted into the bile. Elevations of conjugated (direct) bilirubin are usually specific markers of hepatic dysfunction, acute hepatitis, or cholestasis. Increased unconjugated (indirect) bilirubin is usually caused by hemolytic processes or by a reduction in liver metabolism.