## Biofluol

Biofluol | Biofluol  
Ref. 37.001.00

### FINALIDADE

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### ENGLISH

#### INTENDED USE

Solution for the cleaning of automatic and semi-automatic biochemical analyzers.

#### STORAGE AND HANDLING

- Store at 15 to 30 °C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

#### WORKING PRINCIPLE

The solution acts by removing wastes remaining in pipes, hoses, and continuous flow cuvettes after use.

#### PRODUCT DESCRIPTION

<b>R 1</b>	Buffer Carbonate, stabilizer, detergent and preservative	<b>X</b>
------------	--	----------

#### STABILITY IN USE

- After opening, the product (R1) in use is stable for 3 months, provided that it is followed by recommended storage conditions (15 to 30 °C).
- The stability of the working reagent is 1 week, provided the recommended conditions of preparation and storage are followed (15 to 30 °C).

#### TREATMENT AND HANDLING OF THE PRODUCT

#### REAGENT PREPARATION

**Working Reagent:** Mix in the proportion of: 1 part of R1 + 3 parts of purified water. Homogenize gently.

#### TEST PROCEDURE

##### A) SEMI-AUTOMATIC BIOCHEMICAL ANALYZERS

Aspirate the working reagent according to the device of the equipment, leaving it in the cuvet during the time established in the equipment manual. Then wash thoroughly with purified water.

##### B) AUTOMATIC BIOCHEMICAL ANALYZERS

Use the working reagent as specified in the equipment manual.

### RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix reagents from different lots.
- Do not use the reagent when it has a visual characteristic that is not in accordance with the product's MSDS requirements.
- Avoid leaving the reagents out of the specified storage conditions when they are not in use.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the used material are fundamental factors for the stability of the reagents and obtaining correct results.

### ALERTS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Disposal, Safety and First Aid information are described in the Individual Product Safety Data Sheet (MSDS) for this product, available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or by telephone +55 35 3214 4646.
- Dispose of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS).

### QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: [sac@biotechnicaltda.com.br](mailto:sac@biotechnicaltda.com.br)

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