C4 PLUS

C4 PLUS /C4 PLUS

Ref. 20.023.00

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ENGLISH

INTENDED USE

Kit for the determination of Complement C4 in serum or plasma. In vitro diagnostic use.

STORAGE AND HANDLING

- Store at 2 to 8 °C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE

Immunoturbidimetric method in which anti-human C4 complement antibodies form insoluble complexes when mixed with samples containing Complement C4. The turbidity intensity is proportional to the C4 Complement concentration of the sample.

SAMPLES – TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION

Sample Type: Serum and plasma (EDTA or heparin).

Collection, handling and preparation: Perform the sample collection according to Good Clinical Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.

Preservation: Complement C4 in serum and plasma, is stable for 8 days stored at 4 to 8 °C and 3 months at -20 °C. Do not use lipemic sample or with hemolysis.

PRODUCT DESCRIPTION

R 1 Anti-complement goat antibodies C4; Tris buffer 20 mmol/L pH 8.2; preservative

STABILITY IN USE

- Once opened, the product R1 in use is stable to the expiration date printed on the label, since recommended storage conditions (2 to 8 °C) are followed.
- Reagents should remain outside the specified temperature only for the time required to perform the tests.

TREATMENT AND HANDLING OF THE PRODUCT

A) REAGENT PREPARATION

R1 : Ready-to-use reagent.

B) OPERATING INTERVAL

The operating range of the product is from 1.092 mg / dL to 93.50 mg / dL.

For higher values, dilute the sample with 150 mM NaCl (0.9%), carry out a new dosage and multiply the result obtained by the dilution factor.

QUALITY CONTROL

The use of controls should be routine practice in the laboratory. For Calibration and Internal Quality Control Laboratory it is recommended to use the calibrator and control below:

Plus Multicalibrador	REF	21.011.00
Medium Multiparameter Control		21.003.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

REACTION PARAMETERS	
TYPE OF REACTION	Final Point/Growing
Wave length	340 nm
Temperature	37 °C
Volume of sample/CAL*	5 µL
Volume of R1*	200 µL
Incubation time	2 minutes
Calibration	Calibration Curve**

* Sample volumes / CAL and R1 may be modified since established ratio isn't maintained.

**CALIBRATION CURVE:

Dilute the Plus Multicalibrador in NaCl 9g / L as follows:

Dilution	1	2	3	4	5	6
Calibrador (µL)	-	10	25	50	75	100
NaCl 9g/L (µL)	100	90	75	50	25	-
Fator	0	0,1	0,25	0,5	0,75	1

Multiply the concentration of Complement C4 indicated in the instruction of use of the Plus Multicalibrador by the corresponding factor to obtain the concentration of Complement C4 of each dilution.

B) CALCULATIONS

The equipment automatically calculates the analyte concentration in each sample.

C) INTERPRETATION

C4 is a component of the essential complement for pathway activation of this system. It is a protein synthesized by the liver and to a lesser extent by monocytes and other tissues. Increases and decreases in concentration are of great clinical significance. Increased levels of C4 concentration are related to acute phase responses (trauma, inflammation). Decreases in concentration are related to deficiencies of genetic origin (autoimmune or collagen vascular diseases, particularly Systemic Lupus Erythematosus), or acquired deficiency as a consequence of the formation of immune complexes, autoimmune hemolytic anemia, and autoimmune.

INTERFERING AND LIMITATIONS

Hemolysis, Jaundice and Lipemia: Hemoglobin > 180 mg / dL / Bilirubin > 20 mg / dL and Lipemia interfere with the assay.

Medications: see recommended literature reference (Young, 2000).

PERFORMANCE CHARACTERISTICS

Sensitivity: Limit of detection: 0.036 mg / dL / Limit of quantification: 1,092 mg / dL.

Prozone High Dose Effect: Prozone effect does not occur up to 130 mg / dL Complement C4 .

Analytical Specificity: The product specifically determines Complement C4 in the presence of other interfering substances in the sample up to the concentrations reported above.

Accuracy: The method was compared with a similar method by the determination of 40 samples in duplicate. The regression equation $y = 0,978x + 0,562$ and correlation coefficient $r = 0,9915$ was obtained. Using this equation the estimated total systematic error of 1.546% for a level of 15.00 mg / dL and 1.497% for a level of 80.00 mg / dL.

Accuracy: It was determined using samples at 03 decimal levels, two runs in duplicate per day for 20 days, being obtained:

Samples (mg/dL)	Repetitions	In-run Precision		Total Precision	
		SD (mg/dL)	%CV	SD (mg/dL)	%CV
17,494	80	0,167	1,0	0,182	1,0
43,762	80	0,635	1,5	0,687	1,6
73,241	80	1,315	1,8	1,896	2,6

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange reagent bottle caps to avoid cross-contamination.
- Do not use the reagent when it has a visual characteristic that is not in accordance with the product's MSDS requirements.
- Avoid leaving the reagents out of the specified storage conditions when they are not in use.
- Use specific glass pipettes and disposable tips for each sample, control, standard / calibrator and reagent.
- The water level of the water bath should be higher than that of the test tubes containing the reactions.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the used material are fundamental factors for the stability of the reagents and obtaining correct results.

REFERENCE RANGES

Adults (20 to 60 years)	10 - 40 mg/dL
-------------------------	---------------

These values are for guidance only and it is recommended that each laboratory establishes its own reference range.

MATERIALS REQUIRED TO CARRY OUT THE TEST

- Spectrophotometer or photometer for reading in 340 nm.
- Glass and / or automatic pipettes.
- Test tubes.
- Biotécnica Plus Multicalibrador - CAT BT 21.011.00

ALERTS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Disposal, Safety and First Aid information are described in the Individual Product Safety Data Sheet (MSDS) for this product, available at www.biotecnica.ind.br or by telephone (35) -3214-4646.
- Dispose of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS).

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use Biotécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiry date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@biotecnicaltda.com.br

AUTOMATION

This product is compatible to the most types of biochemical automatic analysers. The applications are available at www.biotecnicaltda.ind.br

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- SAHU A, LAMBRIS, J.D., *Immunology Reviews*. v. 180, p. 35-48, 2001.
- PRICE, C. P.; et al. *Ann Clin Biochem*. v. 20, p. 1-14, 1983.
- DATI, F. et al; *Eur J Clin Chem Clin Biochem*. v. 34, p. 517-520, 1996.
- TIEZ, Textbook of Clinical Chemistry. 3 ed. Burtis CA, Ashwood ER. **WB Saunders** Co, 1999.
- FRIEDMAN AND YOUNG, Effects of the disease on clinical laboratory tests, 3 ed., **AACC Press**, 1997.
- YOUNG, D.S. *Effects of drugs on clinical laboratory tests - vol. 2*, 5 ed. Washington DC: **AACC Press**, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. *Clin. Chem.* v.27 p.493-501, 1981.

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