**CKNAC**  
CKNAC / CKNAC  
BT1100200 / FD3805000 / SD3905900

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**ENGLISH**

**BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.**

**INTENDED USE**


Kit intended to determine the activity of total Creatine Kinase in serum samples. Diagnostic use only.

**STORAGE AND HANDLING**

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- Do not use reagents whose shelf life has expired.
- **Work Reagent:** mix in the proportion of 4 parts of R1 + 1 part of R2. Homogenize gently. The work reagent is stable for 14 days, as long as the preparation instructions and the recommended storage conditions are followed (2 to 8 °C).
- The reagents can be used in a bi-reactant mode, as long as the established proportions are followed. Consult the Scientific Advisory.
- The Dedicated and Semi Dedicated Reagents presentations have onboard stability of 28 days.

**WORKING PRINCIPLE**

**Method:** Kinetic UV  
Creatine Kinase (CK) catalyzes the phosphorylation of ADP by creatine-phosphate to form creatine and ATP. The enzymatic activity is determined by coupling the reactions of Hexokinase (HK) and Glucose-6-phosphate Dehydrogenase (G6P-DH) and is proportional to the rate of formation of NADPH, which is spectrophotometrically measured at 340 nm.

**SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY**

**Sample Type:** serum

**Collection and handling:** collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

**Preservation:**

	Temperature	Sample Stability
Serum	4 to 8 °C	1 month
	20 a 25 °C	4 hours
	-20 °C	1 month

**PRODUCT DESCRIPTION**

**R 1** Imidazole buffer ≥ 50 mmol/L, NADP ≥ 1 mmol/L, D-glucose ≥ 20 mmol/L, N-Acetyl Cysteine ≥ 20 mmol/L, hexokinase ≥ 2000 U/L, activators, preservative.

**R 2** Creatine phosphate ≥ 100 mmol/L, AMP ≥ 5 mmol/L, Di-adenosina-5-pentaphosphate ≥ 20 μmol/L, ADP ≥ 5 mmol/L, glucose-6-phosphate dehydrogenase ≥ 4000 U/L, activators, preservatives.

The determination of CKNAC activity is traceable to the reference method IFCC.

**QUALITY CONTROL**

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and controls below:

Calibrador - Autocal H	<b>REF</b>	BT1300200
Normal Control – Quantinorm		BT1300300
Pathological Control – Quantialt		BT1300400

**NECESSARY EQUIPMENT FOR TESTING**

- Spectrophotometer or photometer for reading at 340 nm.
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

**TEST PROCEDURE, CALCULATION AND INTERPRETATION**

**A) TEST PROCEDURE**

1. Heat the work reagent at 37 °C for 3 minutes.
2. Zero the reading equipment with purified water filter at 340 nm.
3. Pipette in the assay tube:

	Sample
Sample	40 μL
Work Reagent	1.0 mL

4. Homogenize and insert it immediately in the thermostated cuvette at 37 °C.
5. After 2 minutes, measure the initial absorbance A0 and read it again every minute, for 3 minutes (A1, A2 and A3, respectively).

**B) CALCULATIONS**

Using the measured absorbances, calculate the mean variation of absorbance per minute (ΔA/min):

$$\Delta A/\text{min} = \frac{(A1 - A0) + (A2 - A1) + (A3 - A2)}{3}$$

The CKNAC activity is calculate by multiplying the ΔA/min by the factor below:

$$\text{CKNAC (U/L)} = \Delta A/\text{min} \times 4127$$

**Automation:** this product is compatible to most types of automatic analyzers. Instrument settings are available at [www.biotechnicaltda.ind.br](http://www.biotechnicaltda.ind.br)

Note: to calculate the Factor for automatic analyzers use the absorption coefficient of NADPH = 6.3.

**C) INTERPRETATION**

CK is a dimeric enzyme that occurs in four different forms: a mitochondrial isoenzyme and the cytosolic isoenzymes CK-MM (skeletal muscle), CK-BB (brain) and CK-MB (myocardial). The determination of CK activity and its isoenzymes are used in the diagnosis and monitoring of myocardial infarction and myopathies such as Duchenne muscular dystrophy. In neurogenic muscular diseases there is no change of serum enzyme activity. It is very high in malignant hyperthermia.

**PERFORMANCE CHARACTERISTICS**

Operating range
8.35 to 1000 U/L

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
4.43 U/L	8.35 U/L

Analytical Specificity		
Hemoglobin	Bilirubin	Triglycerides
200 mg/dL	30 mg/dL	2000 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	40 in duplicate
Regression Equation	$y = 1.020x - 3.48$
Correlation Coefficient (R)	0.9951

By applying the regression equation, the total systematic error estimated for decision levels of 40.0 U/L and 180.0 U/L was -6.70% and 0.06%, respectively.

#### Precision:

Determined with two runs in duplicate per day for 20 days.

Samples (U/L)	Replicates	Within-Run Precision		Total Precision	
		SD (U/L)	CV (%)	SD (U/L)	CV (%)
74.39	80	1.756	2.40	2.279	3.10
176.18	80	1.502	0.90	2.302	1.30
269.53	80	5.864	2.20	7.050	2.60

CV: Coefficient of variation; SD: Standard deviation

#### RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- The laboratory should establish chemical, microbiological and particulates requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
- The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

#### REFERENCE RANGES

	37 °C
Women	< 170 U/L
Men	< 190 U/L

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the International System of Units (µKat/L):  
CKNAC (U/L) x 0.017 = CKNAC (µKat/L)

#### WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or calling +55 35 3214 4646

#### QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to [sac@biotechnica.ind.br](mailto:sac@biotechnica.ind.br).
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at [sac@biotechnica.ind.br](mailto:sac@biotechnica.ind.br).

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