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**ENGLISH**

**INTENDED USE**



Kit intended to determine chloride in serum and plasma samples. Veterinary Diagnostic use only.

**STORAGE AND HANDLING**

- Store at 15 to 30 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagent ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (15 to 30°C) are followed.
- Do not use reagents whose shelf life has expired.

**WORKING PRINCIPLE**

**Method:** Mercury Thiocyanate Chloride ions from sample react with mercury thiocyanate to form mercuric chloride and free thiocyanate ions. The latter react with ferric ions to form ferric thiocyanate, whose red color can be spectrophotometrically measured at 500 nm. Color intensity is proportional to the chloride concentration in the sample.

**SAMPLE: TYPE, COLLECTION, HANDLING, PREPARATION AND STABILITY**

**Sample Type:** serum and heparin plasma.

**Collection and handling:** collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

**Precision:**

Serum and plasma	Temperature	Sample Stability
	4 to 8 °C	7 days
	20 a 25 °C	7 days
	-20 °C	1 year

\* Serum samples must be separated within 1 hour after collection.

**PRODUCT DESCRIPTION**

R 1	Mercury thiocyanate $\geq$ 1.0 mmol/L, ferric nitrate $\geq$ 20 mmol/L, nitric acid $\geq$ 0.1% v/v, sulfuric acid $\geq$ 0.1% v/v, mercury nitrate $\geq$ 0.001% v/v.
STD	Chloride ions at a concentration equivalent to 100 mEq/L; preservative. Traceable to the certified reference materials NIST SRM 3182 and NIST SRM 999b.

**QUALITY CONTROL**

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and controls below:

Autocal VET	90.039.00
Normal Control – Quantinorm VET	90.040.00
Pathological Control – Quantialt VET	90.041.00

**NECESSARY EQUIPMENT FOR TESTING**

- Spectrophotometer or photometer for reading at 500 nm (420 – 520 nm).
- Glass pipettes and/or automatic, clock or chronometer.

**TEST PROCEDURE, CALCULATION AND INTERPRETATION**

**A) TEST PROCEDURE**

1. Pipette in the test tubes:

	Blank	Standard	Sample
STD	-	5 $\mu$ L	-
Sample	-	-	5 $\mu$ L
R1	1,0 mL	1,0 mL	1,0 mL

2. Homogenize and incubate at 15 to 30 °C for 5 minutes.

3. Measure the Sample's and Standard's absorbance against the blank at 500 nm (420 – 520 nm). The final reaction is stable for 30 minutes.

**B) CALCULATIONS**

Chloride (mEq/L) =  $\frac{\text{Sample Absorbance}}{\text{STD Absorbance}} \times \text{STD Concentration (mEq/L)}$

Calculations with the Calibration Factor:

Calibration Factor =  $\frac{\text{STD Concentration (mEq/L)}}{\text{STD Absorbance}}$

Chloride (mEq/L) = Sample Absorbance x Calibration Factor

**Automation:** this product is compatible to most types of automatic analyzers. Instrument settings are available at [www.biotechnicaltda.ind.br](http://www.biotechnicaltda.ind.br)

**C) INTERPRETATION**

Chloride is the main extracellular anion. Sodium and chloride represent the majority of the osmotically active constituents of plasma. Chloride ions from food are almost completely absorbed in the intestine. It is significantly involved in: maintaining water distribution, osmotic pressure and anion-cation balance in the extracellular fluid (ECF). Determining plasma chloride concentration is useful in the differential diagnosis of acid-base disorders, and is essential for calculating the anion gap. Chloride fluctuations in serum or plasma have little clinical consequence, but serve as signs of a basic disturbance in fluid homeostasis and acid-base balance, and can help differentiate the cause of these disturbances. Hyperchloremia occurs in dehydration, hypernatremia, renal tubular acidosis (RTA), acute renal failure, metabolic acidosis associated with prolonged diarrhea, loss of sodium bicarbonate, diabetes insipidus, states of adrenocorticoid hyperfunction, overload of treatment with saline solutions and in respiratory alkalosis as a result of renal compensation of HCO<sub>3</sub><sup>-</sup> excretion. Hypochloremia is observed in metabolic acidosis caused by increased production or decreased excretion of organic acids (diabetic ketoacidosis and renal failure), persistent gastric secretion, prolonged vomiting, hypochloremic alkalosis and total body chloride depletion with corresponding retention of HCO<sub>3</sub><sup>-</sup>. Other associated conditions include salt-losing nephritis in Addison's disease, hyperaldosteronism, and bromide poisoning. Excess chloride is excreted in urine and eliminated in sweat.

**PERFORMANCE CHARACTERISTICS**

Operating range	
66.87 to 150.00 mEq/L	

For concentrations above the operating range, dilute the sample with purified water, proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity		
Detection Limit	Quantification Limit	
3.05 mEq/L	66.87 mEq/L	

Analytical Specificity			
Species	Hemoglobin	Bilirubin	Triglycerides
Canines	500 mg/dL	30 mg/dL	1000 mg/dL
Felines	500 mg/dL	40 mg/dL	750 mg/dL
Equines	500 mg/dL	40 mg/dL	1000 mg/dL
Bovines	500 mg/dL	40 mg/dL	1000 mg/dL

Concentrations of interfering substances up to the values presented above do not cause significant changes in the results. For drugs, consult the recommended bibliographic reference (Young, 2000).

Accuracy		
Control	Control 1	Control 2
Interval (mEq/L)	79.1 – 105 (mean 92)	90 – 122 (mean 106)
Average (mEq/L)	95.80	107.13
Replicates	10	10
%DIF	4.1	1.1

**Precision:**

Sample Type	Samples (mEq/L)	Replicates	Within-Run Precision		Total Precision	
			SD (mEq/L)	CV (%)	SD (mEq/L)	CV (%)
Serum	86.14	80	0.769	0.90	0.836	1.00
	95.21	80	0.804	0.80	0.852	0.90
	135.92	80	0.815	0.60	0.893	0.70

CV: Coefficient of variation; SD: Standard deviation

**RESIDUAL RISKS, WARNINGS AND PRECAUTIONS**

- Use protective equipment in accordance with the Good Laboratory Practices.
- The laboratory should establish chemical, microbiological and particulates requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
- The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.
- R1 is toxic because it contains Mercury thiocyanate.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- R1 is toxic because it contains Mercury thiocyanate.

**REFERENCE RANGES**

Canines	105 – 115 mEq/L
Felinos	117 – 123 mEq/L
Equines	99 – 109 mEq/L
Bovinos	97 - 111 mEq/L

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

**Conversion to the International System of Units (mmol/L):**

Chloride (mEq/L) x 1.0 = Chloride (mmol/L)

**WARNINGS AND PRECAUTIONS**

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or calling +55 35 3214 4646

**QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE**

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to [sac@biotecnicaltda.com.br](mailto:sac@biotecnicaltda.com.br)

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