Creatinina Mono

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ENGLISH
INTENDED USE
Kit intended to determination of creatinine in serum, plasma and urine. Diagnostic in vitro use.
STORAGE AND HANDLING
• Store at 2 to 8 ° C.
• Protect from light.
• Expiration date that is printed on the packaging label.
• Do not use reagents whose date has expired.

WORKING PRINCIPLE

The creatinine in the sample reacts with the picrate in alkaline medium to give an orange-red complex. The intensity of the formed color is proportional to the creatinine concentration in the sample.

SAMPLE – TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION

Sample Type: Serum, plasma (EDTA and heparin) and urine.

Collection, handling and preparation: Perform the sample collection according to Good Clinical Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.

Preservation: Serum and plasma creatinine is stable for 7 days if stored at 20 to 25 ° C, 7 days if stored at 4 to 8 ° C and 3 months at -20 ° C. Urine creatinine is stable for 2 days if stored at 20 to 25 ° C, 6 days if stored at 4 to 8 ° C and 6 months if stored at -20 ° C. Collect urine for 24 hours. Homogenize the urine, measure the volume and separate a 20 mL sample. Centrifuge for 10 minutes at 3000 rpm. Dilute an aliquot of the centrifuged urine at the ratio of 1:25 with purified water and proceed with the assay. Multiply the result obtained by 25. This dilution should result

in a measurement within the operating range, to higher values perform new dilution by changing the ratio.

PRODUCT DESCRIPTION	
R 1	Picric Acid ≥ 5,0 mmol/L, Sodium Hidroxide ≥ 100 mmol/L, detergents and stabilizers.
STD	Stabilizers, preservative, creatinine concentration equivalent to 2 mg / dL. Traceable to reference material NIST 914a.

STABILITY IN USE

- After opening, the product (R1 and STD) in use is stable to the expiration date printed on the label since recommended storage conditions (2 to 8 ° C) are followed.
- Reagents should remain outside the specified temperature only for the time required to perform the tests.

PRODUCT TREATMENT AND HANDLING

A) REAGENT PREPARATION

R1 and STD: Ready-to-use reagents.

B) OPERATING INTERVAL

Serum / Plasma: The operating range of the product is from 0.67 mg / dL to 20.00 mg / dL. For higher values, dilute the sample with 150 mL/L NaCl (0.9%), carry out a new dosage and multiply the result obtained by the dilution factor.

Urine: The operating range of the product is 2.99 mg / dL to 342.25 mg / dL corrected for dilution factor.

QUALITY CONTROL

The use of controls should be routine practice in the laboratory. It is suggested to use one control in the reference range or at the decision level and another control with value in another range of clinical significance. For Calibration and Internal Quality Control Laboratory it is recommended to use the calibrator serum and the control sera below:

Serum Calibrator - Autocal H	REF	13.002.00
Serum Normal Control - Quantinorm		13.003.00
Serum Pathological Control - Quantialt		13.004.00

TEST PROCEDURE, CALCULATIONS AND INTERPRETATION

A) TEST PROCEDURE

1. Preheat the R1 for 3 minutes at 37 ° C.
2. Zero the equipment at 510nm (500 - 520 nm) with purified water.
3. Pipette into test tubes:

R1	Volume
Sample or STD	100 µL

4. Homogenize and insert into the thermostated cuvette at 37 ° C. Trigger the timer.
5. Record absorbance at 60 seconds (A1) and 180 seconds (A2) at 510 nm of the sample and STD.

B) CALCULATIONS

Serum / Plasma:

Creatinine (mg / dL) = $\frac{(A_2 - A_1 \text{ of the Sample})}{(A_2 - A_1 \text{ of the STD})} \times \text{STD Concentration (mg / dL)}$

With Calibration Factor:

Calibration Factor = $\frac{\text{STD Concentration (mg / dL)}}{(A_2 - A_1 \text{ of the STD})}$

Creatinine (mg / dL) = (A₂ - A₁ of the Sample) x Calibration Factor

Calculation for Urine:

*Determine the creatinine concentration in the urine sample using the same calculation procedure for serum. Multiply the value obtained by 25 to obtain the result in mg / dL. This should be used to calculate the creatinine value in mg / 24h, according to the following equation:

Creatinine (mg / 24 hours) = $\frac{\text{mg / dL} \times \text{urinary volume (in mL)}}{100}$

Creatinine in mg / kg / 24 hours:

Creatinine (mg / kg / 24 hours) = $\frac{\text{Creatinine (mg / 24 hours)}}{\text{Weight (kg)}}$

Endogenous Creatinine Clearance:

Debug (mL / min) = $\frac{U \times VM}{S}$

Where:

U = urine creatinine (mg / dL)

S = serum creatinine (mg / dL)

VM = minute volume (24 hour urine volume in mL divided by 1440 minutes).

Attention: The clearance should be corrected to the patient's body surface.

Calculation of the Body Surface (SC) without use of the Nomogram:

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Urine: Hemoglobina 500 mg / dL / Bilirrubina 6 mg / dL.

Medications: consult recommended literature reference (Young, 2000).

PERFORMANCE CHARACTERISTICS

Sensitivity:

Serum / Plasma: Detection limit: 0,37 mg / dL / Limit of quantification: 0,67 mg / dL.

Urine: Limit of detection: 0,032 mg / dL / Limit of quantification: 2,9 mg / dL, corrected by dilution factor.

Analytical Specificity: The product specifically determines creatinine in the presence of other interfering substances in the sample, up to the concentrations reported above.

Accuracy:

Serum: The method was compared to a similar method by the determination of 40 samples in duplicate. The regression equation $y = 0,992x - 0,088$ and correlation coefficient $r = 0,9966$ was obtained. Using this equation the estimated total systematic error of -5.2% for a level of 2.00 mg / dL and -1.9% for a level of 8.00 mg / dL.

Urine: The method was compared to a similar method by the determination of 20 samples in duplicate. The regression equation $y = 1,007x - 11,172$ and correlation coefficient $r = 0,9996$ was obtained. Using this equation the estimated total methodological error of -0.541% for a level of 900 mg / 24h and 0.079% for a level of 1800 mg / 24h.

Precision:

Serum / Plasma: It was determined using samples at 03 decision levels, two runs in duplicate per day for 20 days, obtaining:

Samples (mg/dL)	Repetitions	Intra-run accuracy		Total accuracy	
		SD (mg/dL)	%CV	SD (mg/dL)	%CV
0,9494	80	0,0090	1,0	0,0121	1,3
4,2231	80	0,0440	1,0	0,0796	1,9
10,8520	80	0,0707	0,7	0,1366	1,3

% CV: Coefficient of variation expressed as a percentage; SD: Standard deviation.

Urine: It was determined using samples at 02 decision levels, two runs in duplicate per day for 20 days, being obtained:

Samples (mg/dL)	Repetitions	Intra-run accuracy		Total accuracy	
		SD (mg/dL)	%CV	SD (mg/dL)	%CV
81,331	80	0,646	0,8	1,049	1,3
181,728	80	1,052	0,6	1,315	0,7

% CV: Coefficient of variation expressed as a percentage; SD: Standard deviation.

RESIDUAL RISKS, CARE AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix reagents from different batches.
- Do not exchange reagent bottle caps to avoid cross-contamination.
- Do not use the reagent when it has a visual characteristic that is not in accordance with the product's MSDS requirements.
- Use specific glass pipettes and disposable tips for each sample, control, standard / calibrator, and reagent.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the used material are fundamental factors for the stability of the reagents and obtaining correct results.

REFERENCE INTERVAL

Serum / Plasma	Men		Women	
	mg/dL	µmol/L	mg/dL	µmol/L
Urine	0,9 - 1,3	80 - 115	0,6 - 1,1	53 - 97
	mg/kg/day	µmol/kg/day	mg/kg/day	µmol/kg/day
	14 - 26	124 - 230	11 - 20	97 - 177

Clearance (mL/min)	
Adults	52,1 – 110

These values are for guidance only and it is recommended that each laboratory establishes its own reference range.

Conversion to International System Unit (SI): µmol / L

Creatinine (mg / dL) x 88.4 = Creatinine (µmol / L)

MATERIALS REQUIRED TO CARRY OUT THE TEST

- Spectrophotometer or photometer for reading at 510 nm (500 - 520).
- Water bath, thermostated at 37 °C.
- Glass and / or automatic pipettes.
- Clock or Stopwatch.
- Test tubes

ALERTS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Disposal, Safety and First Aid information are described in the Individual Product Safety Data Sheet (MSDS) for this product, available at www.biotechnica.ind.br or by telephone (35) -3214-4646.
- Dispose reactions leftover in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS).

QUALITY / SAC GUARANTEE - CUSTOMER SERVICE

Before being released for consumption, all Biotechnique reagents are tested by the Quality Control Department. The quality of the reagents is ensured by the expiry date mentioned on the presentation packaging, provided that they are stored and transported under specified conditions. The data related to the Quality Control of this product (lot printed on the labels of the reagent bottles) or any doubt on the use of this kit, contact the Scientific Advisory Board of Biotechnica Ltda, by phone +55 35 3214 4646 or by email sac@biotechnicaltda.com.br

AUTOMATION

This procedure is automated in most analyzers.

The protocols are available at www.biotechnica.ind.br

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