Fator Reumatóide Turbilátex

Rheumatoid Factor Turbilátex | Factor Reumatóide Turbilátex
Ref. 20.007.00

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ENGLISH

INTENDED USE

Kit intended to determination of the Rheumatoid Factor (RF) in serum.

PRINCIPLE OF THE METHOD

Latex particles covered with Immunoglobulin G are agglutinated by the rheumatoid factors (RF) present in the sample. The agglutination produces increment of the absorbance and is proportional to the RF concentration in the sample and may be determined by comparison with a calibrator of known concentration at 650nm.

REAGENT COMPOSITION

R DIL	TRIS buffer 20 mmol/L pH 8.2, preservatives	
R LAX	Latex particles covered with Human Immunoglobulin G, pH 7.4, preservatives	
CAL	Human serum containing rheumatoid factors; preservatives	

All the components of human origin show negative results to HBs antigen, and anti-HIV. However, it may be treated as potentially infectious material.

STORAGE AND STABILITY

- After first opening the reagents are stable until the expiration date printed on the label if stored 2-8 °C tightly closed
- Do not use reagents over the expiration date.
- Store at 2 to 8 °C.
- Reagents must remain outside the specified temperature only the time required for testing.
- Do not freeze and protect from light.

ADDITIONAL EQUIPMENT

- Spectrophotometer or photometer with a thermostated cuvet able to read at 650 (600-650) nm.
- Thermostatic water bath at 37 °C.
- Pipettes and micropipettes.
- Clock or chronometer.
- Assay tubes.

WARNINGS AND PRECAUTIONS

- Perform the sample collection according to Good Laboratory Practices
- The kit is intended only for in vitro diagnostic use.
- The samples must be considered as potentially infective.
- Wear the individual protection equipment accordingly to Good Laboratory Practice.
- Discard the reactions surplus according to Good Laboratory Practice in a proper place for potentially infective material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, also is available on the site www.biotecnica.ind.br or calling for 55-35-3214-4646.
- Do not combine reagents from different batches
- Do not exchange the caps from different reagents in order to avoid cross contamination.
- Do not use the reagent if it displays any sign in disagreement with the one specified in MSDS of the product
- Use glass pipettes and disposable tips for each sample, control, standard/calibrator and reagent.
- The level of water in the thermostatic bath must be over the level of mix reaction within the assay tubes.
- The laboratories must establish the chemical requirements, microbiological and particulate matter for water prior to and define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
- The appropriate cleaning and drying of the lab material are fundamental factors for the stability of the reagents to reach correct results.

SAMPLE – PREPARATION AND STABILITY

Serum.

Stable 7 days between 2 to 8 °C or 3 months at -20 °C.

Discard lipemic and/or hemolyzed samples.

Samples containing fibrin must be centrifuged for 10 minutes at 3,000 rpm. Separate the supernatant and carry out the test.

TECHNICAL PROCEDURE

A) REAGENT PREPARATION

The reagents (RDIL e RLAX) are ready to use.

Do not freeze the reagent. It must affect its stability and turn it unavailable for use.

RF Calibrator

The calibrator is lyophilized.

- Let the calibrator at room temperature (10 - 30°C) for 10 minutes.
- Reconstitute with 2.0 mL of purified water.
- Dissolve with gentle agitation. Let rest for 10 minutes, avoiding foam.
- Mix gently by inversion until the complete dissolution. The calibrator is stable for 30 days at 2 to 8 °C and 3 months if stored at -20 °C

B) PROCEDURE

B1 – CURVE CALIBRATION

Proceed with calibrator dilution of rheumatoid factor with saline solution (NaCl 0.9%), as follows

Dilution	1	2	3	4	5	6
Calibrator (µL)	-	25	50	100	200	400
NaCl 9g/L (µL)	400	375	350	300	200	-
Factor	0	0,0625	0,125	0,25	0,5	1,0

For the concentration of each dilution of rheumatoid factor, multiplying the concentration of the calibrator by the corresponding factor indicated in the table.

B2 – TEST

- Heat the reagents and the equipment cuvet at 37 °C.
- Adjust the blank (zero) of the equipment with purified water.
- Pipette in the assay tube

	Blank	CAL/Sample
RDIL (mL)	0,8	0,8
RLAX (mL)	0,2	0,2
NaCl 0,9% (µL)	7,0	-
CAL ou amostra (µL)	-	7,0

- Homogenize and insert in the thermostated cuvet at 37 °C.
- Measure the Absorbance at 650nm after 2 minutes (600-650), Blank (Ab) and CAL (Ac) and Sample (Aa).

CALIBRATION/CALCULATION PROCEDURES

- Calculate the difference of the absorbencies $A_c - A_b$ of each point of the calibration curve.
- Construct the calibration curve by plotting in the X axis the concentration of each calibrator dilution and in the Y axis value of the absorbance difference obtained respectively.
- Calculate the absorbance differences ($A_a - A_b$) of each sample.
- Determine the sample concentrations through interpolation of the obtained values in the calibration curve

SENSITIVITY, LINEARITY AND PROZONE

Analytical sensitivity: 7,58 IU/mL.
Linearity: 160 UI/mL.
For values higher than 150 UI/mL, dilute the sample with NaCl solution 150 mM (0.9%), repeat the assay and multiply the obtained result by the dilution factor.
Prozone: The reaction has no prozone effect < 800 UI/mL.

TECHNICAL LIMITATIONS

Anticoagulants:
Anticoagulants are not recommended, they may interfere in the test.
Interferences: Bilirubin > 20 mg/dL ; Hemoglobin > 10 g/L; Triglycerides > 1000 mg/dL

QUALITY CONTROL

The use of controls to evaluate the imprecision and inaccuracy of the results must be a routine practice in the lab. It's suggested to use a control within the reference range and other control within the clinical significance range.

For the internal quality control of the lab it's indicated the use of the control sera below:

RF Calibrator (included in the kit)	20.007.00
Rheumatic Control Serum N1	21.007.00
Rheumatic Control Serum N2	21.008.00

REFERENCE RANGES

Normal Values: up to 20 IU/mL
These values are intended for orientation only. It's recommended that each lab establish its own reference range.

PERFORMANCE CHARACTERISTICS

Intra-Assay
The realization of 20 determinations of the same sample at the same day showed a Coefficient of Variation of 3,58%.

Inter-Assay
The realization of 10 determinations of the same sample at different days showed a Coefficient of Variation of 3,69%.

Analytical Specificity

A comparison with a reference method showed a correlation coefficient (r) of 0.9881 obtained from ambulatory samples. The result equation of the linear regression is $Y = 0,945 X - 1,21$.

CLINICAL SIGNIFICANCE

Rheumatoid factors are a group of antibodies that specifically bind to the Fc portion of Immunoglobulin G. Even being present in a large number of rheumatic disorders such as SLE and Sjögren's syndrome, its main clinical interest is associated with the diagnosis Rheumatoid arthritis (AR)

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@biotecnicaltda.com.br

AUTOMATION

This product is compatible to the most types of biochemical automatic analysers. The applications are available at www.biotecnicaltda.ind.br

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