Fosfatase Alcalina – IFCC VET

Alkaline Phosphatase – IFCC Vet /
Fosfatasa Alcalina – IFCC VET

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| Espécies | Intervalo de Referência |
|----------|-------------------------|
| Caninos | 20,0 – 156,0 U/L |
| Felinos | 25,0 – 93,0 U/L |
| Equinos | 143,0 – 395,0 U/L |
| Bovinos | 0,0 – 488,0 U/L |

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ENGLISH

INTENDED USE

Kit for the quantitative determination of the enzymatic activity of Alkaline Phosphatase in serum and plasma. Veterinary Diagnostic use only.


STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- Do not use reagents whose shelf life has expired.
- Reagents ready for use. The reagents can also be used as a Work Reagent.
- Work Ready:** mix in the proportion of 1 part of R1 + 1 part of R2. Homogenize gently. The work reagent is stable for 14 days, as long as the preparation instructions and the recommended storage conditions are followed (2 to 8 °C).

WORKING PRINCIPLE

Method: Kinetic / adapted IFCC

The determination is based on a kinetic and colorimetric method, in which the sample alkaline phosphatase, in an alkaline medium, hydrolyzes the p-nitrophenylphosphate present in the reagent, releasing p-nitrophenol and inorganic phosphate. The concentration of p-nitrophenol, measured at 405 nm, is directly proportional to the enzymatic activity of alkaline phosphatase in the sample.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum and heparin plasma.

Collection, handling and preparation: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

| | Temperature | Stability Period |
|------------------|-------------|------------------|
| Serum and Plasma | 4 a 8 °C | 7 days |
| | -20 °C | 2 months |

PRODUCT DESCRIPTION

| | | |
|------------|--|----------|
| R 1 | Buffer, ativadores, 2-Amino 2-Methyl-1-Propanol < 800 mmol/L | |
| R 2 | p-Nitrophenylphosphate < 50 mmol/L, preservative | X |

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratorial quality control, it's recommended the use of the calibrator and controls below:

| | | |
|--------------------------------------|------------|-----------|
| Calibrator - Autocal Vet | REF | 90.039.00 |
| Normal Control – Quantinorm Vet | | 90.040.00 |
| Pathological Control – Quantialt Vet | | 90.041.00 |

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 405 nm (390 – 425 nm).
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Pipette in the test tubes:

| | Bi-reactant protocol | | | Monoreagent Protocol | | |
|----------------------------------|----------------------|------------|--------|----------------------|------------|---------|
| | Blank | Calibrator | Sample | Blank | Calibrator | Sample |
| R1 | 500 µL | 500 µL | 500 µL | RT | 1000 µL | 1000 µL |
| Incubate at 37°C for 180 seconds | | | | | | |
| R2 | 500 µL | 500 µL | 500 µL | CAL | - | 25 µL |
| CAL | - | 25 µL | - | Sample | - | 25 µL |
| Sample | - | - | 25 µL | | | |

- Homogenize and incubate the reaction tubes at 37 °C for 54 seconds;
- Immediately after incubation, measure the initial absorbance (A₀) of each reaction at 405 nm;
- Take new readings every 60 seconds (A₁ and A₂).

B) CALCULATIONS

Alkaline phosphatase (U/L) = Calibration Factor * [(ΔA_{sample}/min) – (ΔA_{blank}/min)]

In what:

$$\Delta A / \text{min} = \frac{(A_1 - A_0) + (A_2 - A_1)}{2}$$

$$\text{Calibration Factor} = \frac{\text{Calibrator Concentration (U/L)}}{\Delta A_{\text{calibrador}} / \text{min} - \Delta A_{\text{blank}} / \text{min}}$$

Automation: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotecnicaltda.ind.br

C) INTERPRETATION

Alkaline Phosphatase is present in the body in membranes and cell surfaces located in the mucosa of the small intestine and in the proximal convoluted tubules of the kidney, in bone (osteoblasts), in the liver and in the placenta, being important for the detection of hepatobiliary diseases. Cats have a lower amount of the enzyme, and it is quickly eliminated by the kidneys. In dogs, liver disease, which causes an increase in alkaline phosphatase, courses with cholestasis. Corticosteroid induction may increase AF levels. In horses and ruminants, due to the wide reference ranges of FA in these species, its variation is of little importance in liver diseases. It may also be increased in cases of osteomalacia, fracture healing, vitamin D deficiency, hyperparathyroidism, bone tumor, rickets, hyperadrenocorticoidism, pregnancy, and retained placenta. In addition to corticosteroids, drugs such as barbiturates, cephalosporins, phenobarbital, phenothiazines, phenylbutazone, tetracyclines, tiabendazole and halothane induce an increase in alkaline phosphatase. AF of bone origin may be increased in young animals, or in cases of hyperparathyroidism, osteosarcoma, in fracture healing, osteomalacia or vitamin D deficiency. Castrated animals have greater enzyme activity than intact animals.

PERFORMANCE CHARACTERISTICS

| Operating range |
|----------------------|
| 2.45 to 1823.56 U/L |
| Sensitivity |
| Quantification Limit |

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0,9%), proceed with a new dosage and multiply the result by the dilution factor.

2.45 U/L

| Analytical Specificity | | | |
|------------------------|------------|-----------|---------------|
| Species | Hemoglobin | Bilirubin | Triglycerides |
| Canines | 300 mg/dL | 40 mg/dL | 750 mg/dL |
| Felinos | 500 mg/dL | 40 mg/dL | 750 mg/dL |
| Equinos | 200 mg/dL | 30 mg/dL | 2000 mg/dL |
| Bovinos | 500 mg/dL | 40 mg/dL | 2000 mg/dL |

Interfering substances up to the values presented above do not cause significant alterations in results. For drugs, consult the recommended reference (Young, 2000).

| Accuracy - Canines | |
|-----------------------------|---------------------|
| Number of Samples | 40 in duplicate |
| Regression Equation | y = 0.998x - 0.255 |
| Correlation Coefficient (R) | 1.0000 |
| Accuracy - Felinos | |
| Number of Samples | 20 in duplicate |
| Regression Equation | y = 0.990x + 0.705 |
| Correlation Coefficient (R) | 0.9999 |
| Accuracy - Equinos | |
| Number of Samples | 11 in duplicate |
| Regression Equation | y = 1.002x - 3.087 |
| Correlation Coefficient (R) | 0.9972 |
| Accuracy - Bovinos | |
| Number of Samples | 10 in duplicate |
| Regression Equation | y = 1.032x - 0.4286 |
| Correlation Coefficient (R) | 0.9993 |

Precision:

Determined with two runs and in 20 replicates

| Samples (U/L) | Number of samples | Within-Run Precision | | Samples (U/L) | Number of samples | Between-Run Precision | |
|---------------|-------------------|----------------------|--------|---------------|-------------------|-----------------------|--------|
| | | SD (U/L) | CV (%) | | | SD (U/L) | CV (%) |
| 146.24 | 20 | 6.79 | 4.60 | 143.40 | 40 | 6.70 | 4.70 |
| 482.19 | 20 | 12.34 | 2.60 | 482.48 | 40 | 11.98 | 2.50 |

%CV: Coefficient of Variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

| Species | 37 °C |
|---------|-------------------|
| Canine | 20.0 – 156.0 U/L |
| Feline | 25.0 – 93.0 U/L |
| Equine | 143.0 – 395.0 U/L |
| Bovine | 0.0 – 488.0 U/L |

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to Good Laboratory Practice in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotechnica.ind.br or calling for +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnicaltda.com.br.

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| | |
|---------|-------------------|
| Caninos | 20,0 – 156,0 U/L |
| Felinos | 25,0 – 93,0 U/L |
| Equinos | 143,0 – 395,0 U/L |
| Bovinos | 0,0 – 488,0 U/L |

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- YOUNG, D.S. *Effects of drugs on clinical laboratory tests* - vol. 2, 5 ed. Washington DC: AACCC Press, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. *Clin. Chem.* v.27 p.493-501, 1981.
- KANEKO, Jiro J; HARVEY, John W; BRUSS, Michael L. *Clinical biochemistry of domestic animals*. 5th. ed. San Diego: Academic Press, c1997 932p.
- THRALL, Mary Anna. *Hematologia e bioquímica clínica veterinária*. São Paulo: Roca, 2007. 582 p.

TABELA DE SÍMBOLOS INTERNACIONAIS / TABLE OF INTERNATIONAL SYMBOLS / TABLA DE SÍMBOLOS INTERNACIONALES

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