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ENGLISH

INTENDED USE
 Kit for the quantitative determination of the enzymatic activity of γ -glutamyl transferase (GT Gamma) in serum and plasma. Use *in vitro* diagnostic.

STORAGE AND HANDLING

- Store at 2 to 8 °C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE

The GT Gamma catalyzes the transfer of the glutamyl group of L- α -glutamyl-3-carboxy-4-nitroanilide to glycyglycine, originating L- α -glutamylglycyglycine and 5-amino-2-nitrobenzoate. The catalytic concentration of the enzyme present in the sample is determined from the speed of formation of 5-amino-2-nitrobenzoate, measured at 405nm.

SAMPLES: TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION

Sample Type: Serum and Plasma (EDTA).
Collection, handling and preparation: Perform the sample collection according to Good Clinical Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.
Preservation: The GT Gamma in serum and plasma is stable for 30 days if stored at 4 to 8 °C and 1 year at -20 °C.

PRODUCT DESCRIPTION

R 1	TRIS buffer \geq 50 mmol/L; Glycylglycine \geq 50 mmol / L; preservative.	
R 2	L- γ -Glutamyl-3-carboxy-4-nitroaniline \geq 10 mmol / L; MÉS \geq 1 mmol / L; preservative.	

Traceable to the specific absorbance of the chromogen 5-amino-2-nitrobenzoate using calibrated pipettes and handheld spectrophotometer that provides absolute values, based on modified methodology from IFCC (International Federation of Clinical Chemistry and Laboratory Medicine).

STABILITY IN USE

- The stability of the product (R1 and R2) in use is 15 months, since that recommended storage conditions are followed (2 to 8 °C).
- The stability of the working reagent is 3 weeks, since under recommended conditions of storage and preparation (2 to 8 °C).
- Reagents should remain outside the specified temperature only for the time required to perform the tests.

TREATMENT AND HANDLING OF THE PRODUCT

A) REAGENT PREPARATION
Working Reagent (WR)
 Mix in the ratio: 4 parts of R1 + 1 part of R2. Homogenize gently.

B) OPERATING INTERVAL
 The operating range of the product is 4.00 U/L to 270 U/L.
 For higher values, dilute the sample with 150 mM NaCl (0.9%), carry out a new dosage and multiply the result obtained by the dilution factor.

QUALITY CONTROL

The use of controls should be routine practice in the laboratory. It is suggested to use one control in the reference range or at the decision level and another control with value in another range of clinical significance. For Calibration and Internal Quality Control Laboratory it is recommended to use the calibrator serum and the control sera below:

Calibrator Serum - Autocal H	13.002.00
Normal Control Serum - Quantinorm	13.003.00
Pathological Control Serum - Quantialt	13.004.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Preheat the Working Reagent for 3 minutes at 37 °C.
2. Zero the equipment at 405nm with purified water.
3. Pipette into test tube:

	Volume
Work Reagent	1.0 mL
Sample	100 μ L

4. Homogenize and insert it in the thermostated cuvette. Start the timer.
5. After 60 seconds, note the initial absorbance (A₀) and read again after exactly 1, 2 and 3 minutes (A₁, A₂ and A₃, respectively).

B) CALCULATIONS
 Using the absorbance readings, calculate the mean absorbance change per minute ($\Delta A/\text{min}$):

$$\Delta A/\text{min} = \frac{(A_1 - A_0) + (A_2 - A_1) + (A_3 - A_2)}{3}$$

The catalytic activity of the GT Gama in the sample is calculated by multiplying the $\Delta A/\text{min}$ by the corresponding factor.

$$\text{GT Gama (U/L)} = \Delta A/\text{min} \times 1158$$

Example:
 A₀ = 1,160 A₁ = 1,189
 A₂ = 1,220 A₃ = 1,250

$$\Delta A/\text{min} = \frac{(1,189 - 1,160) + (1,220 - 1,189) + (1,250 - 1,220)}{3}$$

$\Delta A/\text{min} = 0,030$
 GT Gama (U/L) = 0,030 x 1158 = 34,7 U/L

C) INTERPRETATION
 The GT Gamma is present in the proximal renal tubule, liver, pancreas and intestine. The enzyme present in the serum originates primarily from the hepatobiliary system, however, increased GT Gamma activities occur concomitantly with the elevation of other biochemical parameters. Therefore, its analysis together with alkaline phosphatase, transaminases and bilirubins, is useful in the differential diagnosis of primary and secondary liver diseases. It is elevated in chronic alcoholics and individuals who make continuous use of anticonvulsants, such as phenytoin and phenobarbital.

INTERFERING AND LIMITATIONS

Hemolysis, Jaundice and Lipemia: Hemoglobin > 200 mg/dL / Bilirubin > 40 mg/dL / Triglycerides > 900 mg / dL.
Medications: consult recommended literature reference (Young, 2000).

PERFORMANCE CHARACTERISTICS

Sensitivity: Detection limit: 3.46 U/L / Limit of quantification: 4.00 U/L.