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
ENGLISH
INTENDED USE
Kit intended to determination of glucose in serum and plasma. Diagnostic use only.

STORAGE AND HANDLING

- Store at 2 to 8 °C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE
The glucose present in the sample is oxidized by the enzyme glucose oxidase forming gluconic acid and hydrogen peroxide. It oxidizes the phenolic reactant in the presence of 4-aminoantipyrine, under the catalytic action of peroxidase, producing a pink compound (quinoneimine), which exhibits a maximum absorption at 505 nm. The color intensity is proportional to the glucose concentration in the sample.

SAMPLE – TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION
Sample Type: Serum and plasma (fluoreto).
Collection, handling and preparation: Make the collection of the sample in accordance with Good Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material. Serum samples should be separated from the cell elements immediately.
Preservation: The glucose is stable in serum for 3 days at a temperature of 4 to 8 °C, 8 hours at a temperature of 20 to 25 °C and 24 hours at a temperature of -20 °C. Plasma glucose is stable for 3 days at a temperature of 20 to 25 °C.

PRODUCT DESCRIPTION

R 1	Phosphate buffer ≥ 40 mmol/L; 4-aminoantipyrine ≥ 0,2 mmol/L; phenol ≥ 1,0 mmol/L; Peroxidase ≥ 500 U/L; Glucose Oxidase ≥ 5000 U/L; activators; detergents and stabilizer.	X
STD	Glucose 100mg/dL; stabilizer and preservative. Traceable to NIST reference material 917c.	X

STABILITY IN USE

- The stability of the product (R1 e STD) in use is 18 months as long as followed by the recommended storage conditions (2 to 8 °C).
- Reagents must remain outside the specified temperature only the time required for testing

TECHNICAL PROCEDURE
A) REAGENT PREPARATION
Reagent 1 (R1)
Ready-to-use Reagent

B) OPERATING RANGE
The product operating range is from 3,05 mg/dL to 400 mg/dL. For higher values, dilute the sample with NaCl 150 mM (0,9%), repeat the assay and multiply the result by the dilution factor.

QUALITY CONTROL
Use of controls should be a routine practice in the lab. It is suggested to use a control within the reference range or decision level and another control with value in another range of clinical significance. For the internal quality control of the lab it's indicated the use of the calibrator serum and control sera below:
Calibrator Serum - Autocal H **REF** 13.002.00
Normal Control Serum – Quantinorm 13.003.00
Pathological Control Serum – Quantialt 13.004.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE
1. Pipette in the assay tubes:

	Blank	Standard	Sample
STD	-	10 µL	-
Sample	-	-	10 µL
R1	1,0 mL	1,0 mL	1,0 mL

2. Homogenize and incubate 10 minutes at 37 °C.
3. Measure the absorbance of Standard and Sample against the blank at 505 nm. The final reaction is stable for 30 minutes.

B) CALCULATIONS
Glucose (mg/dL) = $\frac{\text{Sample Absorbance}}{\text{Standard Absorbance}} \times \text{STD (mg/dL)}$

Example:
STD concentration = 100 mg / dL
Sample Absorbance = 0,219
Absorbance of STD = 0,302
Glucose (mg / dL) = $\frac{0,219}{0,302} \times 100 = 72,5$ mg / dL

With Calibration Factor (CF):
CF = $\frac{\text{STD (mg/dL)}}{\text{STD Absorbance}}$
Glucose (mg/dL) = Sample Absorbance x CF
Calibration Factor = $\frac{100}{0,302} = 331$
Glucose (mg/dL) = 0,219 x 331 = 72,5 mg/dL

C) INTERPRETATION
Glucose is the predominant carbohydrate in the peripheral blood, being the main source of cellular energy, stored under the forms of glycogen in the liver or fatty acids in adipose tissue. Their amount in the blood is controlled in narrow limits by several hormones, which have the function of reducing (insulin) or increasing (glucagon, adrenalina, cortisol and hormone of the growth). Hyperglycemia occurs in Diabetes mellitus caused by deficiency in the secretion or action of insulin. Other secondary factors can also elevate blood glucose levels, which include: pancreatitis, thyroid