HDL Colesterol Direto

HDL - Cholesterol (Direct) |

HDL Colesterol Directo

Ref. 10.006.00

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ENGLISH

INTENDED USE

Kit for the determination of HDL cholesterol in serum and plasma. In vitro diagnostic use.

STORAGE AND HANDLING

- Store at 2 to 8 °C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

PRINCIPIO DE FUNCIONAMENTO

In the first stage of the reaction, cholesterol oxidase catalyzes the decomposition of non-HDL cholesterol generating hydrogen peroxide that is consumed by peroxidase with DSBmT. In the second step, HDL cholesterol is solubilized by a specific detergent and quantified in the presence of enzymes and a chromogenic reagent. The intensity of the color produced is proportional to the amount of HDL cholesterol present in the sample.

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Sample Type: Serum or plasma (heparin and EDTA).

Collection, handling and preparation: Perform the sample collection according to Good Clinical Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.

Preservation: HDL in serum and plasma is stable for 7 days at a temperature of 4-8 ° C and 2 months at a temperature of -20 ° C.

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R 1	Buffer; Cholesterol oxidase <1000 U/L; Peroxidase <1300 ppg U/L; N-bis (4-sulfobutyl)-m-toluidine-dissodium N-bis (DSBmT) <1 mM; Acelerador <1 mM; Preservante <0.06%; Ascorbate oxidase <3000 U/L.	X
R 2	Buffer; Cholesterol-esterase <1500 U / L; 4-aminoantipirine (4-AAP) <1mM; Detergent <2%; Preservative.	X

CAL	Human serum, preservatives. Value printed on the bottle label. Traceable to CDC reference method for HDL Cholesterol.	X
------------	---	----------

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- **R1 and R2:** After opened, the reagents in use are stable until the expiration date printed on the label, SINCE recommended storage conditions are followed (2 to 8 ° C).
- **CAL:** After reconstitution, CAL is stable for 14 days if stored at 2 to 8 ° C and 3 months at -20 ° C.
- Reagents should remain outside the specified temperature only for the time required to perform the tests.

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R1 and R2: Ready-to-use reagents.

CAL: Freeze-dried. Reconstitute with 1 mL of purified water. Leave to stand for 30 minutes. Homogenize gently by inversion avoiding foaming.

B) INTERVALO OPERACIONAL

The operating range of the product is from 4.211 mg / dL to 200.0 mg / dL.

For higher values, dilute the sample with 150 mM NaCl (0.9%), carry out a new dosage and multiply the result obtained by the dilution factor.

QUALITY CONTROL

The use of controls should be routine practice in the laboratory. It is suggested to use one control in the reference range or at the decision level and another control with value in another range of clinical significance. For Calibration and Internal Quality Control Laboratory it is recommended to use the calibrator serum and control sera below:

HDL Cholesterol Direct CAL (included in kit) Normal	10.006.00
Control Serum- Quantinorm	13.003.00
Serum Pathological Control - Quantial	13.004.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

REACTION PARAMETERS	
Types of reaction	Final/Growing point
Bichromatic Wave length	600 nm/700nm
Temperature	37 °C
Purified water (white) /sample/CAL*	3 µL
R1*	300 µL
Reading 1 (A ₁)	5 minutes
R2*	100 µL
Reading 2 (A ₂)	5 minutes
Calibration	Linear

* Volumes can be modified as long as the established ratio is maintained.

B) CALCULATIONS

HDL Cholesterol (mg/dL) = $\frac{\Delta A_{\text{Sample}} - \Delta A_{\text{White}}}{\Delta A_{\text{Calibrator}} - \Delta A_{\text{White}}} \times \text{CAL Value (mg/dL)}$

Where:

$\Delta A = A_2 - A_1$

C) INTERPRETATION

High-density lipoproteins (HDL) are responsible for the reverse transport of cholesterol from peripheral tissues to the liver. In the liver, cholesterol is transformed into bile acids that are excreted in the gut via the bile duct. Monitoring of HDL cholesterol is of great clinical importance, as there is an inverse correlation between serum HDL concentration and the risk of atherosclerotic disease. Reduced concentrations of HDL cholesterol, along with high triglycerides, increase cardiovascular risk.

INTERFERING AND LIMITATIONS

Hemolysis, Jaundice and Lipemia: Hemoglobin > 500 mg / dL / Bilirubin > 40 mg / dL / Triglycerides > 2000 mg / dL, interfere with the dosage.

Medications: consult recommended literature reference (Young, 2000).

PERFORMANCE CHARACTERISTICS

Sensitivity: Limit of detection: 0.714 mg / dL / Limit of quantification: 4.211 mg / dL.

Analytical Specificity: The product specifically determines HDL cholesterol in the presence of other interfering substances in the sample up to the concentrations reported above.

Accuracy: The method was compared with a similar method by the determination of 40 samples in duplicate. The regression equation $Y = 0,998x + 0,30$ and correlation coefficient $r = 0,9991$ were obtained. Using this equation the estimated total systematic error of 1,00% for a level of 25 mg / dL and 0,05% for a level of 120 mg / dL.

Accuracy: It was determined using samples at 03 decision levels, two runs in duplicate per day for 20 days, being obtained.

Samples (mg/dL)	Repetitions	In-run Precision		Total Precision	
		SD (mg/dL)	%CV	SD (mg/dL)	%CV
61,630	80	0,407	0,70	0,477	0,80

101,479	80	0,327	0,30	0,555	0,50
171,738	80	0,377	0,20	0,424	0,20

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange reagent bottle caps to avoid cross-contamination.
- Do not use the reagent when it has a visual characteristic that is not in accordance with the product's MSDS requirements.
- Avoid leaving the reagents out of the specified storage conditions when they are not in use.
- Use specific glass pipettes and disposable tips for each sample, control, standard / calibrator and reagent.
- The water level of the water bath should be higher than that of the test tubes containing the reactions.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the used material are fundamental factors for the stability of the reagents and obtaining correct results.

REFERENCE RANGES

Age Range	Desirable reference values (mg/dL)
Adults (> 20 years)	> 40
Children and teenagers	> 45

NOTE: Reference values are desirable according to the Brazilian Consensus for the Normalization of Laboratory Determination of the Lipid Profile.

MATERIALS REQUIRED TO CARRY OUT THE TEST

- Spectrophotometer or photometer for reading at 600 nm/700nm.
- Glass and / or automatic pipettes.
- Clock or Stopwatch.

ALERTS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Disposal, Safety and First Aid information are described in the Individual Product Safety Data Sheet (MSDS) for this product, available at www.biotechnica.ind.br or by telephone (35) -3214-4646.
- Dispose of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS).

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use, all BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@biotechnicaltda.com.br

AUTOMATION

This product is compatible to the most types of biochemical automatic analysers. The applications are available at www.biotechnicaltda.ind.br

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- BADIMON JJ, BADIMON L, FUESTER V. Regression of atherosclerotic lesions by high density lipoprotein plasma fraction in the cholesterol-fed rabbit. *Journal of Clinical Investigation*. 85; 1234-41 (1990).
- BARR DP, RUSS EM, EDER HA. Protein-lipid relationships in human plasma. *Am. J. Med.*, 11;480-493 (1951).
- CASTELLI WP et al. HDL Cholesterol and other lipids in coronary heart disease. *Circulation ahaajournals.org*. 55;767-772 (1977).
- CROUSE JR et al. Studies of low density lipoprotein molecular weight in human beings with coronary artery disease. *J. Lipid Res.* 26, 566-574 (1985).
- GORDON T et al. High density lipoprotein as a protective factor against coronary heart disease. *Am. J. Med.*, 62; 707-714 (1977).
- FALUDI AA, IZAR MCO, SARAVA JF, CHACRA APM, BIANCO HT, AFUNSE NETO A et al. Atualização da Diretriz Brasileira de Dislipidemias e Prevenção da Aterosclerose - 2017. *Arq Bras Cardiol* 2017; 109(2Supl.1):1-76.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. *Clin. Chem.* v.27 p.493-501, 1981.
- YOUNG, D. S. *Effects of drugs on clinical laboratory tests - vol. 2*, 5 ed. Washington DC: AACC Press, 2000.
- WARNICK, GR; WOOD, PD. National cholesterol education program recommendations for measurement of high-density lipoprotein cholesterol: executive summary. *Clin. Chem.* V.41. 10, 1427-1433 (1995).

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