HbA1C Plus	Responsável Técnico: Dr. Gilson Sérgio Pizzo CRF MG - 5310
HbA1C Plus HbA1C Plus Ref. 20.025.00 (KIT) 20.025.00 (KIT COM CONTROL)	MS 80027310273

FINALIDADE

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ENGLISH
INTENDED USE Kit intended for the quantitative determination of human hemoglobin A1c (HbA1c) in whole blood. Use in <i>in vitro</i> diagnostic.
STORAGE AND HANDLING <ul style="list-style-type: none">• Store at 2 to 8 ° C.• Protect from light.• Do not freeze.• Stable until the kit expiration date that is printed on the label.• Do not use reagents whose date has expired.
WORKING PRINCIPLE Immunoturbidimetric method for the direct determination of HbA1c in whole blood. Total hemoglobin and HbA1c have the same non-specific absorption rate as the latex particles present in R1. When the human anti-HbA1c mouse monoclonal antibody is added (R2), the latex-HbA1c-human anti-HbA1c monoclonal antibody complex is formed. Agglutination occurs when the goat polyclonal antibody interacts with the mouse monoclonal antibody. The turbidity intensity is proportional to the concentration of HbA1c in the sample.
SAMPLES: TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION Sample Type: Total Blood (EDTA). Collection, handling and preparation: Make the collection of the sample in accordance with Good Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.
Preservation: HbA1c in whole blood is stable 7 days if stored at 4 to 8 ° C.

CONTROL DESCRIPTION		
R 1	Latex 0.13%, buffer, stabilizer, preservative.	
R 2	Mouse anti-human HbA1c monoclonal antibody 0.05 mg / ml, goat anti-mouse 0.08mg / dl polyclonal IgG antibody, buffer, stabilizers, preservative.	
R HML	Purified Water, preservative.	
CONTROL	Whole human blood, glucose and preservative. Traceable to the National Glycohemoglobin Standardization Program (NGSP). (Included in the Kit with CONTROL)	

STABILITY IN USE

- R1, R2, RHML: Once opened, reagents in use are stable to the expiration date printed on the label, provided that they are followed at recommended storage conditions (2 to 8 ° C).
- CONTROL after reconstitution is stable for 30 days under the recommended storage conditions (2 to 8 ° C).
- Reagents should remain outside the specified temperature only for the time required to perform the tests

TECHNICAL PROCEDURE

A) REAGENT PREPARATION

R1, R2 and RHML: Ready for use. Gently homogenize before use.

CONTROL: Lyophilized. Kepp at 10-30 ° C for 10 minutes. Reconstitute with 0.5 mL of purified water. Gently homogenize by inversion avoiding foaming until complete dissolution of the control.

B) OPERATING INTERVAL

The operating range of the product is 3% to 15%.

QUALITY CONTROL

The use of controls should be routine practice in the laboratory. For Calibration and Internal Quality Control Laboratory it is recommended to use the calibration set and control below:
HbA1c Calibration Set **REF** 21.009.00
HbA1C Plus CONTROL (Included in the Kit with CONTROL) 20.025.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

HEMOLYSIS

The procedure of hemolysis must be applied on samples, control and calibrators, as described below:

- 1- Hemolyse 20 µL of sample, control and calibrators previously homogenized, with 1.0 mL of RHML, in different test tubes.
 - 2 – Let the tubes rest for 5 minutes to complete hemolysis.
- Hemolyse are stable for up to 10 days if stored at 2 to 8 ° C.

TEST

REACTION PARAMETERS	
Type of reaction	End Point / Ascending
Wave-length	600 - 660 nm
Temperature	37 °C
Sample Volume / CAL *	5 µL
R1* Volume	180 µL
R2* Volume	60 µL
Incubation Time	5 min.
Calibration	Calibration curve**

* Sample volumes, CAL, R1 and R2 can be modified as long as the established ratio is maintained.

****CALIBRATION CURVE:**

Use the Calibration Set HbA1C- CAT BT 21.009.00.

Point 0	Purified water
Point 1	CAL 1
Point 2	CAL 2
Point 3	CAL 3
Point 4	CAL 4

B) CALCULATIONS

The equipment automatically calculates the analyte concentration in each sample.

C) INTERPRETATION

There are several subtypes of HbA1, such as HbA1a1, HbA1a2, HbA1b and HbA1c. Of these, the HbA1c fraction is the glycated hemoglobin itself and corresponds to 80% of the glycated fraction. HbA1c is formed by the binding of glucose to the valine amino acid of the N-terminal portion of the beta chain of hemoglobin by means of a stable and irreversible binding. HbA1c is found in non-diabetic adults. It reflects the average glycemia of the last 2 to 3 months, which corresponds to the half-life of the red blood cells. The higher the glycemia, the higher the HbA1c concentration. The diagnostic criteria for diabetes mellitus underwent several modifications, precisely due to the appearance of new evidence regarding the association of increasingly low blood glucose levels with micro and macrovascular complications, which result in retinopathy, nephropathy, neuropathy, coronary disease, cerebrovascular disease and peripheral arterial disease.

INTERFERING AND LIMITATIONS

Hemolysis, Jaundice and Lipemia: Triglycerides> 2000 mg / dL / Bilirubin> 50 mg / dL interfere with the dosage.

Medications: consult recommended literature reference (Young, 2000).

PERFORMANCE CHARACTERISTICS

Analytical Sensitivity: Variation of 0.087 absorbance for each 1% HbA1c, to a 4.8% level of HbA1c. **Analytical Specificity:** The product specifically determines HbA1c in the presence of other interfering substances in the sample, up to the concentrations reported above.

Accuracy: The method was compared with a similar method by the determination of 34 samples in duplicate. The regression equation Y = 0.996 X + 0.013 and correlation coefficient r = 0.995 was obtained. Using this equation the estimated total systematic error of - 0.16% for a level of 5.5% and - 0.28% for a level of 10.5%.

Precision: It was determined using samples at 03 decision levels, two runs in duplicate per day for 20 days, obtaining:

Sample (%)	Repetitions	In-run Precision		Total Precision	
		SD (%)	%CV	SD (%)	%CV
5,81	80	0,057	1,0	0,068	1,2
6,98	80	0,072	1,0	0,076	1,1
13,3	80	0,107	0,8	0,156	1,2

% CV: Coefficient of variation expressed as a percentage; SD: Standard deviation.

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with Good Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange the caps from different reagents in order to avoid cross contamination.
- Do not use the reagent if it displays any sign in disagreement with the one specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions, when they are not in use.
- Use glass pipettes and disposable tips for each sample, control, standard / calibrator and reagent.
- The laboratory should establish chemical, microbiological and particulates requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
- The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.

REFERENCE RANGES

Diagnostic criteria for diabetes mellitus recommended by the Brazilian Society of Diabetes			
Test	Normal	Pre-diabetes	Diabetes
Fasting glycemia (mg / dL)	< 100	100 a 125	≥ 126
Glucose 2 hours after oral tolerance test with 75 g of glucose (mg / dL)	< 140	140 a 199	≥ 200
HbA1c (%)	< 5,7	5,7 a 6,4	≥ 6

These values are for guidance only and it is recommended that each laboratory establishes its own reference range.

ADDITIONAL EQUIPMENT

- Spectrophotometer or photometer for reading at 600 - 660 nm.
- Glass and / or automatic pipettes.
- Test tubes.
- Biotécnica Calibration Set - CAT BT 21.009.00

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to Good Laboratory Practice in a proper place for potentially infective material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotecnica.ind.br or calling for 55-35-3214-4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use Biotécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@biotecnicaltda.com.br.

AUTOMATION

This product is compatible to the most types of biochemical automatic analysers. The applications are available at www.biotecnicaltda.ind.br

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