## HbA1c Aut Plus

HbA1c Aut Plus | HbA1c Aut Plus  
Ref. 20.031.00 (KIT) 20.031.00 (KIT COM CONTROL)  
/ FD 38.079.00 / SD 39.081.00

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#### ENGLISH

 **BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.**

#### INTENDED USE

Kit intended for the quantitative determination of human glycosylated hemoglobin (A1c) in whole blood. Use in *in vitro* diagnostic.

#### STORAGE AND HANDLING

- Store at 2 to 8 °C and do not freeze. This product must remain out of the specified temperature only the time required for testing. Protect from light.
- Reagents ready for use.
- Once opened, the product in use is stable for 4 weeks when stored at 2 to 8 °C.
- Do not use reagents whose shelf life has expired.

#### WORKING PRINCIPLE

**Method:** Immunoturbidimetric.

The determination of HbA1c in whole blood is based in antibody and antigen interaction. Total hemoglobin and HbA1c have the same non-specific absorption rate as the latex particles present in R1. When the human anti-HbA1c mouse monoclonal antibody is added (R2), the latex-HbA1c-human anti-HbA1c monoclonal antibody complex is formed. The mouse anti-IgG goat polyclonal antibody cause the agglutination of the already formed complexes, altering the turbidity of the medium. The turbidity intensity is proportional to the concentration of HbA1c in the sample and can be measured spectrophotometrically in 660 nm.

#### SAMPLES: TYPE, COLLECTION, HANDLING AND PRESERVATION

**Sample Type:** Whole Blood (EDTA).

**Collection and handling:** collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

#### Preservation:

	Temperature	Stability Period
Whole blood	4 to 8 °C	7 days

#### PRODUCT DESCRIPTION

<b>R 1</b>	Latex 0.1%, buffer.	
<b>R 2</b>	Mouse anti-human HbA1c monoclonal antibody > 0.03 mg/mL, buffer.	
<b>R HML</b>	Buffer, preservative.	
<b>CONTROL 1</b>	Whole human blood and preservative. Traceable to the National Glycohemoglobin Standardization Program (NGSP).	
<b>CONTROL 2</b>	Whole human blood and preservative. Traceable to the National Glycohemoglobin Standardization Program (NGSP).	
<b>CAL</b>	HbA1c Aut Plus CAL 1, HbA1c Aut Plus CAL 2, HbA1c Aut Plus CAL 3, HbA1c Aut Plus CAL 4: Lyophilized whole human blood and preservative. The HbA1c calibrators' concentrations are printed on the labels of the vials. Traceable to the National Glycohemoglobin Standardization Program (NGSP).	

#### QUALITY CONTROL

The use of controls should be routine practice in the laboratory. For Internal Quality Control Laboratory it is recommended to use the calibration set and control below:

HbA1c Aut Plus CAL	21.013.00
HbA1c Aut Plus CONTROL 1	20.031.00
HbA1c Aut Plus CONTROL 2	20.031.00

#### TECHNICAL PROCEDURE

##### A) REAGENT PREPARATION

**R1:** Ready for use. Gently homogenize before use.

**R2:** Ready for use. Gently homogenize before use.

**RHML:** Ready for use. Gently homogenize before use.

**CONTROL 1 e 2:** Lyophilized. Reconstitute with 0.5 mL of purified water. Gently homogenize by inversion avoiding foaming until complete dissolution of the control. Leave the control at a temperature of 10 – 30°C for 15 minutes before use.

**CALIBRATORS 1 - 4:** Lyophilized. Reconstitute with 0.5 mL of purified water. Gently homogenize by inversion avoiding foaming until complete dissolution of the calibrator. Leave the calibrator at a temperature of 10 – 30°C for 15 minutes before use.

#### NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 660 nm.
- Glass and/or automatic pipettes.
- Test tubes.

#### TEST PROCEDURE, CALCULATION AND INTERPRETATION

##### A) TEST PROCEDURE

###### HEMOLYSIS

The procedure of hemolysis must be applied only on samples.

1- Hemolysate 10 µL of sample previously homogenized, with 500 µL of RHML

Hemolysates are stable for up to 7 days if stored at 2 to 8 °C.

###### TEST

Type of reaction	Fixed Time
<b>Primary Wavelength</b>	660 nm
<b>Secondary Wavelength (optional)</b>	800 nm
<b>Temperature</b>	37 °C
<b>Sample/CAL volume*</b>	4 µL
<b>R1 volume*</b>	150 µL
<b>Incubation time of R1 + A</b>	300 seconds
<b>R2 volume*</b>	50 µL
<b>First reading</b>	30 seconds after addition of R2
<b>Second reading</b>	300 seconds after addition of R2
<b>Calibration**</b>	Calibration curve

\* Sample, CAL, R1 and R2 volumes can be modified as long as the established ratio is maintained.

\*\*Calibration Curve:

Use the HbA1c Aut Plus CAL - CAT BT XX.XXX.XX.

Points	Calibrators	Concentration (%)
Point 0	Purified water	0.00%
Point 1	CAL 1	Available on the product label
Point 2	CAL 2	
Point 3	CAL 3	
Point 4	CAL 4	

It is suggested to use a Spline/Exponential mathematical model for the calibration curve.

#### B) CALCULATIONS

The equipment automatically calculates the analyte concentration in each sample.

#### C) INTERPRETATION

Of the different types of glycated hemoglobin - HbA1a1, HbA1a2, HbA1b and HbA1c – HbA1c corresponds to 80% of the glycated fraction, being the glycated hemoglobin itself. This molecule is formed by the addition, in a stable and irreversible manner, of glucose to the valine amino acid of the N-terminal portion of the beta chain of hemoglobin. HbA1c reflects the average glycaemia of the last 2 to 3 months before its dosage, which corresponds to the half-life of the red blood cells. The higher the glycaemia, the higher the HbA1c concentration in blood. The diagnostic criteria for diabetes mellitus underwent several modifications, precisely due to the publication of new evidence regarding the association of increasingly low blood glucose levels with micro and macrovascular complications, which result in retinopathy, nephropathy, neuropathy, coronary disease, cerebrovascular disease and peripheral arterial disease.

#### PERFORMANCE CHARACTERISTICS

<b>Operating range</b>
2.0 up to 16.0 %

<b>Sensitivity</b>
<b>Quantification limit</b>
0.40 %

Analytical Specificity		
Bilirubin	Lipids	Ascorbic Acid
≤ 50 mg/dL	≤ 700 mg/dL	≤ 50 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy was determined by comparison with a similar method by the analysis of 40 samples. The regression equation  $Y = 1.0572 X - 0.4156$  and correlation coefficient  $r = 0.9909$  were obtained.

#### Precision:

Intra-assay precision studies were performed with one batch and 20 replicates in an analytical run; inter-assay precision studies were performed with one batch and two replicates in two analytical runs.

Sample (%)	Repetitions	Within-run Precision		Total Precision	
		SD (%)	%CV	SD (%)	%CV
5.56	40	0.09	1.11	0.07	1.99
11.98	40	0.10	0.59	0.20	1.49

% CV: Coefficient of variation expressed as a percentage; SD: Standard deviation.

#### RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with Good Laboratory Practice.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents, controls and calibrator outside the specified storage conditions.

#### REFERENCE RANGES

Normal	< 5.7%
Pre-diabetes	5.7 – 6.4%
Diabetes	≥ 6%

These values are for orientation only and it is recommended that each laboratory establishes its own reference range.

#### WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to Good Laboratory Practice, in a proper place for potentially infective material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or calling for 55-35-3214-4646

#### QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to [sac@biotechnica.ind.br](mailto:sac@biotechnica.ind.br).
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at [sac@biotechnica.ind.br](mailto:sac@biotechnica.ind.br).

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- Goldstein DE, Little RR, Lorenz RA, et al. Tests of glycemia in diabetes. Diabetes Care 1995;18:896-909.
- The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes on the development and progression of long-term complications in insulin-dependent diabetes mellitus. N Engl J Med 1993;329:977-986.
- American Diabetes Association, "Standards of Medical Care for Patients with Diabetes Mellitus", (Position statement), Diabetes Care, 21 (Suppl.1): S23 S31 (1998).
- Pantheghini M, John WG on behalf of the IFCC Scientific Division. Implementation of hemoglobin A1c results traceable to the IFCC reference system: the way forward. Clin Chem Lab Med 2007; 45(8): 942-4.

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