originates mainly from the skeletal muscle, brain, skin, kidney marrow and erythrocytes. The concentration of blood lactate depends on the rate of production in these tissues and the metabolic rate in the liver and kidneys. The liver uses approximately 65% (75 g/day) of total basal lactate produced predominantly in gluconeogenesis. Anaerobic glycolysis markedly increases lactate in the blood, especially in prolonged exercise. The anoxia resulting from shock, pneumonia and congestive heart failure; Thiamine deficiency and diabetic ketoacidosis, among other causes, cause increased lactate and blood pyruvate. High blood lactate concentration results in lactic acidosis, which can occur in renal failure and leukemia.

INTERFERES OR LIMITATIONS

Hemolysis, Jaundice and Lipemia: Hemoglobin > 500 mg/dL / Bilirubin > 40 mg/dL / Triglycerides > 1000 mg/dL.

Medications: consult recommended reference (Young, 2000).

PERFORMANCE CHARACTERISTICS

Sensitivity: Detection limit: 0.84 mg/dL / Limit of quantification: 2.32 mg/dL.

Analytical Specificity: The product specifically determines lactate in the presence of other interfering substances in the sample up to the concentrations reported above.

Accuracy: The method was compared with a similar method by the determination of 40 samples in duplicate. The regression equation $y = 0.990x - 0.16$ and correlation coefficient $r = 0.9983$ was obtained. Using this equation the estimated total systematic error of -1.32% for a level of 50 mg/dL and -1.16% for a level of 100 mg/dL.

Accuracy: It was determined using samples at 03 decision levels, two runs in duplicate per day for 20 days, obtaining:

Samples (mg/dL)	Repetition s	Within-Run Precision		Total Precision	
		SD (mg/dL)	%CV	SD (mg/dL)	%CV
25,281	80	0,244	1,0	0,568	2,2
35,246	80	0,247	0,7	0,524	1,5
105,411	80	0,230	0,2	0,631	0,6

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange reagent bottle caps to avoid cross-contamination.
- Do not use the reagent when it has a visual characteristic that is not in accordance with the product's MSDS requirements.
- Avoid leaving the reagents out of the specified storage conditions when they are not in use.
- Use specific glass pipettes and disposable tips for each sample, control, standard / calibrator, and reagent.
- The water level of the water bath should be higher than the test tubes containing the reactions.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to its use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the used material are key factors for the stability of the reagents and obtaining correct results.

REFERENCE RANGES

Sample	mg/dL	mmol/L	Specification
Plasma (Fluoride)	4,5-19,8	0,5-2,2	Venoso
	4,5-14,4	0,5-1,6	Arterial

These values are for guidance only and it is recommended that each laboratory establishes its own reference range.

Conversion to International System Unit (SI): Lactate (mg/dL) x 0.111 = Lactate (mmol/L)

MATERIAL NEEDED TO PERFORM THE TEST

- Spectrophotometer or photometer for reading at 546 nm (530 - 570).
- Water bath, thermostated at 37 ° C.
- Glass and / or automatic pipettes.
- Clock or Stopwatch.
- Assay tubes.

WARNINGS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Dispose, Safety and First Aid information are described in the Individual Safety Data Sheet for this product, available at www.biotechnica.ind.br or by telephone (35) -3214-4646.
- Dispose of the leftovers of the reactions according to the Good Clinical Laboratory Practices (GCLP) and Health Service Waste Management Program (HSWMP).

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being released for consumption, all Biotécnica reagents are tested by the Quality Control Department. The quality of the reagents is ensured until the expiration date mentioned on the presentation packaging since they are stored and transported under specified conditions. The data related to the Quality Control of this product (lot printed on the labels of the reagent bottles) or any doubt in the use of this kit, contact the Scientific Advisory of Biotécnica Ltda, by phone +55 35 3214 4646 or by email sac@Biotechnicaltda.com.br

AUTOMATION

This procedure is compatible to the most types of analysers. The applications are available at www.biotechnicaltda.ind.br

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- HUTCHESSON, A. et al. Measurement of lactate in cerebrospinal fluid in investigation of inherited metabolic disease. **Clin. Chem.** v.43, n.1, p.158-161, 1997.
- GUYTON, A. C.; HALL, J. E. **Fisiologia humana e mecanismos das doenças**. 6ª ed. Rio de Janeiro: Guanabara Koogan, 1998. 580p.
- FOSSATI, P.; PRENICEPE, L.; BERT G. Use of 3,5-dichloro-2-hydroxybenzenesulfonic acid/4-aminophenazone chromogenic system in direct enzymic assay of uric acid in serum and urine. **Clin. Chem.** v.26, p.227-231, 1980.
- BARHAM D.; TRINDER P. An improved colour reagent for the determination of blood glucose by oxidase system. **Analyst** v.27, p.142-145, 1972.
- SHINOJO, N. et al. Test strip method for measuring lactate in whole blood. **Clin. Chem.** v.35, n.9, p.1992-1994, 1989.
- BURTIS, C. A.; ASHWOOD, E. R.; BRUNS, D. E. **Tietz Fundamentos de Química Clínica**, Saunders Elsevier, 6 ed, 2008.
- BURTIS, C. A.; ASHWOOD, E. R. **Tietz Fundamentos de Química Clínica**, Guanabara Koogan, 4 ed, 1998.
- YOUNG, D.S. **Effects of drugs on clinical laboratory tests** - vol. 2, 5 ed. Washington DC: AACCP Press, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. **Clin. Chem.** v.27 p.493-501, 1981.
- WESTGARD, James O.; LAHMEYER, Brenda L.; BIRNBAUM, Marvin L. Use of the Du Pont "Automatic Clinical Analyzer" in Direct Determination of Lactic Acid in Plasma Stabilized with Sodium Fluoride. **Clin. Chem., Cincinnati**, v. 18, nº. 11, p.1334-1338, nov. 1972.

TABELA DE SÍMBOLOS INTERNACIONAIS / TABLE OF INTERNATIONAL SYMBOLS / TABLA DE SÍMBOLOS INTERNACIONALES

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