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ENGLISH

INTENDED USE
 Kit for the quantitative determination of lipase enzymatic activity in serum and plasma. *In vitro* diagnostic use.

STORAGE AND HANDLING

- Store at 2 to 8 °C
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired

WORKING PRINCIPLE

Lipase catalyzes the cleavage of the 1-2-O dilauryl-rac-glycerd-3-gliceric- (6-methylresorufin)-ester substrate in alkaline medium forming 1-2-O dilauryl-rac-glycerd and gliceric-6'-methylresorufin ester which spontaneously decomposes into gliceric acid and methylresorufin. The velocity of formation of methylresorufin is proportional to the catalytic activity of the lipase in the sample.

1-2-O dilauryl-rac-glycerd-3-gliceric- (6-methylresorufin)-ester $\xrightarrow{\text{Lípase}}$ 1-2-O dilauryl-rac-glycerd + Gliceric-6'-methylresorufin-ester (unstable) \rightarrow Gliceric Acid + Methylresorufin

SAMPLES – TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION

Sample Type: Serum and plasma (heparin).
Collection, handling and preparation: Perform the sample collection according to Good Clinical Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.

Preservation: Serum and plasma lipase is stable for 7 days at 20-25 °C, 3 weeks at 4-8 °C and 1 year at -20 °C

PRODUCT DESCRIPTION

Buffer >= 50 mmol/L; Cdpase >= 1 mg/L; sodiumtaurodeoxycholate >= 5 mmol/L, sodiumdeoxycholate >= 1 mmol/L, activator, detergent, preservative
 Tartarate buffer >= 2 mmol/L; sodiumtaurodeoxycholate >= 5 mmol/L, Cdpase >= 0,1 mg/L; 1-2-O dilauryl-rac-glycerd-3-gliceric- (6-methylresorufin)-ester >= 2 mmol/L; Stabilizers; activators, detergents, preservative

Traceable to the specific absorbance of the chromogen Methylresorufin using calibrated pipettes and handheld spectrophotometer that provide absolute values.

STABILITY IN USE

- After opening, the product (R1 and R2) in use is stable to the expiration date printed on the label, recommended storage conditions are followed (2 to 8 °C).
- Reagents should remain outside the specified temperature only for the time required to perform the tests.

TREATMENT AND HANDLING OF THE PRODUCT

A) REAGENT PREPARATION

R1 and R2: Ready-to-use reagents.

B) OPERATIONAL INTERVAL

The operating range of the product is 3.78 U/L to 300 U/L. For higher values, dilute the sample with 150 mM NaCl (0.9%), carry out a new dosage and multiply the result obtained by the dilution factor.

QUALITY CONTROL

The use of controls should be routine practice in the laboratory. It is suggested to use one control in the reference range or at the decision level and another control with value in another range of clinical significance. For Calibration and Internal Quality Control Laboratory it is recommended to use the calibrator serum and control sera below

Serum Calibrador - Autocal H	REF	13.002.00
Normal Control Serum Quantinorm		13.003.00
Serum Pathológico Control - Quantidade		13.004.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE
 1- Adjust the temperature of the photometer to 37 °C, select the wavelength at 580nm and read with purified water.
 2- Pipette into test tube.

	Volum
R1	800 μL
Purified Water / Sample / CAL	10 μL
R2	200 μL

- 3- Homogenize and insert into the thermostated cuvette. Set the stopwatch
- 4- Note the absorbance at 90 seconds (A₁) and 180 seconds (A₂), write, sample and calibrator.

B) CALCULATION

Lipase (U/L) = $\frac{\Delta A \text{ Sample} - \Delta A \text{ White}}{\Delta A \text{ Calibrador} - \Delta A \text{ White}}$ x Concentration Calibrator (U/L)

Where:
 ΔA White = A₂ - A₁ of purified water.
 ΔA Sample = A₂ - A₁ of the sample
 ΔA Calibrator = A₂ - A₁ of the CAL

C) INTERPRETATION

Lipase is an enzyme that hydrolyzes esters of glycerol with long chain fatty acids. Most of the enzyme in the serum is derived from the pancreas. Its determination is used to investigate pancreatic disorders, such as acute pancreatitis, chronic pancreatitis, and pancreatic canal obstruction. After an acute pancreatitis event, serum lipase activity increases within 4 to 8 hours, peaking at around 24 hours and decreasing after 8 to 14 days. However, elevation of serum lipase activity is not proportional to the damage caused by the event. Values are also elevated in acute and chronic renal disease and in acute cholecystitis.

INTERFERING AND LIMITATIONS

Hemolysis, Icteric and Lipemia: Hemoglobina > 150 mg/dL / Bilirubin > 40 mg/dL / Triglicérides > 900 mg/dL

Medications: consult recommended literature reference (Young, 2000).