	
<b>Magnésio Mono</b> Magnesium Mono / Magnesium Mono Ref. 12.007.00	
<b>Responsável Técnico:</b> Dr. Gilson Sérgio Pizzo CRF MG – 5310 <b>MS 80027310168</b>	

**FINALIDADE**

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**ENGLISH**

**INTENDED USE**  
 Kit intended to determination of magnesium in serum, plasma and urine. IVD Diagnostic use only.

**STORAGE AND HANDLING CONDITIONS**

- Store at 15 to 30 ° C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

**WORKING PRINCIPLE**

Magnesium forms with sulfonated Magon (xylydil blue I) in alkaline medium a stained complex that can be measured at 505 nm. The intensity of the formed color is proportional to the amount of magnesium present in the sample.

**SAMPLE – TYPE, COLLECTION, HANDLING, PREPARATION AND STABILITY**

**Sample Type:** Serum, plasma (heparina) and urine.  
**Collection, handling and preparation:** Make the collection of the sample in accordance with Good Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.

**Preservation:** The magnesium in serum and plasma (heparina) is stable for 7 days if stored at a temperature of 4 to 8 °C and 6 months if stored at a temperature of - 20 °C . In urine the magnesium is stable for 3 days if stored at a temperature of 4 to 8 °C, acidified in pH < 2,0.

Urine: Use a sample collected within 24 hours in a bottle containing 20 to 30 mL of 50% HCl. Mix well, measure the volume, check the pH and, if necessary, adjust it to approximately 1.0. Dilute an urine aliquot of 1:5 with purified water (urine dilution should be changed, if necessary, to achieve a result within the operating range of the method). Use diluted urine for testing. Multiply the result obtained by 5 (or by the factor that has been applied).

**PRODUCT DESCRIPTION**

<b>R 1</b>	Tris buffer ≥ 100 mmol / L; Sulphonated magon 0.1 mmol / L; Activators; Preservative; Stabilizers.	<b>X</b>
<b>STD</b>	Magnesium ions in concentration equivalent to 2,0 mg/dL; conservative. Traceable to the reference material NIST 3131a.	<b>X</b>

**STABILITY IN USE**

- The stability of the product (R1 and STD) in use is 24 months as long as followed by the recommended storage conditions (15 to 30 °C).
- Reagents should remain outside the specified temperature only for the time required to perform the tests.

**TECHNICAL PROCEDURE**

**A) REAGENT PREPARATION**  
 R1 and STD: reagents are ready to use.

**B) OPERATING RANGE**  
 The product operating range is from 0,40 mg/dL to 4,50 mg/dL. For higher values, dilute the sample with NaCl 150 mM (0.9%), repeat the assay and multiply the result by the dilution factor.

**QUALITY CONTROL**

The use of controls should be routine practice in the laboratory. It is suggested to use one control in the reference range or at the decision level and another control with value in another range of clinical significance. For Calibration and Internal Control of Laboratory Quality it is recommended to use the calibrator serum and the control sera below:

Calibrator Serum - Autocal H		13.002.00
Normal Control Serum – Quantinorm	<b>REF</b>	13.003.00
Pathological Controle Serum – QuantiAlt		13.004.00

**TEST PROCEDURE, CALCULATION AND INTERPRETATION**

**A) TEST PROCEDURE**

1. Pipette in the assay tubes:

	Blank	Standard	Sample
STD	-	10 µL	-
Sample	-	-	10 µL
R1	1.0 mL	1.0 mL	1.0 mL

2. Homogenize and incubate 2 minutes at 15 to 30 °C.
3. Measure the absorbance of Standard (Ast) and Sample (Asa) against the blank at 505 nm. The color is stable for 15 min.

**B) CALCULATIONS**

**Calculation of serum/plasma**  
 Magnesium (mg/dL) =  $\frac{\text{Sample Absorbance}}{\text{Standard Absorbance}} \times \text{STD Concentration}$

**Example:**  
 STD concentration = 2,0 mg/dL  
 Sample absorbance = 0,175  
 Standard absorbance = 0,191  
 Magnesium (mg/dt) =  $\frac{0,175}{0,191} \times 2,0 = 1,83$  mg/dL

**With Calibration Factor (CF):**  
 CF =  $\frac{\text{STD Concentration (mg/dL)}}{\text{STD Absorbance}}$

Magnesium (mg/dL) = Sample Absorbance x CF

**Example:**  
 Calibration Factor =  $\frac{2,0}{0,191} = 10,47$   
 Magnesium (mg/dL) = 0,175 x 10,47 = 1,83 mg/dL

**Urine 24 hours:**  
 Magnesium urine (mg/dL) =  $\frac{\text{Sample absorbance}}{\text{Standard Absorbance}} \times \text{STD concentration (mg/dL)} \times 5^*$

\*5 = Dilution Factor applied

Magnesium urine (mg/24 hours) =  $\frac{\text{Magnesium urine (mg/dL)} \times \text{urinary volume in 24 hours (ml)}}{100}$

**C) INTERPRETATION**

Aproximadamente 55% of the human Mg<sup>2+</sup> cation is found in the skeleton, another part inside tissue cells and a small proportion in the extracellular fluid (about 1%), being the fourth most abundant cation in the body. Magnesium ion acts in more than 300 enzymatic processes as cofactor, allosteric activator and / or in the formation of enzymatic substrates, such as MgATP. Ingested magnesium is absorbed primarily in the small intestine. Its regulation occurs in the kidneys, in the ascending branch of Henle. Hypomagnesemia can be observed in several clinical conditions, such as:

alcoholism, acute pancreatitis, malabsorption syndromes, metabolic acidosis, renal diseases and protein-calorie malnutrition. Hypermagnesemia may occur in some situations, for example: excessive use of antacids, enemas and parenteral fluids containing magnesium, rhabdomyolysis, acute and chronic renal failure. The dosage of urinary magnesium is mainly used to evaluate hydroelectrolytic

#### INTERFERING AND LIMITATIONS

Hemolysis, Jaundice and Lipemia: Hemoglobin> 500 mg/dL / Bilirubin> 40 mg/dL / Triglycerides> 2000 mg/dL interferes with the dosage.

Medications: consult recommended literature reference (Young, 2000).

#### PERFORMANCE CHARACTERISTICS

Sensitivity: Detection limit: 0.07 mg/dL / Limit of quantification: 0.40 mg/dL.

Analytical Specificity: The product specifically determines magnesium in the presence of other interfering substances in the sample, up to the concentrations reported above.

Accuracy: The method was compared with a similar method by the determination of 40 samples in duplicate. The regression equation  $y = 1.019x - 0.034$  and correlation coefficient  $r = 0.9940$  was obtained. Using this equation the estimated total systematic error of 0.20% for a level of 2.00 mg/dL and 0.77% for a level of 3.00 mg/dL.

Accuracy: It was determined using samples at 03 decision levels, two runs in duplicate per day for 20 days, obtaining:

Samples (mg/dL)	Repetitive ns	Within-Run Precision		Total Precision	
		SD (mg/dL)	%CV	SD (mg/dL)	%CV
1,3169	80	0,0084	0,6	0,0090	0,7
3,4550	80	0,0340	1,0	0,0393	1,1
4,2952	80	0,0491	1,1	0,0562	1,3

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

#### RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange reagent bottle caps to avoid cross-contamination.
- Do not use the reagent when it has a visual characteristic that is not in accordance with the product's MSDS requirements.
- Use specific glass pipettes and disposable tips for each sample, control, standard / calibrator, and reagent.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to its use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the used material are key factors for the stability of the reagents and obtaining correct results

#### REFERENCE RANGES

Serum and Plasma	1,6 - 2,6 mg/dL
Urine 24h	72,9 – 121,5 mg/24 hours

These values are intended for orientation only. It's recommended that each lab establishes its own reference range.

Conversion for International System of Units (SI): mmol/L

Magnesium (mg/dL) X 0.25 = Magnesium (mmol/L)

#### ADDITIONAL EQUIPMENT

- Spectrophotometer or photometer for readings at 505 nm (500-510 nm).
- Glass and / or automatic pipettes.
- Clock or stopwatch.
- Test tubes.

#### WARNINGS AND PRECAUTIONS

- Disposal, Safety and First Aid information are described in the Individual Safety Data Sheet for this product, available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or by telephone (35) -3214-4646.
- Dispose of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS).

#### QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being released for consumption, all Biotechnique reagents are tested by the Quality Control Department. The quality of the reagents is ensured until the expiration date mentioned on the presentation packaging once they are stored and transported under specified conditions. The data related to the Quality Control of this product (lot printed on the labels of the reagent bottles) or any doubt in the use of this kit, contact the Scientific Advisory of Biotécnica Ltda, by phone +55 35 3214 4646 or by email [sac@biotecnicalltda.com.br](mailto:sac@biotecnicalltda.com.br).

#### AUTOMATION

This product is compatible to the most types of biochemical automatic analyzers. The applications are available at [www.biotecnicalltda.ind.br](http://www.biotecnicalltda.ind.br)

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