the purification system must be validated and it is important to ensure that the water obtained continues to meet the requirements by means of periodic checks.

- Proper cleaning and drying of the used material are fundamental factors for the stability of the reagents and obtaining correct results.

#### REFERENCE RANGES

Table with 3 columns: Range, Albumin (mg/24hs), Albumin / Creatinine (mg/g). Rows include Normal, Moderately increased, and Severely increased.

These values are for guidance only and it is recommended that each laboratory establishes its own reference range.

#### MATERIALS REQUIRED TO CARRY OUT THE TEST

- Spectrophotometer or photometer for reading at 540 (520 - 560) nm
• Glass and / or automatic pipettes.
• Test tubes

#### ALERTS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Disposal, Safety and First Aid information are described in the Individual Product Safety Data Sheet (MSDS) for this product...
• Disposal of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS).

#### QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging...

#### AUTOMATION

This product is compatible to the most types of biochemical automatic analyzers. The applications are available at www.biotechind.com.br

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Table listing reagents: DIL, RLAX, CAL, CONTROL with their descriptions and hazard symbols.

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- BAR, J. et al., Microalbuminuria: Prognostic and Therapeutic Implications in Diabetic and Hypertensive Pregnancy. Diabetic Medicine. v.12, p. 649-656, 1995.
• CAMBIASSO, C.L. et al., Immunoassay of Low Concentrations of Albumin in Urine by Latex Particle Counting. Clin. Chem. v. 34, p. 416-418, 1988.

- Kidney Disease: Improving Global Outcomes (KDIGO) CKD Work Group. KDIGO clinical practice guideline for the evaluation and management of chronic kidney disease. Kidney Int. Suppl. 2013, 31-150.
• RASMUSSEN, B.F; et al., Microalbuminuria: An Important Diagnostic Tool. J Diab Comp. v. 8; p. 137-145, 1994.
• MEDCALF, E.A; et al., Rapid, Robust Method for Measuring Low Concentrations of Albumin in Urine. Clin. Chem. v. 36, p. 446-449, 1990.

#### TABELA DE SÍMBOLOS INTERNACIONAIS / TABLE OF INTERNATIONAL SYMBOLS / TABLA DE SÍMBOLOS INTERNACIONALES

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