**PCR TURBI LÁTEX**

CRP Turbilátex / PCR Turbilátex  
Ref. 20.015.00

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

**ENGLISH**  
**INTENDED USE**  
Kit intended to determine the C Reactive protein (CRP) in serum

**PRINCIPLE OF THE METHOD**  
Latex particles covered with antibodies anti-human CRP reacts with CRP present in the sample. The agglutination produces increment of the absorbance and is proportional to the CRP concentration in the sample and may be determined at 540nm by comparison with a calibrator of known concentration.

**REAGENT COMPOSITION**

<b>R DIL</b>	TRS buffer 20 mmol/L pH 8,2, preservatives	
<b>R LAX</b>	Latex particles covered with goat IgG anti C Reactive Protein pH 7,3, preservatives	 
<b>CAL</b>	Human serum containing CRP, preservatives	 

All the components of human origin show negative results to Hbs antigen, and anti-HV. However, it may be treated as potentially infective material.

- STORAGE AND STABILITY**
- After first opening the reagents are stable until the expiration date printed on the label if stored 2-8 °C tightly closed
  - Do not use reagents over the expiration date
  - Store at 2 to 8 °C
  - Reagents must remain outside the specified temperature only the time required for testing
  - Do not freeze and protect from light.

**ADDITIONAL EQUIPMENT**

- Spectrophotometer or photometer with thermostatic cuvet able to read at 540 (530-550) nm
- Thermostatic water bath at 37 °C
- Pipettes and micropipettes.
- Clock or chronometer.
- Assay tubes.

- WARNINGS AND PRECAUTIONS**
- Perform the sample collection according to Good Laboratory Practices
  - The kit is intended only for in vitro diagnostic use
  - The samples must be considered as potentially infective.
  - Wear the individual protective equipment accordingly to Good Laboratory Practices
  - Discard the reaction surplus according to Good Laboratory Practices in a proper place for potentially infective material.
  - The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product. The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, also available on the site [www.biotecnica.ind.br](http://www.biotecnica.ind.br) or calling for 55-35-3214-4646
  - Do not combine reagents from different batches
  - Do not exchange the caps from different reagents in order to avoid cross contamination
  - Do not use the reagent if it displays any significant disagreement with the one specified in MSDS of the product
  - Use glass pipettes and disposable tips for each sample, control, standard/calibrator and reagent.
  - The level of water in the thermostatic bath must be over the level of mix reaction with the assay tubes.
  - The laboratories must establish the chemical requirements, microbiological and particulate matter for water prior to and define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
  - The appropriate cleaning and drying of the lab material are fundamental factors for the stability of the reagents to reach correct results.

**SAMPLE - PREPARATION AND STABILITY**  
Fresh serum  
Stable 7 days between 2 to 8 °C or 3 months at -20 °C  
Discard lipemic and/or hemolyzed samples.  
Samples containing fibrin must be centrifuged for 10 min at 3000 rpm before assay.

**TECHNICAL PROCEDURE**  
**A) REAGENT PREPARATION**  
Reagents (RD L, RLAX and CAL) ready for use

**B) PROCEDURE**

- Heat the reagent at 37 °C
- Zero the read equipment with purified water.
- Pipette in the cuvet:

RD L	800 µL
RLAX	200 µL
CAL or SAMPLE	5 µL

- Homogenize and insert in the thermostatic cuvet at 37 °C. Trigger the timer.
- Measure the initial absorbance (A1) and 120 seconds after (A2) at 540nm (530-550).

**CALCULATIONS AND PROCEDURES**  
CRP (ng/L) =  $\frac{(A2 - A1) \text{ Sample} \times \text{Calibrator concentration (ng/L)}}{(A2 - A1) \text{ Calibrator}}$

Example  
A1 - sample = 1,825  
A2 - sample = 1,877  
A1 - calibrator = 1,896  
A2 - calibrator = 1,990  
Calibrator concentration = 55 ng/L  
CRP =  $\frac{(1,877 - 1,825) \times 55}{(1,990 - 1,896)} = 0,052 \times 55$   
0,094  
CRP = 30,4 ng/L

**SENSITIVITY AND LINEARITY**  
Analytical sensitivity: 1,92 ng/L  
Linearity: 150 ng/L  
For values higher than 150 U/mL, dilute the sample with NaCl solution 150 mM (0,9%), repeat the assay and multiply the obtained result by the dilution factor.  
Prozone: The reaction has no prozone effect < 800 ng/L

**TECHNICAL LIMITATIONS**

- Anti-coagulants: They cause measurement errors.
- Hemolyzed, Jaundiced and lipemic Sera
- Bilirrubin > 20 mg/dL / Hemoglobina > 5g/L / Lipids > 10 g/L
- Others  
Rheumatoid Factors > 800 U/mL

## QUALITY CONTROL

Any clinical laboratory must keep an internal quality control program which defines the aids, procedures and criteria for the tolerance limits, corrective actions and registration of the activities. Asq, it must be kept a defined system for verification of the analytical variability that occurs in any measuring system.

The use of controls to evaluate the imprecision of the analysis must be a routine practice in the lab. The lab must participate of external quality control programs.

For the internal quality control of the lab it's indicated the use of the control sera below

CRP CAL (induded)	20.015.00
Rheumatic Control Serum N1	21.007.00
Rheumatic Control Serum N2	21.008.00

REF

## REFERENCE RANGES

Normal Values: up to 6 mg/L

These values are intended for orientation only. It's recommended that each lab establish its own reference range.

## PERFORMANCE CHARACTERISTICS

### Intra-Assay

The realization of 20 determinations of the same sample at the same day showed a Coefficient of Variation of 2,75%

### Inter-Assay

The realization of 10 determinations of the same sample at different days showed a Coefficient of Variation of 2,81%

### Analytical Specificity

A comparison with a reference method showed a correlation coefficient (r) of 0,95 obtained from a random sample. The result equation of the linear regression is  $Y = 0,955 X - 0,513$ .

## CLINICAL SIGNIFICANCE

The C-reactive protein (CRP) is an acute phase protein present in the serum of healthy individuals, which can significantly increase the majority of bacterial and viral infectious processes, tissue damage, inflammation and malignant neoplasms. The elevation of CRP concentration occurs after a few hours of the onset of the inflammatory process, reaching levels of 300 mg/L in 12-24 hours.

## QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use, Biotécnic reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: [sac@biotecnicatda.com.br](mailto:sac@biotecnicatda.com.br)

## AUTOMATION

This product is compatible to the most types of biochemical automatic analysers. The applications are available at: [www.biotecnicatda.ind.br](http://www.biotecnicatda.ind.br)

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