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ENGLISH

INTENDED USE

Kit for the ultra-sensitive determination of C-reactive protein (CRP) in serum. Use in vitro diagnostic.

STORAGE AND HANDLING

- Store at 2 to 8 °C.
- Protect from light.
- Kit expiration date is printed on the label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE

Immunoturbidimetric method in latex particles coated with anti-human CRP antibodies are agglutinated by the CRP present in the sample. The agglutination causes a turbidity whose intensity is proportional to the concentration of CRP.

SAMPLES – TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION

Sample Type: Serum.

Collection, handling and preparation: Perform the sample collection according to Good Clinical Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.

Preservation: CRP in Serum is stable for 7 days at a temperature of 2 to 8 °C or 3 months at a temperature of -20 °C.

PRODUCT DESCRIPTION

R DIL	Tris buffer 20 mmol/L; preservative.	
R LAX	Tris buffer 20 mmol/L; Latex particles coated with goat anti-C-reactive protein; preservative.	
CAL	Human sérum containing CRP; Stabilizer, Preservative. The CRP value is traceable to Reference Material ERM-DA 472/IFCC.	

CONTROL	Human serum containing CRP; buffer Stabilizer; preservative. (INCLUDED IN THE KIT WITH CONTROL)	
----------------	---	--

STABILITY IN USE

- **RDIL, RLAX e CAL:** Once opened, the product (RDIL, RLAX and CAL) in use is stable up to the expiration date printed on the label, provided that it is followed by recommended storage conditions (2 to 8 °C).
- **CONTROL** after reconstitution is stable for 7 days under the recommended storage conditions (2 to 8 °C).
- Reagents should remain outside the specified temperature only for the time required to perform the tests.

TREATMENT AND HANDLING OF THE PRODUCT

A) REAGENT PREPARATION

RDIL, RLAX and CAL: ready for use.

CONTROL: Lyophilized. Leave at a temperature of 10 to 30 °C for 10 minutes. Reconstitute with 1.0 mL of purified water. Gently homogenize by inversion avoiding foaming until complete dissolution of the control.

B) OPERATING INTERVAL

The operating range of the product is from 0.51 mg / L to 5.82 mg / L. For higher values, dilute the sample with 150 mM NaCl (0.9%), carry out a new dosage and multiply the result obtained by the dilution factor.

QUALITY CONTROL

The use of controls should be routine practice in the laboratory. For Calibration and Internal Quality Control Laboratory it is recommended to use the calibrator and control below:

ULTRA-SENSITIVE CRP CAL (kit INCLUDED)	REF	20.017.00
ULTRA-SENSITIVE CRP CONTROL (INCLUDED IN THE KIT WITH CONTROL)		

TEST PROCEDURE, CALCULATION AND INTERPRETATION

REACTION PARAMETERS	
Type of reaction	Fixed time
Wave-length	540 nm
Temperature	37 °C
RDIL*	240 µL
RLAX*	60 µL
Reading 1 (A ₁)	Immediately
NaCl 9 g / L (Blank) / sample / CAL *	3 µL
Reading 2 (A ₂)	4 minutes
Calibration	Calibration curve**

* Volumes can be modified as long as the established ratio is maintained.

**CALIBRATION CURVE:

Dilute the CAL in 9g / L NaCl as follows:

Dilution	1	2	3	4	5	6
CAL (µL)	-	5	10	25	50	100
NaCl 9g / L (µL)	100	95	90	75	50	-
Fator	0	0,05	0,1	0,25	0,5	1

Multiply the concentration of the Ultra-Sensitive CRP CAL, indicated on the label, by the corresponding factor to obtain the Ultra-Sensitive CRP concentration of each dilution.

B) CALCULATIONS

The equipment automatically calculates the analyte concentration in each sample.

C) INTERPRETATION

C-Reactive Protein is an acute phase protein that can increase significantly in most bacterial and viral infectious processes, inflammation and malignant neoplasms. When excluded from inflammatory, infectious or immune causes of elevation of CRP, it can be used in the detection of atherosclerotic process in patients with asymptomatic cardiac diseases.

INTERFERING AND LIMITATIONS

Hemolysis, Jaundice and Lipemia: Hemoglobin > 1000 mg / dL and Bilirubin > 40 mg / dL may interfere with the determination. Triglycerides > 1250 mg / dL interfere in the determination.

Medications: consult recommended literature reference (Young, 2000).

PERFORMANCE CHARACTERISTICS

Sensitivity: Limit of detection: 0,055 mg/L / Limit of quantification: 0,49 mg/L

Prozone High Dose Effect: No prozone effect is observed up to 15,8 mg/L of CRP.

Analytical Specificity: The product specifically determines CRP in the presence of other interfering substances in the sample, up to the concentrations reported above.

Accuracy: The method was compared with a similar method by the determination of 20 samples in duplicate. The regression equation $Y = 1,019X - 0,018$ and correlation coefficient $r = 0,9935$ was obtained. Using this equation the estimated total systematic error of 1.30% for a level of 3.0 mg / L and 1.54% for a level of 5.0 mg / L.

Precision: It was determined using samples at 03 decision levels, two runs in duplicate per day for 20 days, obtaining:

Samples (mg/L)	Repetitions	In-run Precision		Total Precision	
		SD (mg/L)	%CV	SD (mg/L)	%CV
1,54	80	0,010	0,6	0,032	2,1
3,01	80	0,006	0,2	0,067	2,2
6,48	80	0,012	0,2	0,119	1,8

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange reagent bottle caps to avoid cross-contamination.
- Do not use the reagent when it exhibits a visual characteristic in disagreement with that specified in the product's MSDS.
- Avoid leaving reagents out of the specified storage conditions when they are not in use.
- Use specific glass pipettes and disposable tips for each sample, control, standard / calibrator, and reagent.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the requirements by means of periodic checks.
- Proper cleaning and drying of the used material are fundamental factors for the stability of the reagents and obtaining correct results.

REFERENCE RANGES

Interpretation of values, when excluded from inflammatory, infectious or immune causes of ultra-sensitive CRP elevation, according to the Brazilian Directive on Dyslipidemias and Prevention of Atherosclerosis - 2017.:

RISK	mg/L
Low	<1,0
Medium	1,0 a 2,0
High	> 2,0
Very high	> 10

MATERIALS REQUIRED TO CARRY OUT THE TEST

- Spectrophotometer or photometer for reading at 540 nm
- Glass and / or automatic pipettes.
- Test tubes

ALERTS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Disposal, Safety and First Aid information are described in the Individual Product Safety Data Sheet (MSDS) for this product, available at www.biotechnica.ind.br or by telephone (35) -3214-4646.
- Dispose of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS).

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@biotechnicaltda.com.br

AUTOMATION

This product is compatible to the most types of biochemical automatic analyzers. The applications are available at www.biotechnicaltda.ind.br

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- ADAM, A. *et al.*, Analytical Evaluation of a New Latex Agglutination Test for Quantitative Determination of C-Reactive Protein by Laser Nephelometry. **J. Clin. Chem. Clin. Biochem.** v. 23, p. 787-789, 1985.
- HOKAMA, Y.; NAKAMURA, R.M. C- Reative Protein: Current Status and Future Perspectives. **Journal of Clinical Laboratory Analysis.** v. 1, p. 15-27, 1987.
- FALUDI AA, IZAR MCO, SARAIVA JFK, CHACRA APM, BIANCO HT, AFIUNE NETO A *et al.* Atualização da Diretriz Brasileira de Dislipidemias e Prevenção da Aterosclerose - 2017. **Arq Bras Cardiol** 2017; 109(25supl.1):1-76.
- WESTGARD, J. O. *et al.* A multi-rule shewhart chart quality control in clinical chemistry. **Clin. Chem.** v.27, p.493-501, 1981
- YOUNG, D.S. Effects of drugs on clinical laboratory tests - vol. 2, 5 ed. Washington DC: AACC Press, 2000.

TABELA DE SÍMBOLOS INTERNACIONAIS / TABLE OF INTERNATIONAL SYMBOLS / TABLA DE SÍMBOLOS INTERNACIONALES

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