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Ultra Sensitive CRP / PCR Ultrasensível
BT2001700 (KIT)
BT2001700 (KIT COM CONTROLE)

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ENGLISH

BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE

Kit for the ultra-sensitive determination of C-reactive protein (CRP) in serum. *In vitro* diagnostic use.

STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- RDIL, RLAX and CAL:** Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed. Ready-to-use reagents.
- Do not use reagents whose shelf life has expired.
- CONTROL:** Leave at a temperature of 10 to 30 °C for 10 minutes. Reconstitute with 1.0 mL of purified water. Gently homogenize by inversion avoiding foaming. After reconstitution is stable for 7 days if stored at 2 to 8 °C

WORKING PRINCIPLE

Method: Immunoturbidimetric

Latex particles coated with anti-human CRP antibodies are agglutinated by the CRP present in the sample. The agglutination causes a turbidity whose intensity is proportional to the concentration of CRP.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum

Collection and handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Stability Period
Serum	2 to 8 °C	7 days
	-20 °C	3 months

PRODUCT DESCRIPTION

R DIL	Tris buffer 20 mmol/L; preservative.	 
R LAX	Tris buffer 20 mmol/L; Latex particles coated with goat anti-C-reactive protein; preservative.	 
CAL	Human sérum containing CRP; Stabilizer, Preservative. The CRP value is traceable to Referência Material ERM-DA 472/IFCC.	 
CONTROL	Human serum containing CRP; buffer Stabilizer; preservative. (INCLUDED IN THE KIT WITH CONTROL)	 

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratorial quality control, it is recommended the use of the calibrator and control below: Ultra-Sensitive CRP CAL (kit INCLUDED)
Ultra-Sensitive CRP CONTROL (INCLUDED IN THE KIT WITH **REF** BT2001700 CONTROL)

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 540 nm.
- Glass and/or automatic pipettes.
- Test tubes.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

Automation: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotecnica.ind.br

REACTION PARAMETERS	
Type of reaction	Fixed time
Wave length	540 nm
Temperature	37 ° C
RDIL*	240 µL
RLAX*	60 µL
Reading 1 (A1)	Immediately
NaCl 9g/L (Blank) / Sample / CAL*	3 µL
Reading 2 (A2)	4 minutes
Calibration	Calibration curve**

* Sample volumes/CAL, RDIL and RLAX be modified as long as the established ratio is maintained.

**Calibration Curve:

Dilute the CAL in NaCl 9g / L according to follow:

Dilution	1	2	3	4	5	6
CAL (µL)	-	5	10	25	50	100
NaCl 9g/L (µL)	100	95	90	75	50	-
Factor	0	0.05	0.1	0.25	0.5	1

Multiply the concentration of Rheumatoid Factor Turbilatex Turbilatex CAL, indicated on the label, by the corresponding factor to obtain the RF concentration of each dilution.

B) CALCULATIONS

The equipment automatically calculates the analyte concentration in each sample.

C) INTERPRETATION

C-Reactive Protein is an acute phase protein that can increase significantly in most bacterial and viral infectious processes, inflammation and malignant neoplasms. When excluded from inflammatory, infectious or immune causes of elevation of CRP, it can be used in the detection of atherosclerotic process in patients with asymptomatic cardiac diseases.

PERFORMANCE CHARACTERISTICS

Operating range
0.510 to 5.82 mg/L

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
0.055 mg/L	0.490 mg/L

Analytical Specificity		
Hemoglobin	Bilirubin	Triglycerides
1000 mg/dL	40 mg/dL	1250 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	20 in duplicate
Regression Equation	y = 1.019x - 0.018
Correlation Coefficient (R)	0.9935

By applying the regression equation, the total systematic error estimated for decision levels of 3,0 mg/L and 5,0 mg/L was 1.30% and 1.54%, respectively.

Prozone High Dose Effect:

There is no prozone effect up to 15.8 mg/L of CRP.

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CV: Coeficiente de variación; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs of disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.

REFERENCE RANGES

Interpretation of values, when excluded from inflammatory, infectious or immune cause of ultra-sensitive CRP elevation, according to the Brazilian Directive on Dyslipidemias and Prevention of Atherosclerosis - 2017:

RISK	mg/L
Low	< 1,0
Medium	1,0 a 2,0
High	> 2,0
Very high	> 10

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotecnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotecnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotecnica.ind.br.

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