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ENGLISH

INTENDED USE

Kit for internal quality control in Prothrombin Time (PT) and Activated Partial Thromboplastin Time (APTT) trials. Use in *in vitro* diagnostic.

STORAGE AND HANDLING

- Store at 2 to 8 °C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE

Human lyophilized plasma used to verify the performance of determination of Prothrombin Time (PT) and Activated Partial Thromboplastin Time (APTT).

PRODUCT DESCRIPTION

CONTROL 1

Plasma derivative lyophilized with normal PT, APTT, plus buffer and preservatives. (Included in Plasmas Coagulation Control and Plasma Coagulation Control - C1)

CONTROL 2

Plasma derivative lyophilized with normal PT, APTT, plus buffer and preservatives. (Included in Plasmas Coagulation Control and Plasma Coagulation Control – C2)

The Controls present negative results for the HBs antigen for the HIV and HCV antibody. However, they should be treated as

potentially infectious.

STABILITY IN USE

- After reconstitution, the product (CONTROL 1 and CONTROL 2) is stable for up to 6 hours under the recommended storage conditions (2 to 8 °C).
- Controls should remain outside the specified temperature only for the time required to perform the tests.

TREATMENT AND HANDLING OF THE PRODUCT

A) CONTROL PREPARATION

Reconstitute with 1.0 mL of purified water. Homogenize gently by inversion, avoiding foaming. Wait 30 minutes at room temperature until complete dissolution.

QUALITY CONTROL

The use of controls should be routine practice in the laboratory. It is suggested to use one control in the reference range or at the decision level and another control with value in another range of clinical significance. For Internal Quality Control Laboratory it is recommended to use the following::

Plasma Coagulation Control Level 1	REF	36.003.00
Plasma Coagulation Control Level 2		36.003.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

The controls should be used in the same way as the patient samples according to the procedure established for the Prothrombin (TP) and Activated Partial Thromboplastin Time (APTT) tests; and in accordance with the equipment protocol.

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix controls from different lots.
- Do not exchange controls bottle caps to avoid cross-contamination.
- Do not use the controls when it has a visual characteristic that is not in accordance with the product's MSDS requirements.
- Avoid leaving the controls out of the specified storage conditions when they are not in use.
- Use specific disposable tips for each sample and control.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.

- Proper cleaning and drying of the used material are fundamental factors for the stability of the reagents and obtaining correct results.

MATERIALS REQUIRED TO CARRY OUT THE TEST

- Glass and / or automatic pipettes.
- Coagulometer.
- Prothrombin Kit (PT) (CAT BT 36.001.00).
- Activated Partial Thromboplastin Kit (APTT) (CAT BT 36.002.00)
- Purified Water

ALERTS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Disposal, Safety and First Aid information are described in the Individual Product Safety Data Sheet (MSDS) for this product, available at www.biotechnica.ind.br or by telephone (35) -3214-4646.
- Dispose of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS).

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@biotecnicaltda.com.br

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Plasmas Controle de Coagulação	CONTROL 1	1 x 1,0 mL
	CONTROL 2	1 x 1,0 mL

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- FERNÁNDEZ-ROMERO, J.M.; CASTRO, M.D. Luque de. Continuous flow system for the evaluation of the extrinsic coagulation pathway. **Talanta**, v. 43, p.1531-1537, 1996.
- CLSI. One-Stage Prothrombin Time (PT) Test and Activated Partial Thromboplastin Time (APTT) Test; Approved Guideline-Second Edition; Approved Guideline. **CLSI document H47-A2**. Wayne, PA: Clinical and Laboratory Standards Institute; 2008.
- WHO Expert Committee on Biological Standardization. Guidelines for thromboplastins and plasma used to control oral anticoagulant therapy. World Health Organization, 1999: 7-8 (WHO Technical Report Series, N° 889).
- WANG, Huimin et al. Study on hemostatic mechanism of fully soluble hemostatic fiber. **Blood Coagulation And Fibrinolysis**, v. 18, n. 6, p.555-558, 2007.
- VAN HERREWEGEN, Femke et al. The bleeding child. Part II: Disorders of secondary hemostasis and fibrinolysis. **Eur J Pediatr**. V. 171, p. 207-214, 2012.
- Brasil. Ministério da Saúde. Secretaria de Atenção à Saúde. Departamento de Atenção Especializada e Temática. **Manual de Diagnóstico Laboratorial das Coagulopatias Hereditárias e Plaquetopatias**. Brasília: Ministério da Saúde, 2016. 140 p.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. **Clin. Chem**. v.27 p.493-501, 1981.
- YOUNG, D.S. **Effects of drugs on clinical laboratory tests** - vol. 2, 5 ed. Washington DC: AACPS Press, 2000.

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