# Potássio Enzimático

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### ENGLISH

**BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.**

### INTENDED USE

Kit intended to determine potassium in serum samples. Diagnostic use only.

### STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- The Dedicated and Semi Dedicated Reagents presentations have onboard stability of 30 days.
- Reagent ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- Do not use reagents whose shelf life has expired.

### WORKING PRINCIPLE

**Method:** Enzymatic

Potassium is determined enzymatically through the activity of potassium-dependent pyruvate kinase, which reacts with the substrate phosphoenolpyruvate to form pyruvate. The latter reacts with NADH, in the presence of LDH, to form Lactate and NAD. The consumption rate of NADH, which can be spectrophotometrically measured at 340 nm, is proportional to the potassium concentration in the sample.

### SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

**Sample Type:** serum.

**Collection and handling:** collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

**Preservation:**

	Temperature	Stability Period
Serum	4 to 8 °C	6 weeks
	-20 °C	1 year

### PRODUCT DESCRIPTION

**R 1** Tris buffer 250 mmol/L pH 8,2; cryptand 12 mmol/L / L; PEP ≥ 3,3 mmol/L; ADP ≥ 3,15 mmol/L; α-oxoglutarate ≥ 1,2 mmol/L; NADH ≥ 0,35 mmol/L; GLDH ≥ 11 U/mL; PK ≥ 1,2 U/mL.

<b>R 2</b>	LDH ≥ 65 U/mL.	
<b>STD 1</b>	KCl solution with potassium concentration indicated on the label.	
<b>STD 2</b>	KCl solution with potassium concentration indicated on the label.	

Traceability of the standards STD1 and STD2 was obtained with the GBW09152 secondary standard of concentration defined by ID-MS.

### QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the controls below:

Normal Control – Quantinorm	1300300
Pathological Control – Quantialt	1300400

### NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 340 nm.
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

### TEST PROCEDURE, CALCULATION AND INTERPRETATION

#### A) TEST PROCEDURE

1. Pipette in the test tubes:

	STD1	STD2	Sample
<b>R1</b>	720 µL	720 µL	720 µL
<b>Sample</b>	20 µL	20 µL	20 µL

2. Homogenize and incubate at 37 °C for 5 minutes.

3. Add:

<b>R2</b>	240 µL	240 µL	240 µL
-----------	--------	--------	--------

4. Homogenize and transfer immediately to a thermostated cuvette. Set the stopwatch. Note the absorbance at 120 seconds (A1) and 240 seconds (A2) at 340 nm by zeroing with purified water.

### B) CALCULATIONS

Potassium (mmol/L) = (Sample's ΔAbs x Calibration Factor) - Intersection

Calibration Factor =  $\frac{\text{STD2 Concentration} - \text{STD1 Concentration (mmol/L)}}{\Delta\text{Abs STD2} - \Delta\text{Abs STD1}}$

Intersection = Calibration Factor x ΔAbs STD1 - STD1 Concentration

In which:

Sample's or STD's ΔAbs = Sample's or STD's (A1 - A2)

**Automation:** this product is compatible to most types of automatic analyzers. Instrument settings are available at [www.biotechnicalda.ind.br](http://www.biotechnicalda.ind.br)

### C) INTERPRETATION

The function of electrolytes in the human body is multiple. These include maintenance of osmotic pressure and water distribution in the body's various fluid compartments, maintenance of adequate pH, regulation of heart and other muscle function, participation in oxidation-reduction reactions, and catalysis as enzymatic cofactors. Abnormal levels of electrolytes may be the cause or the consequence of various disorders. Potassium is the major cation of intracellular fluid. Its intracellular concentration is maintained by the Na-K-ATPase pump, which continuously carries potassium into the cell. This pump is responsible for the adjustments of the ionic gradients on which the nerve impulses and contractility of the muscle depend. Elevated levels of potassium (hyperkalemia) are present in different clinical conditions such as: metabolic acidosis, dehydration, obstructive nephropathy, use of potassium-sparing diuretics, etc. Low levels of potassium (hypokalemia) occur in alkalosis/metabolic acidosis, extra-renal loss (diarrhea, excessive sweating), poor diet in potassium, among others.

### PERFORMANCE CHARACTERISTICS

Operating range
0,37 to 10 mmol/L

For concentrations above the operating range, dilute the sample with purified water, proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
0,22 mmol/L	0,37 mmol/L

Analytical Specificity		
Hemoglobin	Bilirubin	Triglycerides
200 mg/dL	40 mg/dL	1250 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	40 in duplicate
Regression Equation	$y = 0.997x - 0.001$
Correlation Coefficient (R)	0.9905

By applying the regression equation, the total systematic error estimated for decision levels of 3,80 mmol/L and 8,0 mmol/L was -0.33% and -0.31%, respectively.

#### Precision:

Determined with two runs in duplicate per day for 20 days.

Samples (mmol/L)	Replicates	Within-Run Precision		Total Precision	
		SD (mmol/L)	CV (%)	SD (mmol/L)	CV (%)
3.67	80	0.05	1.50	0.06	1.60
5.71	80	0.08	1.40	0.09	1.50
6.83	80	0.10	1.40	0.10	1.50

CV: Coefficient of variation; SD: Standard deviation

#### RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.

#### REFERENCE RANGES

Serum	3.5 - 5.1 mmol/L
-------	------------------

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

#### Conversion to the International System of Units (mmol/L):

Potassium (mmol/L) x 1.0 = Potassium (mmol/L)

#### WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or calling +55 35 3214 4646

#### QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to [sac@biotechnica.ind.br](mailto:sac@biotechnica.ind.br).
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at [sac@biotechnica.ind.br](mailto:sac@biotechnica.ind.br).

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