Proteína Total VET

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ENGLISH

INTENDED USE

Kit intended to determine total protein in serum samples. Veterinary diagnostic use only.

STORAGE AND HANDLING

- Store at 15 to 30 °C and protect from light.
- Reagent ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (15 to 30 °C) are followed.
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE

Method: Colorimetric

The Cu²⁺ ions, in alkaline medium, react with the peptide bonds of the proteins forming a colored complex that can be spectrophotometrically measured at 550 nm. Color intensity is proportional to the protein concentration in the sample.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum.

Collection and handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Stability Period
Serum	4 to 8 °C -20 °C	3 days 7 days

PRODUCT DESCRIPTION

R 1

Copper sulphate ≥ 5 mmol/L, sodium potassium tartrate ≥ 20 mmol/L, potassium iodide ≥ 10 mmol/L, sodium hydroxide ≥ 0.1 mol/L, detergent.

STD

Phosphate buffer ≥ 20 mmol/L; bovine albumin at a concentration equivalent to 5.0 g/dL; preservative. Traceable to the reference material NIST 927d.

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and controls below:

Calibrador - Autocal Vet	90.039.00
Normal Control – Quantinorm Vet	90.040.00
Pathological Control – Quantialt Vet	90.041.00

REF

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 550 nm (540 – 560 nm).
- Test tubes, glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Pipette in the test tubes:

	Blank	Standard	Sample
STD	-	10 µL	-
Sample	-	-	10 µL
R1	1,0 mL	1,0 mL	1,0 mL

2. Homogenize and incubate at room temperature for 10 minutes.
3. Measure the Sample's and Standard's absorbance against the Blank at 550 nm (540 - 560 nm). The final reaction is stable for 2 hours.

B) CALCULATIONS

Protein (g/dL) = $\frac{\text{Sample Absorbance}}{\text{STD Absorbance}} \times \text{STD Concentration (g/dL)}$

Calculations with the Calibration Factor:

Calibration Factor = $\frac{\text{STD Concentration (g/dL)}}{\text{STD Absorbance}}$

Protein (g/dL) = Sample Absorbance x Calibration Factor

Automation: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotecnica.com.br

C) INTERPRETATION

Proteins are predominantly synthesized in the liver, plasma cells, lymph nodes, spleen and bone marrow. In the course of a disease, the total protein concentration and its fractions may deviate significantly from the normal value. Measurement of the total protein can be used in the diagnosis and treatment of a variety of diseases involving the liver, kidney or bone marrow as well as other metabolic or nutritional disorders. Hypoproteinemia can be caused by diseases and disorders such as: blood loss, nephrotic syndrome, severe burns, salt retention syndrome and Kwarshiorok (acute protein deficiency). Hyperproteinemia can be observed in cases of severe dehydration and in diseases such as multiple myeloma. Alterations in the relative percentage of plasma proteins can occur without change in the amount of total protein, being caused by variation in the percentage of one of its fractions. The A/G ratio is commonly used as an index of distribution of the albumin and globulin fractions. Changes in this ratio can be observed in hepatic diseases, glomerulonephritis, nephrotic syndrome, acute hepatitis, lupus erythematosus as well as in certain acute and chronic inflammations.

PERFORMANCE CHARACTERISTICS

Operating range
0.23 to 10.46 g/dL

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity
Quantification Limit
0.22 g/dL

Analytical Specificity			
Species	Hemoglobin	Bilirubin	Triglycerides
Canines	500 mg/dL	40 mg/dL	500 mg/dL
Felines	300 mg/dL	40 mg/dL	500 mg/dL
Equines	200 mg/dL	40 mg/dL	500 mg/dL
Bovines	200 mg/dL	40 mg/dL	250 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy - Canines	
Number of Samples	40 in duplicate
Regression Equation	y = 1.0056x - 0.0273
Correlation Coefficient (R)	0.9991
Accuracy - Felines	
Number of Samples	20 in duplicate
Regression Equation	y = 0.9747x + 0.2443
Correlation Coefficient (R)	0.9830
Accuracy - Equines	
Number of Samples	11 in duplicate
Regression Equation	y = 0.9564x + 0.3505
Correlation Coefficient (R)	0.9968
Accuracy - Bovines	

Number of Samples	10 in duplicate
Regression Equation	$y = 0.9639x + 0.2468$
Correlation Coefficient (R)	0.9954

Precision:

Determined with two runs in 20 replicates.

Samples (U/L)	Number of samples	Within-Run Precision		Samples (U/L)	Number of Samples	Total Precision	
		SD (U/L)	CV (%)			SD (U/L)	CV (%)
3.88	20	0.07	1.90	3.89	40	0.07	1.80
11.54	20	0.20	1.70	11.54	40	0.24	2.10

CV: Coefficient of variation; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.

REFERENCE RANGES

Canine	5.4 – 7.1 g/dL
Feline	5.4 – 7.8 g/dL
Bovine	6.74 – 7.46 g/dL
Equine	5.2 – 7.9 g/dL

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the International System of Units (g/L):

Total Protein (g/dL) X 10 = Total Protein (g/L)

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotechnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotecnicaitda.com.br.

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- YOUNG, D.S. Effects of drugs on clinical laboratory tests - vol. 2, 5 ed. Washington DC: AACC Press, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. Clin. Chem. v.27 p.493-501, 1981.
- BRASIL. Farmacopéia Brasileira. 5ª ed. Brasília: Agência Nacional de Vigilância Sanitária, 2010.
- KANEKO, Jiro J; HARVEY, John W; BRUSS, Michael L. Clinical biochemistry of domestic animals. 5th. ed. San Diego: Academic Press, c1997 932p.
- THRALL, Mary Anna. Hematologia e bioquímica clínica veterinária. São Paulo: Roca, 2007. 582 p.
- KERR, Morag G. Exames laboratoriais em medicina veterinária: bioquímica clínica e hematologia. 2. ed. São Paulo: Roca, 2003. 436 p.
- DUNCAN, J. Robert; PRASSE, Keith W. Patologia clínica veterinária. Rio de Janeiro: Guanabara Koogan, 1982. 217p.
- BUSH, B. M. Manual del laboratorio veterinario de analisis clinicos. Zaragoza: Acribia, 1982. 467p.
- Stockham, Steven L, Scott, Michael A. Fundamentals of veterinary clinical pathology. Blackwell Publishing, 2008.
- WILLARD, Michael D; TVEDTEN, Harold; TURNWALD, Grant H. Small animal clinical diagnosis by laboratory methods. 2nd ed. Philadelphia: W.B. Saunders, 1994. 377p.
- COLES, Embert H; GOMES E NASCIMENTO, Sonia Cardoso de Aguiar; NASCIMENTO, Fernando Gomes do. Patologia clínica veterinária. 3. ed. São Paulo: Manole, 1984. 566p.

TABELA DE SÍMBOLOS INTERNACIONAIS / TABLE OF INTERNATIONAL SYMBOLS / TABLA DE SÍMBOLOS INTERNACIONALES

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