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ENGLISH
BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE
Kit intended for the quantitative determination of protein in urine and cerebrospinal fluid (CSF) samples. *In vitro* diagnostic use.

STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagent ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8°C) are followed.
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE
Method: Pyrogallol Red
The protein reacts with pyrogallol red and sodium molybdate, in acidic medium, forming a colored complex that can be spectrophotometrically measured at 600 nm. The intensity of the color is proportional to the concentration of protein in the sample.

SAMPLE: TYPE, COLLECTION, HANDLING, PREPARATION AND STABILITY
Sample Type: urine and cerebrospinal fluid (CSF).
Collection and Handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.
Preservation:

Urine	Temperature	Stability Period
	4 to 8 °C	7 days
Cerebrospinal Fluid	-20 °C	1 month
	4 a 8 °C	6 days
	-20 °C	1 year

Preparation:

Urine: homogenize the sample and separate 20 mL aliquot. Centrifuge for 10 minutes at 3000 rpm and utilize the supernatant to proceed with the assay.
Cerebrospinal fluid: centrifuge for 10 minutes at 3000 rpm and use the supernatant in the assay.

R

1

STD

CONTROL

Succinate buffer > 10 mmol / L; sodium molybdate > 10 µmol / L; pyrogallol red > 10 µmol / L; preservative; detergent; stabilizers. Pipes buffer; preservatives; stabilizer; Bovine albumin at a concentration corresponding to 50,0 mg / dL protein. Traceable to NIST reference material 927d.

Phosphate buffer; stabilizers; detergent; preservatives and bovine albumin. (INCLUDED IN THE KIT WITH CONTROL)

QUALITY CONTROL
The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the control below:

SENSIBLE PROTEIN CONTROL (INCLUDED IN THE KIT WITH CONTROL)	REF	BT1001600
---	------------	-----------

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 600 nm (580 – 620 nm).
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE
1. Pipette in the test tubes:

	Blank	Sample	STD
STD	---	---	50 µL
Sample	---	50 µL	---
R1	1,0 mL	1,0 mL	1,0 mL

- Homogenize and incubate for 10 minutes at 37 °C. Remove it from the water bath and let it rest in temperature of 15 to 30 °C before reading.
- Measure the Standard's and Sample's absorbance against the blank at 600 nm (580 – 620 nm). The final reaction is stable for 15 minutes.

B) CALCULATIONS
Protein (mg/dL) = $\frac{\text{Sample Absorbance} - \text{Blank x STD Concentration (mg/dL)}}{\text{STD Absorbance} - \text{Blank}}$

Calculations with the Calibration Factor:
Calibration Factor = $\frac{\text{STD Concentration (mg/dL)}}{\text{STD Absorbance} - \text{Blank}}$
Protein (mg/dL) = (Sample Absorbance - Blank) x Calibration Factor

Calculations for urine samples:
Protein (mg/24 hours) = $\frac{\text{Protein (mg/dL)} \times \text{Urinary Volume (mL)}}{100}$

Automation: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotechnicaltda.ind.br

C) INTERPRETATION
Urine:
Proteinuria is defined as the presence of excess protein in the urine, which occurs in many types of kidney disease, after intense exercise or dehydration. The causes may be due to increased glomerular permeability, altered tubular reabsorption, elevated low-molecular-weight plasma proteins, abnormal secretion of proteins in the urinary tract and post-renal bleeding.
Cerebrospinal fluid:
Examination of CSF total proteins is mainly used to detect increased permeability of the blood-brain barrier or increased intrathecal production of immunoglobulins. The permeability of the blood-brain barrier to plasma proteins is increased by intracranial pressure due to brain tumor, intracerebral hemorrhage, traumatic injury, inflammation due to bacterial or viral meningitis, encephalitis or poliomyelitis.

Sensitivity		
Sample type	Limit of Detection	Limite of Quantification
Urine	0.50 mg/dL	4.80 mg/dL
Cerebrospinal fluid	-	14.39 mg/dL

Operating range	
Urina	4.80 a 97.00 mg/dL
Cerebrospinal fluid	14.28 a 68.55 mg/dL

BIOTÉCNICA IND.COM. LTDA - CNPJ: 025340690001/20, Avenida Washington Ribeiro, 200 - Industrial Miguel de Luca
CEP 37072-030 Varginha-MG BRASIL Tel/fax: +55 35 3214 4646 www.biotechnica.ind.br

Revisão/Revision/Revisión: 04 – 04/04/2024

For concentrations above the operating range, dilute the sample with purified water, proceed with a new dosage and multiply the result by the dilution factor.

Analytical Specificity		
Sample type	Hemoglobin	Bilirubin
Urine	Any concentration	> 1.5 mg/dL
cerebrospinal fluid	Any concentration	NA

Interfering substances up to the values presented above do not cause significant alterations in results. Any concentration of membrane stabilizer (Transfix®) interferes with the determination. For drugs, consult the recommended reference (Young, 2000).

Sample type		Accuracy	
Urine	Number of samples	40 em duplicata	
	Regression Equation	y = 0,992x + 0,028	
	Correlation Coefficient (R²)	0,9975	
Cerebrospinal fluid	Number of samples	20 em duplicata	
	Regression Equation	y = 0,979x + 1,0645	
	Correlation Coefficient (R*)	0,9964	

Urine: by applying the regression equation, the total systematic error estimated for decision levels of 30.0 mg/dL and 70.0 mg/dL was -0.71% and -0.76%, respectively.

Precision:

Sample Type	mg/dL	Repetition	Precision Intra-Assay		Precision Inter-Assay	
			SD (mg/dL)	CV (%)	SD (mg/dL)	CV (%)
Urin	15,20	80	0.276	1.80	0.752	5.00
	47,50	80	0.518	1.10	1.905	4.00
	87,10	80	0.693	0.80	1.730	2.00
Cerebrospi nal Fluid	18,80	40	0.352	1.90	0.384	2.00
	51,16	40	0.557	1.10	0.518	1.00
	81,16	40	0.900	1.10	0.838	1.00

CV: Coefficient of variation; SD: Standard deviation.

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

Urine		Cerebrospinal Fluid
Adults	1 - 14 mg/dL	15 - 40 mg/dL
	< 100 mg/dia	
Pregnant	< 150 mg/dia	

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to International System of Units (g/L):

Protein (mg/dL) x 0,01 = Protein (g/L)

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotechnica.ind.br or calling for +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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- BIRMINGHAM, D. J. et al. Relationship between Albuminuria and Total Proteinuria in Systemic Lupus Erythematosus Nephritis: Diagnostic and Therapeutic Implications. Clinical Journal Of The American Society Of Nephrology, [s.l.], v. 3, n. 4, p.1028-1033, 9 abr. 2008
- BURTS, C. A.; ASHWOOD, E. R. **Tietz textbook of clinical chemistry**, 4 ed. Philadelphia: W.B.Saunders, 1994. 836p.
- GNUTZMANN, L. V.; PLEWKHA, J.; SULDOFSKI, M. T.; FELISBERTO, M.; NESI, V. Análise dos valores de referência do líquido cefalorraquidiano. RBAC, v. 48, n. 3, p. 189-197, 2016.
- MARSHALL, Thomas; WILLIAMS, Katherine M.. Total Protein Determination in Urine: Elimination of a Differential Response between the Coomassie Blue and Pyrogallol Red Protein Dye-binding Assays. Clinical Chemistry, v. 46, n. 3, p.392-398, 2000.
- SOLDIN, J. O.; WONG, E. C.; BRUGNARA, C. Pediatric Reference Intervals, 7ª.edition Washington: AACCPress, 2011.p.193.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. **Clin. Chem.** v.27 p.493-501, 1981.
- YALAMATI, Padma; KARRA, Madhu Latha; BHONGIR, Aparna V.. Comparison of Urinary Total Proteins by Four Different Methods. Indian Journal Of Clinical Biochemistry, [s.l.], v. 31, n. 4, p.463-467, 8 fev. 2016. Springer Nature.
- YOUNG, D.S. **Effects of drugs on clinical laboratory tests - vol. 2**, 5 ed. Washington DC: AACCPress, 2000.
- HRISHI, Ajay Prasad; SETHURAMAN, Manikandan. Cerebrospinal fluid (CSF) analysis and interpretation in neurocritical care for acute neurological conditions. Indian Journal of Critical Care Medicine: peer-reviewed, official publication of Indian Society of Critical Care Medicine, v. 23, n. Suppl 2, p. S115, 2019.

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