TROMBOPLASTINA PARCIAL

ATIVADA (TTPa)

Activated Partial Thromboplastin Time (APTT) | Tromboplastina Parcial Activada (TTPa)

Ref. 36.002.00

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ENGLISH

INTENDED USE

Kit for the determination of time of Partial Activated Thromboplastin Time (APTT) in plasma. In vitro diagnostic use.

STORAGE AND HANDLING

- Store at 2 to 8 °C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE

The activator (ellagic acid) together with the procoagulant phospholipid (cephalin) activate the coagulation factors of the citrated plasma and with the addition of calcium ions the coagulation process begins. The time of formation of the fibrin clot is the time of Activated Partial Thromboplastin

SAMPLE – PREPARATION AND STABILITY

Sample Type: Plasma (citrate).

Collection, handling and preparation: Perform the collection of the sample according to the Good Practices of the Clinical Laboratory. The samples to be analyzed should be treated as potentially infectious material.

- Use citrated plasma in the proportion of 9 parts of venous blood and 1 part of 3.2% sodium citrate.
- Avoid prolonged garroting.
- Do not use the sample obtained by traumatic puncture.
- If clots are present, discard the sample and perform a new collection.
- Use polypropylene or silicized glass tubes.

Preservation: The plasma sample is stable for 4 hours if stored at 15 -30 °C and / or 4 to 8 °C, 2 weeks at -20 °C and 6 months at -70 °C. Preferably samples should be thawed in a water bath at 37°C for approximately 5 minutes.

PRODUCT DESCRIPTION

R 1 MES buffer ≥ 0.2 mmol / L; rabbit brain extract solution; ellagic acid <0.05 g/L; stabilizer; detergent; preservative.

R 2 Calcium chloride solution ≥ 10 mmol / L; preservative.

STABILITY IN USE

- After being opened, the product (R1 and R2) in use is stable to the expiration date printed on the label since recommended storage conditions are followed (2 to 8 °C)..
- Reagents should remain outside the specified temperature only for the time required to perform the tests

TREATMENT AND HANDLING OF THE PRODUCT

A) REAGENT PREPARATION

R1 and R2: Ready-to-use reagents

BEFORE THE TEST, GENTLY HOMOGENIZE THE R1 WITH SUCESSIVE INVERSIONS

QUALITY CONTROL

The use of controls should be routine practice in the laboratory. It is suggested to use one control in the reference range or at the decision level and another control with value in another range of clinical significance. For the Internal Control of Laboratory Quality it is recommended to use the following plasmas:

Plasma Coagulation Control **REF** 36.003.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

It is recommended to carry out the analysis in duplicate for each sample.

1. Preheat the **R2** to 37 °C.
2. Incubate 100 µL of the sample at 37 °C for 2 minutes.
3. Add 100 µL of **R1**, homogenize and incubate the reaction tube at 37 °C for 3 minutes.
4. Add 100 µL of the preheated **R2**, turning the timer on simultaneously.
5. The reaction should occur under constant agitation.
6. Record, in seconds, the reaction time until clot formation.

B) CALCULATIONS

Calculate the average of the times obtained in the duplicates for each sample.

C) INTERPRETATION

APTT is a screening test for the evaluation of intrinsic and common coagulation pathway factors, the presence of lupus anticoagulant and to monitor levels of unfractionated heparin in plasma. The prolongation of APTT may occur in inherited coagulopathies (such as hemophilia) or acquired (presence of autoantibodies inhibiting blood clotting, such as lupus anticoagulant). Other conditions include: use of anticoagulants such as heparin, disseminated intravascular coagulation, and liver diseases.

INTERFERING AND LIMITATIONS

Hemolysis, Jaundice and Lipemia: Hemoglobin> 200 mg dL / Triglycerides> 2000 mg/dL / Bilirubin> 40 mg dL interferes with the dosage.

Medications: see recommended literature reference (Young, 2000).

PERFORMANCE CHARACTERISTICS

Sensitivity: The product was evaluated using deficient plasmas in coagulation factors, obtaining for a percentage of 30% of activity the following results:

	Fator VIII	Fator IX	Fator XI	Fator XII
TTPa (seconds)	27,4	28,2	39,3	28,8

Analytical Specificity: The product specifically determines the activated partial thromboplastin time in the presence of other interfering substances in the sample, up to the concentrations reported above.

Accuracy: The method was compared to a similar method by the determination of 40 samples in duplicate. The regression equation $y = 0,999x - 0,03$ and

correlation coefficient $r = 0.9904$ was obtained. Using this equation the total systematic error estimated was -0.22% for a level of 27.8 seconds and -0.14% for a level of 64.7 seconds.

Precision: It was determined using two samples of control plasma, two daily runs in duplicate for 10 days, obtaining:

Samples (seconds)	Repetitions	Intra-run Precision		Total Precision	
		SD	%CV	SD	%CV
24,54	40	0,76	3,1	1,42	5,8
60,04	40	0,59	1,0	1,35	2,2

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with Good Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange the caps from different reagents in order to avoid cross contamination.
- Do not use the reagent if it displays any sign in disagreement with the one specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions, when they are not in use.
- Use glass pipettes and disposable tips for each sample, control and reagent.
- The level of the water bath must be greater than the test tubes containing the reaction.
- The laboratory should establish chemical, microbiological and particulates requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
- The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.

REFERENCE RANGES

TROMBOPLASTINE TIME PARTIAL ACTIVATED (seconds)
26 to 38

These values are for guidance only. For coagulation tests, the laboratory should establish the reference interval considering the representativeness of the local population, selecting blood donors after clinical and laboratory screening

ADDITIONAL EQUIPMENT

- Coagulometer / Water bath (37 °C).
- Pipettes.
- Clock or Stopwatch.
- Test tubes.

WARNINGS AND PRECAUTIONS

- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotechnica.ind.br or calling for 55-35-3214-4646.
- Dispose of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS)

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@biotechnicaltda.com.br

AUTOMATION

This product is automated. The applications are available at www.biotechnicaltda.ind.br

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