## Triglicérides

Triglycerides | Triglicéridos  
Ref. 10.010.00

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**ENGLISH**

**INTENDED USE**



Kit intended to determination of triglycerides in serum and plasma. Diagnostic use only.

**STORAGE AND HANDLING**

- Store at 2 to 8 °C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

**WORKING PRINCIPLE**

Triglycerides are hydrolyzed by lipoprotein lipase producing free glycerol, which is transformed into glycerol-3-phosphate by the action of glicerolquinase. The glycerol 3-phosphate is oxidized to dihydroxyacetone and hydrogen peroxide by glycerol-P oxidase. Hydrogen peroxide reacts with 4-aminoantipyrine and 4-chlorophenol under the catalytic action of peroxidase, producing a pink compound (quinonimine) having an absorption maximum at 505 nm.

**SAMPLE – PREPARATION AND STABILITY**

**Sample Type:** Serum and plasma (EDTA).  
**Collection, handling and preparation:** Make the collection of the sample in accordance with Good Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.

**Preservation:** The triglycerides is stable in serum and plasma (EDTA) for 3 days from 4 to 8 °C and 30 days from - 20 °C.

**PRODUCT DESCRIPTION**

<b>R</b>	<b>1</b>	Pipes buffer ≥ 20 mmol/L; 4-chlorophenol ≥ 1 mmol/L; ; 4-aminoantipyrine ≥ 0,1 mmol/L; ATP - adenosine triphosphate ≥ 0,5 mmol/L; Glycerol kinase ≥ 500 U/L; Peroxidase ≥ 1000 U/L; Lipoprotein lipase ≥ 1000 U/L; Glycerol-3-phosphate oxidase ≥ 1000 U/L; activators; detergents, stabilizers and preservative.	
<b>STD</b>		Glycerol is equivalent to the triglyceride concentration of 200 mg / dL, preservatives.	
Traceable to NIST reference material 1951b.			

**STABILITY IN USE**

- The stability of the product (R1 e STD) in use is 15 months as long as followed by the recommended storage conditions (2 to 8 °C).
- Reagents must remain outside the specified temperature only the time required for testing

**TECHNICAL PROCEDURE**

**A) REAGENT PREPARATION**  
**Reagent 1 (R1)**  
Ready-to-use Reagent

**B) OPERATING RANGE**  
The product operating range is from 9,00 mg/dL to 800 mg/dL. For higher values, dilute the sample with NaCl 150 mM (0,9%), repeat the assay and multiply the result by the dilution factor.

**QUALITY CONTROL**

Use of controls should be a routine practice in the lab. It is suggested to use a control within the reference range or decision level and another control with value in another range of clinical significance. For the internal quality control of the lab it's indicated the use of the calibrator serum and control sera below:  
Calibrator Serum - Autocal H 13.002.00  
Normal Control Serum – Quantinorm 13.003.00  
Pathological Control Serum – Quantialt 13.004.00

**TEST PROCEDURE, CALCULATION AND INTERPRETATION**

**A) TEST PROCEDURE**

1. Pipette in the assay tubes:

	Blank	Standard	Sample
STD	-	10 µL	-
Sample	-	-	10 µL
R 1	1,0 mL	1,0 mL	1,0 mL

2. Homogenize and incubate 10 minutes at 37 °C.  
3. Measure the absorbance of Standard and Sample against the blank at 505 nm (490-510 nm). The final reaction is stable 25 minutes.

**B) CALCULATIONS**

**Calculations**  
Triglycerides (mg/dL) =  $\frac{\text{Sample Absorbance}}{\text{STD Absorbance}} \times \text{STD concentration (mg/dL)}$

**Example:**  
Sample Absorbance = 0,277  
STD Absorbance = 0,248  
STD concentration = 200 mg/dL  
Triglycerides (mg/dL) =  $\frac{0,277}{0,248} \times 200 = 223 \text{ mg/dL}$