	
URÉIA UV Urea UV Urea UV Ref. 10.012.00	Responsável Técnico Dr. Glison Sérgio Rizzo CRF MG - 5310 MS 80027310263

FINALIDADE

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ENGLISH

INTENDED USE


Kit intended to determine urea in serum, plasma and urine. Diagnostic use only.

STORAGE AND HANDLING

- Store at 2 to 8 °C
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE

The urea from sample is hydrolyzed by the enzyme urease producing carbon dioxide and ammonium ions. The ammonium ions, in presence of the enzyme glutamate dehydrogenase and the substrates NADH and α-ketoglutarate, produces NAD⁺ and glutamate. The velocity of decrease of NADH in the assay may be measured in spectrophotometer at 340 nm being proportional to the urea concentration.

SAMPLE – PREPARATION AND STABILITY


Sample Type: Serum, plasma and urine.

Collection, handling and preparation: Make the collection of the sample in accordance with Good Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.

Preservation: The urea is stable in serum and plasma for 7 days from 2 to 8 °C and 12 months from -20 °C.

The urea is stable in urine for 2 days from 20 to 25 °C, 7 days from 4 to 8 °C and 28 days from -20 °C. Work with samples collected in the period of 24 hours with 2mL HD 50% Centrifuge the sample for 10 minutes at 3,000 rpm. Dilute a small part of the centrifuged urine in a 1:50 proportion with distilled or deionized water. When necessary, the dilution of urine should be changed to obtain a result within the operating range of the method. Use the diluted urine to proceed with the assay. Multiply the result by 50 (or the fact that has been applied).

PRODUCT DESCRIPTION

R 1	TRIS buffer > 100 mmol/L; EDTA > 1 mmol/L; α-ketoglutarate > 2 mmol/L; ADPK > 1 mmol/L; Urease > 5000 U/L; glutamate dehydrogenase > 1000 U/L; activators; detergent; stabilizers; preservative.		
R 2	Carbonate buffer > 2 mmol/L; NADH > 0,5 mmol/L; preservative.		
STD	Phosphate buffer > 2 mmol/L; preservative; Urea in a concentration equivalent to 70 mg/dL.		

The determination of Urea is traceable to reference material NIST 912a.

STABILITY IN USE

- The stability of the product (R1, R2 e STD) in use is 15 months as long as followed by the recommended storage conditions (2 to 8 °C).
- The stability of working reagent is 4 weeks, as long as followed by the conditions of preparation and recommended storage (2 to 8 °C).
- Reagents must remain outside of the specified temperature only the time required for testing.

TECHNICAL PREPARATION

A) REAGENT PREPARATION
Work Reagent (WR)
 Mix in proportion 4 parts of R1 + 1 part of R2. Homogenize it gently. The Work Reagent is stable 4 weeks at 2 to 8 °C.

OBS: The absorbance of the reagent below 0,900 at 340 nm with the instrument adjusted to zero with purified water, indicates deterioration.

B) OPERATING RANGE

The product operating range is from 13,3 mg/dL to 200 mg/dL. For higher values, dilute the sample with NaCl 150 mM (0,9%), repeat the assay and multiply the result by the dilution factor.

QUALITY CONTROL

Use of controls should be a routine practice in the lab. It is suggested to use a control within the reference range or decision level and another control with value in another range of clinical significance. For the internal quality control of the lab it is indicated the use of the calibrator serum and control sera below:

Calibrator Serum - Autocal H	13.002.00
Normal Control Serum - Quantinorm	13.003.00
Pathological Control Serum - Quantit	13.004.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Heat the reagent for 3 minutes at 37°C
2. Zero the read ng equipment with purified water at 340 nm
3. Repeat in the assay tube:

Volumen	
Wbrk Reagent	1.0 mL
Sample	10 µL

3. Homogenize and insert immediately in the thermostatic cuvette at 37°C
4. After 30 sec. note the initial absorbance (A1) and read again after exactly 120 sec. (A2) of sample and standard.

B) CALCULATIONS

Calculation of Serum Urea Nitrogen

$$\text{Urea (mg/dL)} = \frac{(A1 - A2 \text{ sample})}{(A1 - A2 \text{ standard})} \times \text{STD Concentration (mg/dL)}$$

Calibration factor (CF)

$$\text{CF} = \frac{\text{STD Concentration (mg/dL)}}{(A1 - A2 \text{ standard})}$$

$$\text{Urea (mg/dL)} = (A1 - A2 \text{ sample}) \times \text{CF}$$

Calculations of Urea

$$\text{Urea (mg/24 hours)} = \frac{\text{Urea (mg/dL)} \times \text{urinary volume (mL)}}{100}$$

*Value corrected by dilution factor

Conversion to g/24 hours unit

$$\text{Urea (g/24 hours)} = \frac{\text{Urea (mg/24 hours)}}{1000}$$

C) INTERPRETATION

Urea is the major nitrogenous metabolite of protein catabolism. More than 90% of urea is excreted through the kidneys, therefore renal disease is associated with accumulation of urea in the blood. The dosage of blood and urinary urea can provide useful clinical information. It can be elevated due to various factors, such as high protein diet, high protein catabolism, reabsorption of blood proteins after gastritis, renal bleedings, treatment with corticosteroids and androgens, dehydration and reduced perfusion of the kidneys. It is also elevated in obstructive post-renal conditions such as malignant tumors, nephritis and prostatitis.

INTERFERENCES AND LIMITATIONS

Anti-coagulants: Citrate, fluoride and oxalate may interfere in the reaction.
Hemolysis, Jaundice and Lipemic Serum: Hemoglobin > xx mg/dL, Bilirubin > xx mg/dL, Triglycerides > xx mg/dL.

Drugs: see recommended reference (Young, 2000).

PERFORMANCE CHARACTERISTICS

Sensitivity: Detection limit: 3.52 mg/dL / quantification limit: 13.30 mg/dL.

Analytical Specificity: The product determines Urea specifically in the presence of other interfering substances in the sample up to the concentrations stated above.

Accuracy: The method was compared with a standard method for the determination of 40 samples in duplicate. The regression equation obtained was $y = 1.016x - 0.572$ and the correlation coefficient $r = 0.9971$. Using this equation the total systematic error estimated is 0.46% to a level of 50 µg/L and 1.16% to a level of 130 µg/L.

Precision: Determined using samples in 03 levels of decision, with two runs in duplicate per day during 20 days, being obtained:

Sample (µg/L)	Repetitions	Within-Run Precision		Total Precision	
		SD (µg/L)	%CV	SD (µg/L)	%CV
31,726	80	0.346	1.1	0.474	1.5
108,31	80	0.801	0.7	1,654	1,5
198,705	80	1,345	0,7	2,764	1,4

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with Good Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange the caps from different reagents in order to avoid cross contamination.
- Do not use the reagent if it displays any significant disagreement with the one specified in the product MSDS.
- Avoid leaking reagents outside of the specified storage conditions, when they are not in use.
- Use glass pipettes and disposable tips for each sample control, standard / calibrator and reagent.
- The level of the water bath must be greater than the test tubes containing the reaction.
- The level of water in the thermostatic bath must be over the level of mix reaction within the assay tubes.
- The laboratory should establish chemical, microbiological and particulates requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
- The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.

REFERENCE RANGES

	mg/dL	mmol/L
Serum and Plasma	13 - 45	2.17 - 7.51
Urine	26 - 43 g/24 hours	

These values are intended for orientation only. It's recommended that each lab establishes its own reference range.

Conversion for International System of Units (SI): mmol/L

$$\text{Urea (mg/dL)} \times 0.167 = \text{mmol/L}$$

ADDITIONAL EQUIPMENT

- Spectrophotometer or photometer with a thermostatic cuvette able to read at 340 nm
- Thermostatic water bath at 37°C
- Glass Pipettes and/or automatic
- Clock or chronometer.
- Assay tubes

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to Good Laboratory Practice in a proper place for potentially infective material.
- The information for disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.bioteccalindbr.com or calling for 55-35-3214-4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use, Biotecca reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sc@bioteccaltda.com.br

AUTOMATION

This product is compatible to the most types of biochemical automatic analyzers. The applications are available at www.bioteccaltda.com.br

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