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ENGLISH
INTENDED USE
Kit for the determination of zinc in serum and plasma. In vitro diagnostic use.
STORAGE AND HANDLING

- Store at 2 to 8 °C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE
Zinc ions (Zn++) form with nitro-PAPS a colored complex that can be measured at 570-582 nm. The intensity of the formed color is proportional to the amount of zinc present in the sample.

SAMPLES – TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION
Sample Type: Serum and plasma (heparin).
Collection, handling and preparation: Perform the sample collection according to Good Clinical Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.
Preservation: Zinc in serum and plasma is stable for 8 days if stored at 2 to 8 °C and 1 year if stored at -20 °C.

PRODUCT DESCRIPTION

R 1	Borate buffer pH 8,2; Salicylaldoxime; Dimethylglyoxime; Surfactant and Preservative.	X
R 2	Nitro-PAPS, Preservative.	X

STD	Zinc Nitrate in solution. Traceable to reference material NIST-SRM 3168a.	X
------------	---	----------

STABILITY IN USE

- After opening, the product (R1, R2, STD) in use is stable up to the expiration date printed on the label, provided that it is followed by recommended storage conditions (2 to 8 °C).
- Reagents should remain outside the specified temperature only for the time required to perform the tests.

TREATMENT AND HANDLING OF THE PRODUCT

A) REAGENT PREPARATION
R1, R2 and STD: Ready-to-use reagents.

B) OPERATING INTERVAL
The operating range of the product is 36.3 µg / dL at 718 µg / dL. For higher values, dilute the sample with 150 mM NaCl (0.9%), carry out a new dosage and multiply the result obtained by the dilution factor.

QUALITY CONTROL
The use of controls should be routine practice in the laboratory. For Internal Quality Control Laboratory it is recommended to use the control serum below:

Control Serum - Quantialt	REF	13.004.00
---------------------------	------------	-----------

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE
1. Before starting, allow the reagents to reach a temperature of 15 to 30 °C.
2. Pipette into test tubes:

	Blank	Standard	Sample
STD	---	50 µL	---
Sample	---	---	50 µL
Purified water	50 µL	---	---
R1	800 µL	800 µL	800 µL

2. Homogenize gently and incubate the tubes for 5 minutes at 37 °C.

3. Add:

	Blank	Standard	Sample
R2	200 µL	200 µL	200 µL

4. Homogenize gently and incubate the tubes for 5 minutes at 37 °C.
5. Measure the absorbance of the Standard and the Sample against the Blank at 578 nm (570-582). The final color of the reaction is stable for 15 minutes.

B) CALCULATIONS
Zinc (µg/dL) = $\frac{\text{Sample Absorbance}}{\text{Absorbance of STD}} \times \text{Concentration of STD (µg/dL)}$
Example:
Concentration of STD = 200 µg/dL
Sample Absorbance = 0,491
Absorbance of STD = 0,263
Zinc (µg/dL) = $\frac{0,491}{0,263} \times 200 = 373,4 \mu\text{g/dL}$

With Calibration Factor:
Calibration Factor = $\frac{\text{Concentration of STD (µg/dL)}}{\text{Absorbance of STD}}$
Zinc (µg/dL) = Sample Absorbance x Calibration Factor
Example:
Calibration Factor = $\frac{200}{0,263} = 760,5$
Zinc (µg/dL) = 0,491 x 760,5 = 373,4 µg/dL.

C) INTERPRETATION
Zinc is a cofactor in more than 70 enzymatic reactions, involving different metabolic processes such as the synthesis or degradation of carbohydrates, lipids, proteins and nucleic acids. Zinc deficiency can cause anemia, hepatosplenomegaly, slow growth and healing of ulcers and wounds, as well as changes in smell and taste. Low levels of zinc can occur in several physiological situations: pregnancy, contraceptive use and in some pathologies such as myocardial infarction, cirrhosis, lung infection, lymphoma and carcinoma.

INTERFERING AND LIMITATIONS
Hemolysis, Jaundice and Lipemia: Hemoglobin > 25 mg / dL, Bilirubin > 20 mg / dL and Lipemia interfere with the dosage.
Medications: consult recommended literature reference (Young, 2000).

PERFORMANCE CHARACTERISTICS

Sensitivity: Limit of detection: 9.48 µg / dL / Limit of quantification: 36.3 µg / dL
Analytical Specificity: The product specifically determines Zinc in the presence of other interfering substances in the sample up to the concentrations reported above.

Accuracy: The method was compared with a similar method by the determination of 30 samples in duplicate. The regression equation $Y = 1.006 X - 0.87$ and correlation coefficient $r = 0.9993$ were obtained. Using this equation the total systematic error estimated was - 0.27% for a level of 100 µg / dL and 0.45% for a level of 603 µg / dL.

Precision: It was determined using samples at 03 decision levels, two runs in duplicate per day for 20 days, being obtained.

Samples (µg/dL)	Repetitions	In-run Precision		Total Precision	
		SD (µg/dL)	%CV	SD (µg/dL)	%CV
64,5	80	1,06	1,6	1,38	2,1
135	80	1,23	0,9	1,43	1,1
245	80	1,04	0,4	1,27	0,5

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange reagent bottle caps to avoid cross-contamination.
- Do not use the reagent when it has a visual characteristic that is not in accordance with the product's MSDS requirements.
- Avoid leaving the reagents out of the specified storage conditions when they are not in use.
- Use specific glass pipettes and disposable tips for each sample, control, standard / calibrator and reagent.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the used material are fundamental factors for the stability of the reagents and obtaining correct results.

REFERENCE RANGES

80 - 120 µg/dL

These values are for guidance only and it is recommended that each laboratory establishes its own reference range.

MATERIALS REQUIRED TO CARRY OUT THE TEST

- Spectrophotometer or photometer for reading at 578 nm (570-582).
- Water bath, thermostated at 37 °C.
- Glass and / or automatic pipettes.
- Clock or Stopwatch.
- Test tubes

ALERTS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Disposal, Safety and First Aid information are described in the Individual Product Safety Data Sheet (MSDS) for this product, available at www.biotechind.com.br or by telephone (35) -3214-4646.
- Dispose of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS).

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@biotechind.com.br

AUTOMATION

This product is compatible to the most types of biochemical automatic analysers. The applications are available at www.biotechind.com.br

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- Commission Directive 1999/11/EC of 8 March 1999 adapting to technical progress the principles of good laboratory practice as specified in Council Directive 87/18/EEC on the harmonisation of laws, regulations and administrative provisions relating to the application of the principles of good laboratory practice and the verification of their applications for tests on chemical substances (OJ L 77, 23.03.1999).
- MAKINO, T. A sensitive, direct colorimetric assay of serum zinc using Nitro-PAPS and microwell plates. *Clinica Chimica Acta.* n. 197, p. 209-220, 1991.
- GIROUX, E.L.; SCHECHTER, P. A Study of Zinc Distribution in Human Serum. *Bioinorganic Chemistry.* n. 5, p. 211-218, 1976
- WESTGARD, J. O. *et al.* A multi-rule shewhart chart quality control in clinical chemistry. *Clin. Chem.* v. 27, p.493-501, 1981.
- YOUNG, D.S. *Effects of drugs on clinical laboratory tests.* vol. 2, 5 ed. Washington DC: AACC Press, 2000.

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