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**ENGLISH**

**INTENDED USE**

Kit for the determination of C-reactive Protein (CRP) in serum. Use in *in vitro* diagnostic.

**STORAGE AND HANDLING**

- Store at 2 to 8 ° C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

**WORKING PRINCIPLE**

Latex PCR is a plate agglutination method for qualitative and semi-quantitative detection, where the latex particles coated with anti-human PCR antibodies are bound by the PCR present in the sample.

**SAMPLES – TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION**

**Sample Type:** Serum.  
**Collection, handling and preparation:** Perform the sample collection according to Good Clinical Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.

**Preservation:** Serum PCR is stable for 7 days at 2 to 8 °C and for 3 months when stored at -20 °C.

**PRODUCT DESCRIPTION**

	Suspension of Latex particles coated with goat anti-human IgG, pH 8.2. Preservative. The sensitivity of RLAX is traceable to ERM-DA 472/IFCC (International Federation of Clinical Chemistry)	
	Human serum with a concentration of CRP> 20 mg / L. Preservative. (INCLUDED IN THE PCR LATEX KIT)	
	Human serum in buffer Phosphate, stabilizer and preservative. (INCLUDED IN THE PCR LATEX KIT)	

**STABILITY IN USE**

- After opening, the product (RLAX, CONTROL (+) and CONTROL (-)) in use is stable until the expiration date printed on the label, since recommended storage conditions (2 to 8 ° C) are followed.
- Reagents should remain outside the specified temperature only for the time required to perform the tests.

**TREATMENT AND HANDLING OF THE PRODUCT**

**A) REAGENT PREPARATION**

RLAX: Ready for use.  
CONTROL (+) e CONTROL (-): Ready for use.

**B) OPERATING INTERVAL FOR PCR LÁTEX RLAX**

Up to 300 mg/L

**QUALITY CONTROL**

- The use of controls should be routine practice in the laboratory. It is recommended to use Positive and Negative Controls to monitor the performance of the procedure and to serve as a comparative model for a better interpretation of the results. For Internal Quality Control Laboratory it is recommended to use the control sera below:  
CONTROL (+) (included in the kit)      **REF**      30.005.00  
CONTROL (-) (included in the kit)      **REF**      30.005.00

**TEST PROCEDURE, CALCULATION AND INTERPRETATION**

**A) TEST PROCEDURE**

**Qualitative Method**

1. Before starting the tests, allow the RLAX, controls and samples to reach room temperature. Homogenize before use.
2. On the reaction plate use a different circle for each sample or control.
3. Place 20 µL of the sample / controls.
4. In each circle add 20 µL of RLAX.
5. Mix each determination using disposable rods (one for each sample).
6. Homogenize the plate for 2 minutes in circular motions, tilting it back and forth. Observe the result under artificial light: Presence or absence of agglutination.

**Semi-quantitative method**

1. Dilute the sample with saline solution (0.9% NaCl) as follows:

Dilution	Saline (µL)	Sample (µL)	Dilution Factor
1/2	100	100 (pure)	2
1/4	100	100 (sol.1/2)	4
1/8	100	100 (sol.1/4)	8
1/16	100	100 (sol.1/8)	16

2. Perform the test of each dilution according to the qualitative method.

**B) CALCULATIONS**

The approximate concentration in the semi-quantitative test is obtained by dilution factor by the sensitivity of the test (6 mg / L).

Example:  
Agglutination up to 1/8 dilution.  
Dilution Factor = 8  
PCR ≈ 48 mg / L (8 x 6)

**C) INTERPRETATION**

C-Reactive Protein is an acute phase protein that can significantly increase in most bacterial and viral infectious processes, inflammation and malignant neoplasms. Its elevation occurs rapidly and can reach values of 300 mg / L in 12 to 24 hours.

**INTERFERING AND LIMITATIONS**

**Hemolysis, Jaundice and Lipemia:** Hemoglobin > 1000 mg / dL, Bilirubin > 40 mg / dL, Triglycerides > 2000 mg / dL and Rheumatoid Factor > 40 IU / mL may interfere with the determination.

**Medications:** consult recommended literature reference (Young, 2000)

## PERFORMANCE CHARACTERISTICS

**Sensitivity:** Detection limit: 6 mg / L.

**Pro-zone effect:** Pro-zone high-dose effect is not observed up to concentrations of 300 mg / L.

**Analytical Specificity:** The product specifically determines the PCR in the presence of interferents in the sample, up to the concentrations reported above.

**Accuracy:** To estimate the concordance percentages, the method was compared with a similar one by the determination of 100 samples, which were obtained:

Percentage of Positive Agreement: 96%

Negative Agreement Percentage: 100%

**Precision:** In order to estimate the intra-assay precision, the study was performed in two daily runs, in duplicate, for 10 days and for the inter-run one daily run of two replicates, for 5 days. In both, four samples were used, two positive and two negative with 100% of concordance .

## RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange reagent bottle caps to avoid cross-contamination.
- Do not use the reagent when it has a visual characteristic that is not in accordance with the product's MSDS requirements.
- Avoid leaving the reagents out of the specified storage conditions when they are not in use.
- Use specific disposable tips for each sample, control and reagent.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the used material are fundamental factors for the stability of the reagents and obtaining correct results.

## REFERENCE RANGES

Until 6 mg/L

These values are for guidance only and it is recommended that each laboratory establishes its own reference range.

## MATERIALS REQUIRED TO CARRY OUT THE TEST

- Reaction plate.
- Disposable rod.
- Automatic pipette.
- Clock or Stopwatch

## ALERTS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Disposal, Safety and First Aid information are described in the Individual Product Safety Data Sheet (MSDS) for this product, available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or by telephone (35) -3214-4646.
- Dispose of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS).

## QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiry date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: [sac@biotechnicaltda.com.br](mailto:sac@biotechnicaltda.com.br)

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- BURTIS, C.A.; ASHWOOD, E.R.; BRUNS, D.E. Tietz: Fundamentos de Química Clínica. 6. ed. Rio de Janeiro: Saunders Elsevier, 2008. 959 p.
- HOKAMA, Yoshitsugi; NAKAMURA, Robert M. C-Reactive Protein: Current Status and Future Perspectives. *Journal Of Clinical Laboratory Analysis*, v. 1, p.15-27, 1987.
- OTSUIJI, Shogo; SHIBATA, Hideaki; UMEMA, Mamoru. Turbidimetric Immunoassay of Serum C-reactive Protein. *Clin. Chem.* v. 28, p.2121-2124, 1982.
- PEPYS, M. B.; DASH, A. C.; ASHLEY, M. J. Isolation of C-reactive protein by affinity chromatography. *Clin. Exp. Immunol.* v. 30, p.32-37, 1977.
- WADSWORTH, Charles; WADSWORTH, Elsa. Efficacy of latex agglutination and quantification methods for determination of C-reactive protein (CRP) in pediatric sera. *Clinica Chimica Acta*, v. 138, p.309-318, 1984.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. *Clin. Chem.* v.27 p.493-501, 1981.
- YAMAMOTO, S. et al. Avidity of antibody and agglutinability of antibody-sensitized latex in latex agglutination test. *Veterinary Immunology And Immunopathology*, v. 36, p.257-264, 1993.
- YOUNG, D.S. Effects of drugs on clinical laboratory tests - vol. 2, 5 ed. Washington DC: AACCC Press, 2000.

## TABELA DE SÍMBOLOS INTERNACIONAIS / TABLE OF INTERNATIONAL SYMBOLS / TABLA DE SÍMBOLOS INTERNACIONALES

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