VDRL - Sífilis

VDRL Syphilis / VDRL Sifilis

BT3001200 (R1 / CONTROL+ / CONTROL-) BT3001200 (R1)

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Dr. Gilson Serio Pizzo

CRF MG – 5310

Anvisa 80027310320

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ENGLISH

BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE

Kit for the determination of non-treponemal antibodies (reagins) in serum, plasma and cerebrospinal fluid (CSF). *In vitro* diagnostic use.

STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing. Do not freeze.
- Once opened, the product is stable for 8 weeks, as long as the recommended storage conditions (2 to 8 °C) are followed.
- Do not use reagents whose shelf life has expired.
- Ready-to-use reagent and controls.

WORKING PRINCIPLE

Method: Plate flocculation

VDRL is a non-treponemal qualitative or semi-quantitative method of plate flocculation to detect non-treponemal antibodies (reagins). The antigenic suspension, composed by a mixture of complex lipids, is agglutinated in the presence of non-treponemal antibodies from the sample of a patient infected by *Treponema pallidum*.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum, plasma citrate, and cerebrospinal fluid.

Collection and handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Stability Period
Serum and Plasma	2 to 8 °C	7 days
	-20 °C	3 months

Cerebrospinal fluid (CSF) samples can be stored in the refrigerator (2-8°C) for up to 5 days. For longer periods, store at -20°C.

PRODUCT DESCRIPTION

R 1

CONTROL+

CONTROL-

Cardiolipin 0,3 g/L; lecitin 2,1 g/L; cholesterol 9 g/L; phosphate buffer 1,5 mmol/L; preservative 0,95 g/L. Traceable to reference material 05/132 from the National Institute for Biological Standards and Control (NIBSC).

Artificial serum in 170 mmol/L NaCl solution; BSA 10 g/L; metanil yellow 10 g/L; reagin with titer >1/8; preservative 0,95 g/L.

Animal serum in 150 mM glycine buffer; BSA 5 g/L; preservative 0,95 g/L.

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratorial quality control, it is recommended the use of the controls below:

VDRL – Syphilis CONTROL (+)	REF	BT3001200
VDRL – Syphilis CONTROL (-)		BT3001200

NECESSARY EQUIPMENT FOR TESTING

- Transparent glass plate;
- Pipette, clock or stopwatch;
- Orbital shaker adjustable to 180 rpm and optical microscope.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

Qualitative Method

1. Allow the reagent to reach room temperature. Homogenize before use.
2. Pipette 50 µL of the sample into the plate's well and then add 20 µL of the R1.
3. Homogenize the plate in the orbital shaker at 180 rpm for 4 minutes. Immediately after this period, read the results with an optical microscope. False-positive results might happen if the test is read after 4 minutes from the end of the agitation period.

Semi-quantitative Method

1. Perform a serial dilution of the sample with saline solution (NaCl 9 g/L).
2. Proceed to the test determining each dilution as stated in the Qualitative Method.

B) READING AND INTERPRETATION

Evaluate, with an optical microscope, the presence or absence of agglutination immediately after the homogenization period. The results should be interpreted as follows:

Visual Observation	Reading	Result
Large or medium-sized agglutination	P	Positive
Small-sized agglutination	WP	Weakly positive
No agglutination or a “dust”	N	Negative

For the semi-quantitative method, the titer is defined as the highest titration with a positive result.

C) INTERPRETATION

Syphilis is a sexually transmitted infection, of chronic character, caused by the bacteria *Treponema pallidum*. The course of the disease can be categorized in four stages: primary, secondary, latent and tertiary. In the primary stage, between 10 and 90 days after the infection, a sore called chancre is developed. After 6 weeks to 6 months from the cicatrization of the chancre, the secondary stage of the disease begins. It is characterized by a rash, mainly in thorax region, and skin lesions in the mucosa. Later, syphilis reaches the latent stage, in which no symptoms are observed. The tertiary stage happens between 1 and 40 years after the beginning of the infection, in around 15 to 25% of non-treated infections, being characterized by tissue damage, mainly in the nervous and cardiovascular systems. The syphilis infection causes the immunological system to produce treponemal and non-treponemal antibodies. The non-treponemal antibodies (reagins) recognize as antigens the tissue fragments from the lesions caused by the infection. These antibodies can also be generated by other diseases; thus, non-treponemal methods, such VDRL, should be use for the diagnostic screening of syphilis. The treponemal antibodies, on the other hand, are specific for *Treponema pallidum*; thus, treponemal tests, such TPHA, are confirmatory for the syphilis diagnostic.

PERFORMANCE CHARACTERISTICS

Analytical Sensitivity	Serum and plasma	Cerebrospinal fluid
Limit of Detection	0,375 IU/mL	1:32

Analytical Specificity				
Serum and plasma				Cerebrospinal fluid
Hemoglobin	Bilirubin	Triglycerides	Rheumatoid Factor	Hemoglobin
2000 mg/dL	20 mg/dL	1000 mg/dL	500 UI/mL	2000 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	Serum and Plasma	Cerebrospinal fluid
Number of Samples	80 samples	71 samples
Diagnostic Sensitivity	92,5%	100%
Diagnostic Specificity	100%	100%

Prozone High Dose Effect:

Serum and plasma: There is no prozone effect up to titers of 1/128.

Cerebrospinal fluid: There is no prozone effect up to titers of 1/32.

Precaution:

The intra-assay studies were done in two daily runs, with 2 replicates for 20 days; the inter-assay studies were done with 5 lots, with 3 replicates in two analytical runs.

Samples	Intra-assay Precision	Inter-assay Precision
	Concordance (%)	Concordance (%)
Negative	100%	100%
Positive 1/8	100%	100%
Positive 1/16	100%	100%

BIOTÉCNICA IND.COM. LTDA - CNPJ: 025340690001/20, Avenida Washington Ribeiro, 200 - Industrial Miguel de Luca
CEP 37072-030 Varginha-MG BRASIL Tel/fax: +55 35 3214 4646 www.biotecnica.ind.br

Revisão/Revision/Revisión: 03 – 05/11/2024

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- In addition to the undiluted sample, a test with a 1:8 dilution with 0.9% NaCl (for all samples) is recommended in order to avoid the prozone effect.
- VDRL is non-specific for syphilis.All reactive samples should be retested with treponemic methods such as TPHA and FTA-Abs to confirm the results.
- A non-reactive result by itself does not exclude a diagnosis of syphilis.
- False positive results have been reported in diseases such as infectious mononucleosis, viral pneumonia, toxoplasmosis, and autoimmune diseases and in pregnant women.
- CSF samples contaminated with blood can generate false-positive results.
- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotechnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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