daily run of four replicates, for 5 days. In both, four samples were used, two positive and two negative with 100% agreement.

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange reagent bottle caps to avoid cross-contamination.
- Do not use the reagent when it has a visual characteristic that is not in accordance with the product's MSDS requirements.
- Avoid leaving the reagents out of the specified storage conditions when they are not in use.
- Use specific glass pipettes and disposable tips for each sample, control, standard / calibrator and reagent.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the used material are fundamental factors for the stability of the reagents and obtaining correct results.

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MATERIALS REQUIRED TO CARRY OUT THE TEST

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ALERTS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Disposal, Safety and First Aid information are described in the Individual Product Safety Data Sheet (MSDS) for this product, available at www.biotechind.com.br or by telephone (35) -3214-4646.
- Dispose of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS).

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use, all BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, please contact the Scientific Advisory Board of Biotécnica Ltda, by phone +55 35 3214 4646 or by email sac@biotechind.com.br.

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- ASSIMEH, S. N.; JOHNSON, P. M. A Haemagglutination Method for Detection of Rheumatoid Factor Using Preserved Erythrocytes Covalently Coated With Human IgG. *Journal Of Immunological Methods*, v. 34, p.205-215, 1980.
- DORNER, Robert W.; ALEXANDER JUNIOR, Roy L; MOORE, Terry L. Rheumatoid factors. *Clinica Chimica Acta*, v. 167, p.1-21, 1987.
- FULFORD, Karen M.; TAYLOR, Roger N.; PRZYBYCZEWSKI, Vncent A. Reference Preparation to Standardize Results of Serological Tests for Rheumatoid Factor. *Journal Of Clinical Microbiology*, v. 7, n° 5, p.434-441, 1978.
- KORITZ, T. N. et al. Measure of Rheumatoid Factor in Human Sera by Passive Haemagglutination of Human Erythrocytes Carrying Immunoglobulin Linked by Chromic Chloride. *Journal Of Immunological Methods*, v. 32, p.1-9, 1980.
- ROSE, Harry M. et al. Differential Agglutination of Normal and Sensitized Sheep Erythrocytes by Sera of Patients with Rheumatoid Arthritis. *Proc. Soc. Exp. Biol. Med.* v. 68, n. 1, p.1-6, 1948.
- SHMERLING, Robert H.; DELBANCO, Thomas L. The Rheumatoid Factor: An Analysis of Clinical Utility. *The American Journal Of Medicine*, v. 91, p.528-534, 1991.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. *Clin. Chem.* v.27 p.493-501, 1981.
- WOLFE, Frederick; CATHEY, Mary Ann; ROBERTS, F. Kay. The Latex Test Revisited: Rheumatoid Factor Testing in 8,287 Rheumatic Disease Patients. *Arthritis And Rheumatism*, v. 34, n° 8, p.951-960, 1991.
- YOUNG, D.S. Effects of drugs on clinical laboratory tests - vol. 2, 5 ed. Washington DC: AAC Press, 2000.

TABELA DE SÍMBOLOS INTERNACIONAIS / TABLE OF INTERNATIONAL SYMBOLS / TABLA DE SÍMBOLOS INTERNACIONALES

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