NEFA VET

NEFA VET / NEFA VET
Ref. 90.054.00

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ENGLISH

INTENDED USE

Kit intended to determine non-esterified fatty acids (NEFA) in serum and plasma samples. *In vitro* diagnostic use only.




STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagents ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8°C) are followed.
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE

Method: Enzymatic - Trinder

Non-esterified fatty acids (NEFA) from sample react with coenzyme A (CoA), in the presence of Acyl-CoA Synthetase (ACS), to form Acyl-CoA. The latter is oxidized by Acyl-Coenzyme A Oxidase (ACOD), generating hydrogen peroxide. A Trinder compound, in the presence of hydrogen peroxide and under the catalytic action of peroxidase (POD), is converted to a colored product that can be spectrophotometrically measured at 546 nm. The intensity of the color is proportional to the free fatty acids concentration in the sample.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum and plasma (EDTA and citrate).

Collection and handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Stability Period
Serum and Plasma*	4 to 8 °C	1 day
	-20 °C	10 days

* Samples should be collected after, at least, 12 hours fasting. It is recommended to perform the assay immediately after the collection, as the fatty acids concentration is increased by the hydrolytic action of some enzymes. Samples from patients being treated with heparin should not be used as they can promote wrongly elevated results.

PRODUCT DESCRIPTION

R 1 Phosphate buffer ≥ 50 mmol/L; Coenzyme A ≥ 0,9 mmol/L; ATP ≥ 5 mmol/L; Acyl-CoA Synthetase ≥ 0,3 KU/L; MgCl₂ ≥ 2,0 mmol/L; 4-aminoantipyrine ≥ 1,5 0 mmol/L; stabilizers; activators; preservatives. **X**

R 2 Phosphate buffer ≥ 50 mmol/L; Acyl-CoA Oxidase ≥ 10 KU/L; Peroxidase ≥ 7,5 KU/L. **X**

STD Concentration equivalent to ≈ 1,00 mmol/L of NEFA; stabilizers; preservatives. **X**

CONTROL 1 Concentration equivalent to ≈ 0,45 mmol/L of NEFA; stabilizers; preservatives. **X**

CONTROL 2 Concentration equivalent to ≈ 0,95 mmol/L of NEFA; stabilizers; preservatives. **X**

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and controls below:

STD – NEFA STD	90.054.00
Control Level 1 – NEFA CONTROL 1	REF 90.054.00
Control Level 2 – NEFA CONTROL 2	90.054.00

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 546 nm.
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Pipette in the test tubes:

	Blank	STD	Sample
STD	-	20 µL	-
Sample	-	-	20 µL
R1	800 µL	800 µL	800 µL
Homogenize and incubate at 37°C for 5 minutes.			
Measure the Sample's and Calibrator's A1 absorbance at 546-550 nm.			
R2	200 µL	200 µL	200 µL

2. Homogenize and incubate at 37°C for 5 minutes.

3. Measure the Sample's and Calibrator's A2 absorbance at 546 nm.

B) CALCULATIONS

NEFA (mmol/L) = $\frac{\text{Sample's (A2 - A1)}}{\text{Calibrator's (A2 - A1)}}$ x STD Concentration (mmol/L)

Calculations with the Calibration Factor

Calibration Factor = $\frac{\text{STD value (mmol/L)}}{(A2 - A1) \text{ STD}}$

NEFA (mmol/L) = Sample's (A2 - A1) x Calibration Factor

Automação: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotecnicaitda.ind.br

C) INTERPRETATION

Non-esterified fatty acids are used by the organism as a metabolic energy source, as substrate for structures of the cellular membrane and as a precursor for several intracellular signaling molecules, for example, the prostaglandins. Non-esterified fatty acids are usually formed by lipolysis of the adipose tissue, which is affected by diet and fluctuations of the insulin levels. The increase in NEFA concentrations in pathological states such as insulin resistance/type 2 diabetes, adiposity, malignant disease, and metabolic syndrome is associated to the development of cardiovascular diseases. Free fatty acids in the bloodstream can reflect the organism's fat metabolism.

PERFORMANCE CHARACTERISTICS

Operating range
0.02 to 2.78 mmol/L

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0,9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity - Quantification Limit				
0,02 mmol/L				
Analytical Specificity				
Species	Hemoglobin	Bilirubin	Triglycerides	Ascorbic Acid
Caninos	500 mg/dL	40 mg/dL	2000 mg/dL	18 mg/dL
Felinos	500 mg/dL	30 mg/dL	250 mg/dL	0 mg/dL
Equinos	500 mg/dL	10 mg/dL	100 mg/dL	0 mg/dL
Bovinos	500 mg/dL	18 mg/dL	500 mg/dL	10 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy - Canines	
Number of Samples	20 em duplicate
Regression Equation	y = 0.9945x + 0.0439
Correlation Coefficient (R)	0.9940
Accuracy - Felines	
Number of Samples	10 em duplicate
Regression Equation	y = 1.0465x - 0.0082
Correlation Coefficient (R)	0.9946
Accuracy - Equines	
Number of Samples	10 em duplicate
Regression Equation	y = 0.9355x + 0.0191
Correlation Coefficient (R)	0.9994
Accuracy - Bovines	
Number of Samples	10 em duplicate
Regression Equation	y = 0.9902x + 0.0223
Correlation Coefficient (R)	0.9963

Precision:

The studies were performed in two runs per day, in duplicate.

Intra-assay Precision				
Sample	Repetitions	Average mmol/L	Intra-assay Precision	
			DP	%CV
Low	20	0,42	0,02	3,7%
High	20	0,97	0,02	2,5%

Inter-assay Precision				
Sample	Repetitions	Average mmol/L	Inter-assay Precision	
			DP	%CV
Low	40	0,44	0,02	4,6%
High	40	0,99	0,03	3,0%

CV: Coefficient of variation; SD: Standard deviation.

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

Canines	0,13 – 1,25 mmol/L
Felinos	0,15 – 1,05 mmol/L
Equinos	0,02 – 0,43 mmol/L
Bovinos	0,07 – 0,46 mmol/L

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the International System of Units (mmol/L):

NEFA (mg/dL) = NEFA (mmol/L) / 0,0354

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotechnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnicaltda.com.br.

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