FINALIDADE

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ENGLISH

BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE

Kit for the qualitative determination of human chorionic gonadotrophin (HCG) in serum and urine. *In vitro* diagnostic use.

STORAGE AND HANDLING

- Store at 4 to 30 °C.
- The product is stable until the expiration date printed on the label, as long as the recommended storage conditions are followed.
- After opening the individual packaging the strip is stable for 1h under the following conditions: humidity from 20 to 90% and temperature from 10 to 50 ° C
- Ready-to-use strip.

WORKING PRINCIPLE

Method: Immunochromatographic
Nitrocellulose membrane chromatographic immunoassay coated with anti-HGC antibodies in the test region and goat anti-mouse antibodies in the control region. The sample migrates through the membrane by capillarity. In the presence of HCG, a colored band will form in the test region. As a control of the procedure, a colored band should always appear in the control region.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum and urine
Collection and handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material. Preferably use the first morning urine. Samples with visible precipitates must be centrifuged.

	Temperature	Stability Period
Serum	4 to 8 °C	7 days
	-20 °C	3 months
Urine	2 to 8 °C	2 days
	-18 °C	4 months

PRODUCT DESCRIPTION

TIRA	Single-use nitrocellulose strip coated with human anti-HCG antibodies and rat goat anti-HCG antibodies.
	The sensitivity of the strip is traceable to NIBSC 75/551 (National Institute for Biological Standards and Control).

QUALITY CONTROL

The strip has a control line that should be visible for the interpretation of the result.

NECESSARY EQUIPMENT FOR TESTING

- Test tubes;
- Clock or stopwatch.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

Before opening the individual packs, leave the strips to be used at room temperature for 30 minutes (20 to 30 ° C).

Urine

Immerse for 5 to 8 seconds the lower end (indicated by the arrows) in the sample, without exceeding the limit line (MAX). Remove the strip and place it on a dry, clean, flat surface with the reaction area facing up. Turn on the timer and observe the result within 5 minutes.

Serum

Immerse for 10 seconds the lower end (indicated by the arrows) in the sample, without exceeding the limit line (MAX). Remove the strip and place it on a dry, clean, flat surface


with the reaction area facing up. Activate the timer and observe the result within 5 to 7 minutes.

B) ANALYSIS OF RESULTS

Positive: two colored lines appear (red or pink), one line in the control region (C) and another in the test region (T).

Negative: one colored line appears only in the control region (C).

Invalid: no colored line appears in the control region (C).

C) INTERPRETATION
HCG is a glycoprotein hormone synthesized in placental syncytiotrophoblast cells. It stimulates the corpus luteum to synthesize progesterone to maintain the endometrium suitable for fixation and implantation of the embryo. Minimum quantities are also produced in the pituitary gland of men and non-pregnant women..

PERFORMANCE CHARACTERISTICS

Operating range	
Serum and Urine	Up to 50,000 mUI/mL
Sensitivity	
Detection Limit	
Serum and Urine	25 mUI/mL

Analytical Specificity		
LH	FSH	TSH
500 mUI/mL	1000 mUI/mL	1000 µUI/mL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	Serum	Urine
Number of Samples	50	100
Positive Agreement	95%	100%
Negative Agreement	100%	96%

Prozone High Dose Effect:
There is no prozone effect up to 50,000 mUI/mL de HCG.

Precision:

It was determined using a sample with concentration of 25 mUI/mL of HCG in 20 replicates with three different batches. The results obtained were in agreement with those expected

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Only use strips whose packages are perfectly sealed.
- The test reading should not exceed 5 minutes for urine and 7 minutes for serum.
- The result of this test is not confirmatory for verification of pregnancy. As with any diagnostic procedure, the result of this test should always be interpreted by a qualified professional, armed with other available clinical information.
- Concentrations below the detection limit may cause negative results. In this case, a new test must be performed after two days to confirm the result.
- hBhCG values less than 50 mUI/mL may cause weakly positive results, characterized by a faint pinkish to light red coloring line in the test line. The test should be interpreted in conjunction with the patient's signs and symptoms and other confirmatory laboratory tests.
- Patients with non-gestational trophoblastic diseases (carcinoma, hydatidiform mole and others) may have a positive result in the absence of pregnancy.
- This test does not differentiate a normal pregnancy from an ectopic pregnancy.
- Low concentrations of hBhCG may occur in normal non-pregnant women.
- False negative result can be obtained in cases of ectopic pregnancy, imminent abortion, pregnancy toxemia, due to possible decreased excretion of hBhCG.
- The presence of heterophile antibodies can produce false positive results. In the presence of positive results without supporting clinical data, perform a quantitative test.
-

WARNINGS AND PRECAUTIONS
- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material. - The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotechnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE
- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br. - To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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