βHCG*
βHCG* | βHCG*
Ref. 30.010.00

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ENGLISH

INTENDED USE

Kit for the qualitative determination of human chorionic gonadotrophin (HCG) in serum and urine. *In vitro* diagnostic use.

STORAGE AND HANDLING CONDITIONS

- Store at 4 to 30 °C.
- The validity of the kit is printed on the package label.
- Do not use the reaction strip if expire date is over.

WORKING PRINCIPLE

Nitrocellulose membrane chromatographic immunoassay coated with anti-HCG antibodies in the test region and goat anti-mouse antibodies in the control region. The sample migrates through the membrane by capillarity. In the presence of HCG, a colored band will form in the test region. As a control of the procedure, a colored band should always appear in the control region.

SAMPLES: TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION

Sample Type: Serum and Urine.

Collection, handling and preparation: Perform the sample collection according to Good Clinical Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.

Preservation: **Serum:** HCG is stable for 1 day if stored at 20 to 25 °C, 7 days at 4 to 8 °C and 12 months at -20 °C.

Urine: Urine HCG is stable for 2 days if stored at 2 to 8 °C and 4 months at -18 °C.

Preferably use the first morning urine. Samples with visible precipitates must be centrifuged.

PRODUCT DESCRIPTION

Single-use nitrocellulose strip coated with human anti-HCG antibodies and rat goat anti-HCG antibodies.
The sensitivity of the strip is traceable to NIBSC 75/551 (National Institute for Biological Standards and Control).

STABILITY IN USE

- After opening the individual packaging the strip is stable for 1h under the following conditions: humidity from 20 to 90% and temperature from 10 to 50 °C.

PRODUCT TREATMENT AND HANDLING

A) REAGENT PREPARATION

Strip: Ready for use.

B) OPERATING INTERVAL

Serum / Urine: Up to 50.000 mIU/mL.

QUALITY CONTROL

The strip has a control line that should be visible for the interpretation of the result.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Before opening the individual packs, leave the strips to be used at room temperature for 30 minutes (20 to 30 °C).

1.1- **Urine:** Immerse for 5 to 8 seconds the lower end (indicated by the arrows) in the sample, without exceeding the limit line (MAX). Remove the strip and place it on a dry, clean, flat surface with the reaction area facing up. Turn on the timer and observe the result within 5 minutes.


1.2- **Serum:** Immerse for 10 seconds the lower end (indicated by the arrows) in the sample, without exceeding the limit line (MAX). Remove the strip and place it on a dry, clean, flat surface with the reaction area facing up. Activate the timer and observe the result within 5 to 7 minutes.

B) ANALYSIS OF RESULTS

POSITIVE: two colored lines appear (red or pink), one line in the control region (C) and another in the test region (T).

NEGATIVE: A colored line appears only in the control region (C).

INVALID: no colored line appears in the control region (C).

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HCG is a glycoprotein hormone synthesized in placental syncytiotrophoblast cells. It stimulates the corpus luteum to synthesize progesterone to maintain the endometrium suitable for fixation and implantation of the embryo. Minimum quantities are also produced in the pituitary gland of men and non-pregnant women.

INTERFERENTES OR LIMITATIONS

Specificity: LH < 500 mIU/mL, FSH < 1000 mIU/mL, TSH < 1000 µIU/mL did not interfere.

Medications: consult recommended literature reference (Young, 2000).

PERFORMANCE CHARACTERISTICS

Sensitivity: **Serum / Urine:** Detection limit: 25 mIU/mL.

High dose Prozone effect: **Serum / Urine:** No high dose prozone effect at concentrations up to 50,000 mIU / mL of HCG is observed.

Analytical Specificity: The product specifically determines HCG in the presence of interferents in the sample, up to the concentrations reported above.

Accuracy: Serum: To estimate the percentages of concordance, the method was compared with a similar one by the determination of 50 samples, obtaining:

Percentage of Positive Agreement: 95%

Negative Agreement Percentage: 100%

Urine: To estimate the percentages of agreement, the method was compared with a similar one by the determination of 100 samples, which were obtained:

Percentage of Positive Agreement: 100%

Negative Agreement Percentage: 96%

Accuracy: It was determined using a sample with concentration of 25 mIU / mL of HCG in 20 replicates with three different batches. The results obtained were in agreement with those expected.

RESIDUAL RISKS, CARE AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not use the product when it has a visual characteristic that is not in accordance with the product's MSDS.
- Don't leave the product out of the specified storage conditions.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Only use strips whose packages are perfectly sealed.
- The test reading should not exceed 5 minutes for urine and 7 minutes for serum.

MATERIALS REQUIRED TO CARRY OUT THE TEST

- Clock or Stopwatch.
- Test tubes

ALERTS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Disposal, Safety and First Aid information are described in the Individual Product Safety Data Sheet (MSDS) for this product, available at www.biotechnica.ind.br or by telephone (35) -3214-4646.
- Disposal of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS).

QUALITY / SAC GUARANTEE – CUSTOMER SERVICE

Before being released for consumption, all Biotécnica reagents are tested by the Quality Control Department. The quality of the reagents is ensured by the expiry date mentioned on the presentation packaging, since they are stored and transported under specified conditions. The data related to the Quality Control of this product (lot printed on the labels of the reagent bottles) or any doubt on the use of this kit, contact the Scientific Advisory Board of Biotécnica Ltda, by phone +55 35 3214 4646 or by email sac@biotecnica Ltda.com.br.

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