Frutosamina Bi
Fructosamine Bi | Fructosamina Bi
Ref. 10.036.00

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ENGLISH

INTENDED USE

Kit for determination of Fructosamine in serum and plasma. Use in in vitro diagnostic.

STORAGE AND HANDLING

- Store in a temperature of 2 to 8°C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE

Fructosamine is formed by the non-enzymatic glycation of plasma proteins. In alkaline medium, fructosamine's keto-enol balance is shifted, so that the analyte can now reduce the nitrotetrazolium blue dye into a purple formazan which is spectrophotometrically measured at 550 nm. The difference between the absorbances taken within 5 and 10 minutes of incubation at 37°C is proportional to the sample's concentration of fructosamine.

SAMPLES: TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION

Sample Type: Serum and Plasma (Heparin and EDTA).

Collection, handling and preparation: Perform the sample collection according to Good Clinical Laboratory Practices. The samples to be analyzed should be treated as potentially infectious material.

Preservation: Fructosamine in serum and plasma is stable for 3 days if stored at 20 to 25°C, 2 weeks at 4 to 8°C and 2 months if stored at -20°C.

PRODUCT DESCRIPTION

| | | |
|-----|---|---|
| R 1 | Buffer ≥ 120 mmol / L, nitrotetrazolium blue (NBT) ≥ 100 µmol / L, Uricase ≥ 1000 U / L, preservative, surfactants and stabilizers. | X |
| | | |

| | | |
|-----|---|-----------|
| R 2 | Buffer > 1,0 mmol/L | Biohazard |
| | | |
| CAL | Lyophilized serum calibrator. Value printed on the bottle label. Traceable to a Polylysine standard. | |

STABILITY IN USE

- **R1 and R2:** Once opened, the reagents in use are stable to the expiration date printed on the label, provided that it is followed by recommended storage conditions (2 to 8°C).
- **CAL:** After reconstitution, the calibrator is stable for 15 days if stored at 2 to 8°C and 2 months at -20°C.
- Reagents should remain outside the specified temperature only for the time required to perform the tests.

TECHNICAL PREPARURE

A) REAGENT PREPARATION
R1: Reagent ready for use.
R2: Reagent ready for use.

CAL: Lyophilized. Reconstitute with 1.0 mL of purified water. Mix by gentle inversion to avoid foaming.

B) OPERATING RANGE

The operating range of the product is 25,40 a 807,00 µmol/L. For higher values, dilute the sample with 150 mM NaCl (0.9%), carry out a new dosage and multiply the result obtained by the dilution factor.

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. It is suggested to use one control in the reference range or in the decision level and another control with value in another range of clinical significance. For Calibration and Internal Control of Laboratorial Quality it is recommended to use the calibrator serum and the control sera below:

| | | |
|--|------------|-----------|
| Normal Control Serum - Quantinorm | REF | 13.003.00 |
| Serum Pathological Control - Quantialt | | 13.004.00 |

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Pipette in the assay tubes:

| | |
|---------------------|---------------|
| | Volume |
| R1 | 800 µL |
| R2 | 200 µL |
| Sample / CAL | 50 µL |

2. Homogenize and insert into the cuvette holder thermostated at 37 ° C, zeroing the equipment with purified water.
3. Measure the absorbance at 550 nm (520-580) at 5 minutes (A1) and at 10 minutes (A2).

B) CALCULATIONS

Fructosamine (µmol/L) = $\frac{(A_2 - A_1) \text{ Sample}}{(A_2 - A_1) \text{ Calibrator}}$ x Calibrator value (µmol/L)

Example:
A₁ Calibrator = 0.082 / A₂ Calibrator = 0.197 / A₁ Sample = 0.136 / A₂ Sample = 0.212 / Calibrator value = 350 µmol/L

Fructosamine (µmol/L) = $\frac{(0,212 - 0,136)}{(0,197 - 0,082)} \times 350 = 231,3 \mu\text{mol/L}$

With Calibration Factor:

Calibration Factor = $\frac{\text{Calibrator factor} (\mu\text{mol/L})}{(A_2 - A_1) \text{ Calibrator}}$
Fructosamine (µmol/L) = (A₂ - A₁) Sample x Calibration Factor

Example:
Calibration Factor = $\frac{350}{0,197 - 0,082} = \frac{350}{0,115} = 3043,48$

Fructosamine (µmol/L) = (0.212 - 0.136) x 3043.48 = 231.3 µmol/L

C) INTERPRETATION

Fructosamine is the generic name for ketoamines. The name refers to the structure of the rearrangement product of the ketoamine formed by the interaction of glucose with the ε-amino group on the lysine residues of albumin. As glycated hemoglobin determinations, the dosage of fructosamine can be used to monitor the average blood glucose concentration over a prolonged period.

INTERFERING AND LIMITATIONS

Endogenous interferers: Hemoglobin > 200.0 mg/dL / Bilirubin > 7.5 mg/dL / Triglycerides > 586.7 mg/dL / Ascorbic Acid > 8.0 mg/dL / Uric Acid > 20.0 mg/dL / Glucose > 1000.0 mg/dL interfere with the dosage.

Medications: consult recommended literature reference (Young, 2000).

PERFORMANCE CHARACTERISTICS

Sensitivity: Detection limit: 6.46 µmol/L / Limit of quantification: 20.0 µmol/L.
Analytical Specificity: The product specifically determines fructosamine in the presence of other interfering substances in the sample, up to the concentrations reported above.

Accuracy: The method was compared with a similar method by the determination of 40 samples in duplicate. The regression equation $Y = 0.987 + 1.48$ and correlation coefficient $R = 0.9941$ were obtained. Using this equation, the estimated total systematic error is -0.56% for a level of 200 µmol/L and -1.05% for a level of 600 µmol/L.

Precision: It was determined using samples at 03 decision levels, two runs in duplicate per day for 20 days, obtaining:

| Samples (µmol/L) | Repetition | Within-Run Precision | | Total Precision | |
|------------------|------------|----------------------|-----|-----------------|-----|
| | | SD (µmol/L) | %CV | SD (µmol/L) | %CV |
| 136.43 | 80 | 1.19 | 0.9 | 2.21 | 1.6 |
| 259.03 | 80 | 1.96 | 0.8 | 4.04 | 1.6 |
| 478.17 | 80 | 2.60 | 0.5 | 6.82 | 1.4 |

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with Good Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange the caps from different reagents in order to avoid cross contamination.
- Do not use the reagent if it displays any sign in disagreement with the one specified in the product MSDS.
- Use glass pipettes and disposable tips for each sample, control, standard / calibrator and reagent.
- The laboratory should establish chemical, microbiological and particulates requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
- The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.

REFERENCE RANGES

| | |
|------------------|------------------|
| Serum and Plasma | 205 - 285 µmol/L |
|------------------|------------------|

These values are for guidance only and it is recommended that each laboratory establishes its own reference range.

ADDITIONAL EQUIPMENT

- Spectrophotometer or photometer for reading at 550 nm (520-580) with thermostated cuvette at 37°C.
- Glass pipettes and/or automated.
- Clock or chronometer.
- Assay tubes.

WARNINGS AND PRECAUTIONS FOR DISPOSAL OF PRODUCTS

- Disposal, Safety and First Aid information are described in the Material Safety Data Sheet for this product, available at www.biotechnica.ind.br or by telephone (35) 3214-4646.
- Dispose of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS).

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being released for consumption, all Biotechnique reagents are tested by the Quality Control Department. The quality of the reagents is ensured by the expiration date mentioned on the presentation packaging once they are stored and transported under specified conditions. The data related to the Quality Control of this product (lot printed on the labels of the reagent bottles) or any doubt in the use of this kit, contact the Scientific Advisory of Biotechnica Ltda, by phone +55 35 3214 4646 or by email sac@Biotechnicaltda.com.br

AUTOMATION

This procedure is automated in most analyzers. The protocols are available at www.biotechnicaltda.ind.br

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- ARMBRUSTER, David A. Fructosamine: Structure, Analysis, and Clinical Usefulness. **Clinical Chemistry**, v. 33, n. 12, p.2153-2163, 1987.
- BAKER, John R. et al. Use of Protein-Based Standards in Automated Colorimetric Determinations of Fructosamine in Serum. **Clinical Chemistry**, v. 31, n. 9, p.1550-1554, 1985.
- HURST, Paul L. Effect of Anticoagulants on Fructosamine Determination. **Clinical Chemistry**, v. 33, n.10, p.1947-1948, 1987.
- KRUSE-JARRES, J. D. et al. A new colorimetric method for the determination of fructosamine. **Lab.med.**, v. 13, p.245-253, 1989.
- LIM, Yong Seng; STALEY, Michael J. Measurement of Plasma Fructosamine Evaluated for Monitoring Diabetes. **Clin. Chem.**, v. 31, nº. 5, p.731-733, 1985.
- BURTIS, Carl A.; ASHWOOD, Edward R.; BRUNS, David E. **Tietz: Fundamentos de Química Clínica**. 6. ed. Rio de Janeiro: Saunders Elsevier, 2008. 959 p.
- YOUNG, D.S. **Effects of drugs on clinical laboratory tests - vol. 2**. 5 ed. Washington DC: AAC Press, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. **Clin. Chem.** v.27 p.493-501, 1981.

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