Samples (µmol/L)	Replicates	Within-Run Precision		Total Precision	
		SD (µmol/L)	CV (%)	SD (µmol/L)	CV (%)
136,43	80	1.19	0.90	2.21	1.60
259,03	80	1.96	0.80	4.04	1.60
478,17	80	2.60	0.50	6.82	1.40

CV: Coefficient of variation; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

Serum and Plasma	205 - 285 µmol/L
------------------	------------------

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to Good Laboratory Practice in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotechnica.ind.br or calling for +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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- ARMBRUSTER, David A. Fructosamine: Structure, Analysis, and Clinical Usefulness. **Clinical Chemistry**, v. 33, n. 12, p.2153-2163, 1987.
- BAKER, John R. et al. Use of Protein-Based Standards in Automated Colorimetric Determinations of Fructosamine in Serum. **Clinical Chemistry**, v. 31, n. 9, p.1550-1554, 1985.
- HURST, Paul L. Effect of Anticoagulants on Fructosamine Determination. **Clinical Chemistry**, v. 33, n. 10, p.1947-1948, 1987.
- KRUSE-JARRES, J. D. et al. A new colorimetric method for the determination of fructosamine. **Lab.med.**, v. 13, p.245-253, 1989.
- LIU, Yong Seng; STALEY, Michael J. Measurement of Plasma Fructosamine Evaluated for Monitoring Diabetes. **Clin. Chem.**, v. 31, nº. 5, p.731-733, 1985.
- BURTIS, Carl A.; ASHWOOD, Edward R.; BRUNS, David E. **Tietz:** Fundamentos de Quimica Clínica. 6. ed. Rio de Janeiro: Saunders Elsevier, 2008. 959 p.
- YOUNG, D.S. **Effects of drugs on clinical laboratory tests - vol. 2**, 5 ed. Washington DC: AACC Press, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. **Clin. Chem.** v.27 p.493-501, 1981.