Beta Hidroxibutirato

Beta Hydroxybutyrate / Beta Hidroxibutirato
BT1004500

Responsável Técnico:
Dr. Gilson Serio Pizzo
CRF MG – 5310
Anvisa 80027310308

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ENGLISH

BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE

Kit intended to determine Beta Hydroxybutyrate in serum samples. Diagnostic use only.


STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagents ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE

Method: Enzymatic

Beta hydroxybutyrate (β-HB) from sample is oxidized to acetoacetic acid (ACAC) in the presence of NAD⁺, which is reduced to NADH by the catalytic action of the enzyme hydroxybutyrate dehydrogenase (β-HBDH). The concentration of beta hydroxybutyrate in the sample is proportional to the formation rate of NADH, which is spectrophotometrically measured at 340 nm.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum.

Collection and Handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Stability Period
Serum	4 to 8 °C	14 days
	-20 °C	40 days

PRODUCT DESCRIPTION

R 1 Tris 100 mmol/L and beta hydroxybutyrate dehydrogenase 2 KU/L.

R 2 NAD⁺ 2,5 mmol/L and calcium oxalate 20 mmol/L.

CAL Beta hydroxybutyrate in concentration determined in the lable, sodium chloride 100 mmol/L, BSA 20 g/L, sucrose 10% and preservative 0,03%. Traceable to the purified material 3-hydroxybutyric acid – Sigma Aldrich.

QUALITY CONTROL

The use of controls should be routine laboratory practice. It is recommended to perform the internal laboratory quality control to ensure the obtained results. For Calibration, it is recommended to use the calibrator below:

Calibrator – Beta Hydroxybutyrate CAL **REF** BT1004500

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 340 nm.
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

This product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotechnica.ind.br

Reaction Type	Kinetic / Positive
Wavelength	340 nm
Temperature	37 °C
Sample volume*	6 µL
R1 Volume*	180 µL
Sample + R1 Incubation Time	300 seconds
R2 Volume*	60 µL
R1 + Sample + R2 Incubation Time	90 seconds
Measure the absorbance continuously for 120 seconds	
Calibration	2 points

* Volumes can be modified as long as the established proportion is maintained.

B) CALCULATIONS

The equipment automatically calculates the analyte's concentration in each sample.

C) INTERPRETATION

Beta Hydroxybutyrate, also known as 3-hydroxybutyrate, is the main component of ketone bodies in the blood (representing 78%). It is synthesized in the liver from the oxidation of fatty acids. The formation of ketone bodies (ketogenesis) is elevated when blood glucose levels fall, during fasting and in the practice of prolonged exercise. Beta hydroxybutyrate is used in the early diagnosis and monitoring of the treatment for alcoholic or diabetic ketoacidosis, when its levels are abnormally high. In ketoacidosis, patients have an antagonistic effect on insulin, resulting in disturbance of protein synthesis and degradation metabolism. Therefore, the level of beta hydroxybutyrate during fasting in the morning can be measured to understand the mobilization of lipoproteins in patients with severe and acute injuries, severe illness, sepsis and nutritional support. Thus, beta hydroxybutyrate can be used as a marker of hepatic energy metabolism.

PERFORMANCE CHARACTERISTICS

Operating range
0.05 to 4.50 mmol/L

For concentrations above the operating range, dilute the sample with NaCl 150 mmol/L (0,9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity
Detection Limit
0.05 mmol/L

Analytical Specificity			
Hemoglobin	Bilirubin	Triglycerides	Ascorbic Acid
250 mg/dL	25 mg/dL	600 mg/dL	25 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	40 samples
Regression Equation	y = 0.9481x - 0.0085
Correlation Coefficient (R)	0.9980

Precision:

Within-run precision was determined with 3 lots and 10 replicates in an analytical run; between-run precision was determined with 3 lots and 3 replicates in an analytical run.

Samples (mmol/L)	Within-Run Precision	Between-Run Precision
	SD (mmol/L)	CV (%)
	SD (mmol/L)	CV (%)
0.90	0.011	1.23
	0.017	2.59

CV: Coefficient of variation; SD: Standard deviation.

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

Serum	0.02 – 0.27 mmol/L
-------	--------------------

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotechnica.ind.br or calling for +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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