**Proteína Transportadora do Retinol - RBP**

 Retinol Binding Protein / Proteína de Unión al Retinol  
 Ref. 20.029.00

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Kit intended to determine retinol binding protein in serum and plasma samples. Diagnostic use only.

**STORAGE AND HANDLING**

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagents ready for use.
- Once opened, the product in use is stable for 4 weeks when stored *on board*.
- Do not use reagents whose shelf life has expired.

**WORKING PRINCIPLE**
**Method:** Immunoturbidimetric  
 Retinol binding protein from sample binds to the specific anti-retinol binding protein antibody adsorbed in latex particles by an agglutination reaction. The formed antigen-antibody complex increases the turbidity of the solution, which can be spectrophotometrically measured at 570 nm. The intensity of the turbidity is proportional to the retinol binding protein concentration in the sample.

**SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY**
**Sample Type:** serum and plasma (EDTA or heparin).

**Collection and Handling:** collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

**Preservation:**

	Temperature	Stability Period
Serum and Plasma	4 to 8 °C	3 days
	-20 °C	6 months

**PRODUCT DESCRIPTION**

<b>R 1</b>	Buffer 50 mmol/L.	
<b>R 2</b>	Latex particles coated with rabbit anti-retinol binding protein antibody 5 g/L.	
<b>CAL S</b>	Retinol binding protein, buffer 50 mmol/L.	
<b>CONTROL S</b>	Retinol binding protein, buffer 50 mmol/L.	

**QUALITY CONTROL**

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and controls below:

Calibrador – Retinol Binding Protein CAL S	<b>REF</b>	20.029.00
Control – Retinol Binding Protein CONTROL S		20.029.00

**NECESSARY EQUIPMENT FOR TESTING**

- Spectrophotometer or photometer for reading at 570 nm.

**TEST PROCEDURE, CALCULATION AND INTERPRETATION**
**A) TEST PROCEDURE**

 This product is compatible to most types of automatic analyzers. Instrument settings are available at [www.biotecnicaltda.ind.br](http://www.biotecnicaltda.ind.br)

Reaction Type	Fixed Time / Positive
Wavelength	570 nm
Temperature	37 °C
Sample volume*	2 µL
R1 Volume*	180 µL
Sample + R1 Incubation Time	300 seconds
R2 Volume*	60 µL
R1 + Sample + R2 Incubation Time	30 seconds
First Reading	Immediately
Second Reading	270 seconds
Calibration**	Calibration curve

\* Volumes can be modified as long as the established proportion is maintained.

**\*\*CALIBRATION CURVE:**

Dilute the Retinol Binding Protein CAL S in a 9,0 g/L NaCl solution as follows:

Dilution	1	2	3	4	5
Calibrator (µL)	0	25	50	100	200
NaCl 9,0 g/L (µL)	200	175	150	100	0
Factor	0	0,125	0,25	0,5	1

Multiply the Retinol Binding Protein CAL S concentration indicated on the label by the corresponding dilution factor to obtain the retinol binding protein concentration of each dilution.

**B) CALCULATIONS**

The equipment automatically calculates the analyte's concentration in each sample.

**C) INTERPRETATION**

Retinol binding protein (RBP) is synthesized by the liver and its main function is to transport retinol (vitamin A) from the liver to the target-tissues (retina, skin, etc.). Later on, RBP is metabolized, filtered and reabsorbed in the kidney. RBP is a sensible marker for malnutrition, so that monitoring its serum concentration reveals the patient's nutritional state. Reduced serum and plasmatic concentrations of RBP can be observed in hypovitaminosis A, tubular nephropathies, hepatocellular failure and severe or acute inflammation. Increased concentrations of RBP are correlated with glomerular renal failure, type 2 diabetes or steatosis.

**PERFORMANCE CHARACTERISTICS**

Operating range
8.64 to 141.56 mg/L

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0,9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
1.75 mg/L	2.15 mg/L

Analytical Specificity		
Bilirubin	Hemoglobin	Triglycerides
20 mg/dL	500 mg/dL	500 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	120 samples
Regression Equation	y = 0.9961 + 1.6839
Correlation Coefficient (R)	0.9985

**Precision:**

Within-run precision was determined with 20 replicates in an analytical run; between-run precision was determined with 2 lots, in two periods of time, for 20 days.

Samples (mg/L)	Within-Run Precision		Between-Run Precision	
	SD (mg/L)	CV (%)	SD (mg/L)	CV (%)
36.31	0.5	1.3	1.20	3.26

CV: Coefficient of variation; SD: Standard deviation.

**RESIDUAL RISKS, WARNINGS AND PRECAUTIONS**

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.

**REFERENCE RANGES**

	Men	Women
Serum and Plasm	36 - 72 mg/L	22 - 53 mg/L

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

**WARNINGS AND PRECAUTIONS**

- Discard the reactions surplus according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or calling for +55 35 3214 4646

**QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE**

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific

Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotecnicaltda.com.br.

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