IVD

Homocisteína - HCY

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ENGLISH

BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE

Kit intended to determine homocysteine (HCY) in serum and plasma samples. Diagnostic use only.

STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagents ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8°C) are followed.
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE

Method:

Enzymatic

Oxidized homocysteine is reduced to free homocysteine, which reacts with serine, under catalytic action of cystathionine β-synthase (CBS), to form L-cystathionine. This, in the presence of cystathionine β-lyase (CBL) is converted back to homocysteine, releasing pyruvate and ammonia in the process. Pyruvate is converted to lactate by lactate dehydrogenase (LDH) and NADH is oxidized to NAD⁺. The concentration of homocysteine in the sample is determined by the consumption rate of NADH, which is spectrophotometrically measured at 340 nm.

Serine + Homocysteine

^{CBS}

L-cystathionine

^{CBL}

Homocysteine + Pyruvate + NH₃

^{LDH}

Pyruvate + NADH

Lactate + NAD⁺

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type:

serum and plasma (EDTA or heparin).

Collection and Handling:

collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Stability Period
Serum and Plasma	4 to 8 °C	4 weeks
	-20 °C	4 years

PRODUCT DESCRIPTION

R 1

Lactate dehydrogenase 35 kU/L; serine 0,76 mmol/L and NADH 0,47 mmol/L.

R 2

Cystathionine β-synthase 20 kU/L and cystathionine β-lyase 10 kU/L.

CAL

Homocysteine; sodium chloride 100 mmol/L; BSA 20 g/L; sucrose 10% and preservative 0,03%. Traceable to reference material NIST SRM1850.

CONTROL 1

Homocysteine $\approx$ 15 μ mol/L; sodium chloride 100 mmol/L, BSA 20 g/L; sucrose 10% and preservative 0,03%.

CONTROL 2

Homocysteine $\approx$ 37 μ mol/L; sodium chloride 100 mmol/L, BSA 20 g/L; sucrose 10% and preservative 0,03%.

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratorial quality control, it is recommended the use of the control below:

Calibrador – HCY CAL

BT1004200

Level 1 Control – HCY CONTROL 1

BT1004200

Level 2 Control – HCY CONTROL 2

BT1004200

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 340 nm.
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

This product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotechnica.ind.br

Reaction Type	Kinetic / Negative
Wavelength	340 nm
Temperature	37 °C

Sample volume*	16 μ L
R1 Volume*	234 μ L
Sample + R1 Incubation Time	180 seconds
R2 Volume*	36 μ L
R1 + Sample + R2 Incubation Time	210 seconds
Measure the absorbance continuously for 180 seconds	
Calibration	2 points

* Volumes can be modified as long as the established proportion is maintained.

B) CALCULATIONS

The equipment automatically calculates the analyte's concentration in each sample.

C) INTERPRETATION

Homocysteine is a sulfur-containing amino acid that is widely used as a cardiovascular disease marker. High levels of homocysteine are observed in coronary atherosclerosis and in myocardial infarction; the increase in the concentration of homocysteine is proportional to the risk of cardiovascular disease. The increase in the concentration of homocysteine is also related to the mortality risk in diabetic patients, kidney failure, Alzheimer's disease, non-alcoholic fatty liver disease, folic acid and vitamin B12 deficiency.

PERFORMANCE CHARACTERISTICS

Operating range

1.0 to 50.0 μ mol/L

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity

Detection Limit

1.0 μ mol/L

Analytical Specificity

Bilirubin	Hemoglobin	Triglycerides
25 mg/dL	250 mg/dL	500 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	40 samples
Regression Equation	$y = 0.999x - 0.767$
Correlation Coefficient (R)	0.9999

The product's accuracy was evaluated with a homocysteine control, at the decision level of 3,98 μ mol/L, with a total systematic error of 2,18%.

Precision:
Within-run precision was determined with 3 lots and 10 replicates in an analytical run; between-run precision was determined with 3 lots and 3 replicates in an analytical run.

Samples (μ mol/L)	Within-Run Precision		Between-Run Precision	
	SD (μ mol/L)	CV (%)	SD (μ mol/L)	CV (%)
28,0	0,95	3,45	1,17	2,96

CV: Coefficient of variation; SD: Standard deviation.

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

Serum and Plasma

4.0 - 15.4 μ mol/L

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotechnica.ind.br or calling for +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific

BIOTÉCNICA IND.COM. LTDA - CNPJ: 025340690001/20, Avenida Washington Ribeiro, 200 - Industrial Miguel de Luca
CEP 37072-030 Varginha-MG BRASIL Tel/fax: +55 35 3214 4646 www.biotechnica.ind.br

Revisão/Revision/Revisión: 01 – 10/04/2024

Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotecnica.ind.br.

- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotecnica.ind.br.

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