	20-50 ng/mL	Generally considered adequate for bone and overall health in healthy individuals.
	> 50 ng/mL	Linked to potential adverse effects, particularly at > 60 ng/mL.

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotechnica.ind.br or calling for +55 35 3214 4646.

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotecnicatda.com.br
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br

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Vitamin D is an essential steroid hormone for human health that plays a crucial role in regulating calcium and phosphorus metabolism, as well as maintaining bone mineral homeostasis. It exists in two forms: vitamin D2, primarily acquired through diet, and vitamin D3, which is synthesized in the skin through exposure to ultraviolet light. In the liver, vitamin D is converted into 25-hydroxyvitamin D (25-OH vitamin D), the main circulating form, which serves as an accurate indicator of an individual's vitamin D nutritional status. Adequate levels of 25-OH vitamin D are essential for the prevention of diseases such as osteoporosis, rickets, and osteomalacia. Additionally, its deficiency has been associated with a variety of clinical conditions, including cardiovascular and autoimmune diseases. Therefore, regular assessment of 25-OH vitamin D levels is fundamental not only for bone health but also for the prevention and management of other systemic diseases, making it a vital diagnostic tool in clinical practice.

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