Sensitivity		
Detection Limit	Quantification Limit	
10,1 ng/mL	17,97 ng/mL	

Analytical Specificity		
Bilirubin	Hemolysis	Lipemia
40 mg/dL	1000 mg/dL	1000 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in results. For drugs, consult the recommended reference (Young, 2000)

Accuracy	
Number of Samples	41 in duplicate
Regression Equation	$y = 0,9028x + 6,7362$
Correlation Coefficient (R)	0,9913

Serum: by applying the regression equation, the total systematic error estimated for decision levels of 60.0 ng/mL, 230 ng/mL and 500 ng/mL was 1,5%, -6,8% and -8,4%, respectively.

Precision					
Samples (ng/mL)	Within-Run Precision		Samples	Between-Run Precision	
	SD (ng/mL)	CV (%)		SD (ng/mL)	CV (%)
122	3,29	2,71	122	3,49	2,86
303	4,51	1,49	364	10,66	2,93
487	12,02	2,47	-	-	-

CV: Coefficient of variation; SD: Standard deviation.

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices’ instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination.
- The reagent contains preservatives. If it comes into contact with your eyes, mouth, or skin, rinse immediately with clean water and, if necessary, go to the hospital.

REFERENCE RANGES

Serum	60 – 230 ng/mL
-------	----------------

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotecnica.ind.br or calling for +55 35 3214 4646.

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and all current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotecnicaltda.com.br
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotecnica.ind.br

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- TURPEINEN, U.; HÄMÄLÄINEN, E. Determination of cortisol in serum, saliva and urine. Best Practice & Research Clinical Endocrinology & Metabolism, v. 27, n. 6, p. 795-801, 2013.
- GATTI, R.; ANTONELLI, G.; PREARO, M.; et al. Cortisol assays and diagnostic laboratory procedures in human biological fluids. Clinical Biochemistry, v. 42, n. 12, p. 1205-1217, 2009.
- LEVINE, A.; ZAGOORY-SHARON, O.; FELDMAN, R.; et al. Measuring cortisol in human psychobiological studies. Physiology & Behavior, v. 90, n. 1, p. 43-53, 2007.
- NIEMAN, L. K.; BILLER, B. M. K.; FINDLING, J. W.; et al. The Diagnosis of Cushing's Syndrome: An Endocrine Society Clinical Practice Guideline. Journal of Clinical Endocrinology & Metabolism, v. 93, n. 5, p. 1526-1540, 2008.
- MOORE, A.; AITKEN, R.; BURKE, C. W.; et al. Cortisol Assays: Guidelines for the Provision of a Clinical Biochemistry Service. Annals of Clinical Biochemistry, v. 22, n. 5, p. 435-454, 1985.
- WORLD HEALTH ORGANIZATION et al. Use of anticoagulants in diagnostic laboratory investigations. World Health Organization, 2002.

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