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ENGLISH

A rapid test for the diagnosis of Syphilis to detect antibodies (IgG and IgM) to *Treponema pallidum* (TP) qualitatively in human whole blood, Serum or plasma.

For professional *in vitro* diagnostic use Only.

INTENDED USE

The Syphilis Rapid Test Cassette (Whole Blood/Serum/Plasma) is a rapid chromatographic immunoassay for the qualitative detection of antibodies (IgG and IgM) to *Treponema pallidum* (TP) in whole blood, serum or plasma to aid in the diagnosis of Syphilis.

SUMMARY

Treponema pallidum (TP) is the causative agent of the venereal disease Syphilis. TP is a spirochete bacterium with an outer envelope and a cytoplasmic membrane. Relatively little is known about the organism compared to other bacterial pathogens. According to the Center for Disease Control (CDC), the number of cases of Syphilis

infection has increased markedly since 1985. ² Some main factors that contributed to this increase include the cocaine and crack epidemic and the high incidence of prostitution among drug users. drugs. ³ One study reported a significant epidemiological correlation between contamination and transmission of the HIV virus and Syphilis. ⁴ Syphilis features multiple clinical stages, long latency periods and asymptomatic infection. Primary syphilis is defined by the presence of a chancre at the site of inoculation. The antibody response to TP bacteria can be detected 4 to 7 days after the cancer appears. The infection remains detectable until the patient receives appropriate treatment. ⁵

The Syphilis Rapid Test Cassette (whole blood/serum/plasma) uses a dual antigen combination of a particle coated with Syphilis antigen and Syphilis antigen immobilized on the membrane to detect TP antibodies (IgG and IgM) qualitatively and selectively in the whole blood, serum or plasma.

PRINCIPLE

The Syphilis Rapid Test Cassette (whole blood/serum/plasma) is a qualitative membrane-based immunoassay for the detection of TP antibodies (IgG and IgM) in whole blood, serum or plasma. In this testing procedure, recombinant Syphilis antigen is immobilized on the test region. After the sample is added to the sample well of the cassette, it reacts with the Syphilis antigen-coated particles in the test. This mixture migrates through chromatography throughout the test and interacts with the immobilized Syphilis antigen. The dual antigen test format can detect both IgG and IgM in samples. If the sample contains TP antibodies, a colored line will appear in this region, indicating a positive result. If the sample does not contain TP antibodies, the colored line will not appear in the test region, indicating a negative result. To serve as procedural control, a line color will always appear in the control region, indicating that the appropriate volume of sample has been added and absorption by the membrane has occurred.

REAGENTS

The test contains Syphilis antigen-coated particles and Syphilis antigen coated on the membrane.

PRECAUTIONS

Read all information in this leaflet before running the test.

- In vitro* diagnostic use only. Do not use after the expiration date.
- Do not eat, drink or smoke in the area where samples or kits are handled.
- Do not use if packaging is damaged.
- Handle all samples as if they contain infectious agents. Observe established precautions against microbiological risks in all procedures and follow standard procedures for proper disposal of samples.
- Wear protective clothing such as laboratory coats, disposable gloves, and eye protection when samples are analyzed.
- The test used must be disposed of in accordance with local regulations.
- After use, the test must be disposed of in accordance with local regulations.
- Humidity and temperature can negatively affect results.
- This product is guaranteed to perform within specifications until the expiration date indicated on the labels, provided that the precautions for use and storage indicated on the labels and in these instructions are correctly followed.
- Single use product, do not reuse or recycle.
- The cassette, capillary tube/dropper, lancet and alcohol pad are disposable and cannot be reused

STORAGE AND STABILITY

Store as packaged in unsealed containers at room temperature or refrigerated (2-30°C) until the expiration date indicated on the packaging. The test is stable until the expiration date printed on the sealed packaging. The test must remain in the sealed packaging until use. **DO NOT FREEZE.** Do not use after the expiration date.

SAMPLE COLLECTION AND PREPARATION

- The Syphilis Rapid Test Cassette (Whole Blood/Serum/Plasma) can be performed using whole blood (from venipuncture or fingerstick), serum or plasma.
- Collect the sample according to standard procedures.
- If serum or plasma is used, separate the serum or plasma from the blood as soon as possible to avoid hemolysis. Use only clear, non-hemolyzed samples.
- The test must be performed immediately after the sample has been collected. Do not leave samples at room temperature for prolonged periods. Serum and plasma samples can be stored at 2-8°C for up to 3 days. For long-term storage, samples should be kept below -20°C. Whole blood collected by venipuncture should be stored at 2-8°C if the test is to be performed within 2 days of collection. Do not freeze whole blood samples. Whole blood collected by finger stick should be tested immediately.
- Allow samples to reach room temperature before testing. Frozen samples must be completely thawed and mixed well before testing. Samples should not be repeatedly frozen and thawed.
- If samples are transported, they must be packaged in accordance with local regulations for the transport of etiological agents.
- EDTA K2, sodium heparin, sodium citrate and potassium oxalate can be used as anticoagulants for sample collection.

MATERIALS

Materials provided

- Test Cassettes (with silica sachet)
- Dropper
- Buffer
- Instructions for use

Materials can be provided to replace or add as per request

- Capillary tubes
- Sterile lancets

Materials needed but not provided

- Sample collection containers
- Centrifuge
- Alcohol pad
- Timer

GUIDELINES FOR USE


(IMPORTANT NOTE): The instructions for use must be read before executing (performing) the test and the test expiration date must be checked. If the expiration date has passed, use another batch.

Organizing your workspace to perform the test.

- Allow the test, sample, and buffer to reach room temperature (15-30°C) before performing the test.
- Set up your work area on a clean, level surface, having the Provided Materials and the Materials Required but Not Provided to perform the test on hand.
- Open the sealed aluminum wrapper when it is at room temperature, remove the test cassette, and place on your work area. Use it as quickly as possible as the best results will be obtained if the test is carried out in approximately one hour. Identify the cassette with the patient's initials to avoid errors when delivering results;


For Serum or Plasma sample:

- Hold the dropper vertically and transfer **1 drop of serum or plasma (approximately 40 µL)** to the sample well, then add **1 drop of buffer (approximately 40 µL)** and start the timer, see the illustration opposite.

For venipuncture whole blood sample:

- Hold the dropper vertically and transfer **2 drops of whole blood (approximately 80 µL)** to the sample well, then add **1 drop of buffer (approximately 40 µL)** and start the timer. See the illustration below.


For the fingerstick whole blood sample:

1. Clean the puncture site (fingertip) with a cotton swab or cotton pad dipped in alcohol or wash the patient's hand with soap and water. Let it dry;
2. Press the finger to be punctured (ring or middle finger), without touching the area to be punctured, to promote greater blood accumulation in the area;
3. Remove the cap and press the lancet **firmly** onto the puncture site. Wipe away the first sign of blood with a sterile gauze or cotton pad;
4. Continue pressing your finger, with the help of your thumb and index finger, to apply light pressure near the puncture site and increase blood flow, in order to form a large, rounded drop of blood at the puncture site, facilitating collection;

Using a dropper :


- Squeeze the dropper bulb first, then make the open end of the dropper come into contact with the blood, and then slowly release it to let the blood enter the dropper until the blood rises near of the line marked on the dropper.
NOTE: Do the entire process slowly and do not release the bulb too quickly. When the sample reaches the mark, do not relax the bulb, hold it.

- Squeeze the bulb again to release all of the blood (approximately 2 drops) into the sample well (S). If there are many bubbles or the sample volume is not enough, continue pressing your finger several times and with intensity, repeat step 1 and add 1 drop of blood.
- Add **1 drop of buffer** and start the timer.


Using a capillary tube :

- Hold the capillary tube, without pressing the bulb, in a horizontal position and parallel to the finger or tilted at a maximum angle of 10° measured from the horizontal.
- Touch the open end of the capillary tube to the drop of blood. The blood will enter the capillary tube by capillary action.
- Continue pressing your finger several times and with intensity until blood fills the capillary tube. Repeat this procedure if the amount of blood is insufficient and try to avoid bubbles and excessive movement of the capillary tube during collection .
- Squeeze the bulb of the capillary tube to dispense the entire volume of whole blood into the sample well of the test cassette, then vertically add **1 drop of buffer solution** into the sample well, start the timer. It is expected that blood residues will remain on the inner walls of the capillary tube;


Using hanging blood drops:

- Allow 2 suspended drops of fingerstick blood sample to fall into the sample well, then add 1 drop of buffer and begin timing.

- Wait for the lines to form before reading the results.
- **Read the results in 5 minutes** . Do not interpret after 20 minutes.

INTERPRETATION OF RESULTS

(See illustration above)

POSITIVE: **Two lines are displayed.** One colored line must be in the control region (C) and another apparent colored line must be in the test region (T).
***NOTE:** The intensity of the color in the test region (T) will vary depending on the concentration of *TP antibodies* present in the sample. Therefore, any color shade in the test region (T) should be considered positive.

NEGATIVE: **One colored line appears in the control region (C).** No lines appear in the test region (T).

INVALID : **A control line does not appear.** Insufficient sample volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, immediately discontinue use of the test kit and contact your local distributor.

DISPOSAL OF INPUTS

- Kit materials that have had contact with the patient's biological material must have a specific destination:
- Lancets must be disposed of in sharps waste.
- Cassettes (test devices), capillary tubes (optional), droppers (optional), sterile lancets (optional), alcohol pads (not provided) and buffer should be disposed of in biological waste.
- Boxes, instructions for use, aluminum packaging and silica sachet can be disposed of in general trash.
- Instructions must be followed exactly to obtain accurate results.
- To dispose of the materials used, follow the guidelines of local biosafety standards.

QUALITY CONTROL

An internal procedural control is included in the test. One colored line that appears in the control region (C) is considered an internal procedural control. It confirms sufficient sample volume, adequate membrane absorption, and correct procedural technique. Control standards are not provided with this kit; however, it is recommended that positive and negative controls be tested as good laboratory practice to confirm the testing procedure and verify proper test performance.

DIAGNOSTIC MEANING

The diagnostic interval is the period between detection of infection and detection of antibodies to the primary infection. During this period, tests that detect antibodies against TP are unable to obtain a reactive result. Specific Anti-TP antibodies are produced approximately 2 weeks after exposure. Anti-TP IgM and IgG can be detected within 3 days after the appearance of the lesion in primary syphilis.

LIMITATIONS

1. The Syphilis Rapid Test Cassette (Whole Blood/Serum/Plasma) is for *in vitro* diagnostic use only. The test should only be used for the detection of *TP antibodies* in whole blood, serum or plasma. Neither the quantitative value nor the rate of increase of *TP antibodies* can be determined by this qualitative test.
2. The Syphilis Rapid Test Cassette (Whole Blood/Serum/Plasma) will only indicate the presence of *TP antibodies* in the sample and should not be used as the sole criterion for diagnosing *TP infection* .
3. As with all diagnostic tests, all results must be interpreted in conjunction with other clinical information available to the physician.
4. If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not preclude, at any time, the possibility of a *TP infection* .

PERFORMANCE CHARACTERISTICS

Sensitivity and Specificity

The Syphilis Rapid Test Cassette (Whole Blood/Serum/Plasma) identified samples correctly from a performance panel and was compared to a leading commercial TPPA Syphilis test using clinical samples. The test results were tabulated as below.

Method	TPPA		Total Result	
	Results	Positive		Negative
Syphilis Rapid Test Cassette (Whole Blood/Serum/Plasma)	Positive	299	1	300
	Negative	1	419	420
Total Result		300	420	720

Relative Sensitivity: 99.7% (95%CI*: 98.2% -> 99.9%)

Relative Specificity: 99.8% (95%CI*: 98.7% -> 99.9%)

Accuracy: 99.7% (95%CI*: 98.9% -> 99.9%)

*Confidence interval

Precision Intra-Assay

Within-analysis precision was determined using 10 replicates of four samples: one negative, one low positive, one medium positive, and one high positive. Negative, low positive, medium positive, and high positive values were correctly identified >99% of the time.

Inter-Rehearsal

Between-assay precision was determined by 10 independent assays on the same four samples: one negative, one low positive, one medium positive, and one high positive. Three different lots of the Syphilis Rapid Test Cassette (Whole Blood/Serum/Plasma) were tested over a 3-day period using negative, low positive, medium positive, and high positive samples. Samples were identified correctly >99% of the time.

Cross Reactivity

The Syphilis Rapid Test Cassette (Whole Blood/Serum/Plasma) has been tested for positive samples of HAMA, RF, HBsAg, HBsAb, HBeAg, HBeAb, HBcAb, HCV, HIV, *H.Pylori*, MONO, CMV, Rubella and Toxoplasmosis. The results showed no cross-reactivity.

Interfering Substances

The following substances with potential for interference were added to the Syphilis negative and positive samples.

Paracetamol: 20 mg/dL	Caffeine: 20 mg/dL
Acetylsalicylic Acid: 20 mg/dL	Gentisic Acid: 20 mg/dL
Ascorbic Acid: 2 g/dL	Albumin: 2 g/dL
Creatine: 200 mg/dL	Hemoglobin: 1.1 mg/dL
Bilirubin: 10 mg/dL	Oxalic Acid: 600 mg/dL

None of the substances interfered with the test at the concentrations tested.

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1. Claire M. Fraser. *Complete genome sequence of Treponema Pallidum, the Syphilis spirochete*, Science 1998; 281 July: 375-381
2. Center for Disease Control. *Recommendations for diagnosing and treating Syphilis in HIV-infected patients*, MMWR Morb. Mortal Wkly Rep. 1988; 37: 601
3. Aral R. Marx. *Crack, sex and STD*, Sexually Transmitted Diseases, 1991; 18:92-101
4. J.N. Wasserheit. *Epidemiological Synergy: Interrelationships between human immunodeficiency virus infection and other sexually transmitted diseases*, Sexually Transmitted Diseases 1992; 19:61-77
5. Johnson Phillip C. *Testing for Syphilis*, Dermatologic Clinic 1994; 12 Jan: 9-17

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Fabricante de lanceta estéril:

PROMISEMED HANGZHOU MEDITECH CO., LTD.

No. 1388 Cangxing Street, Cangqian Community, Yuhang District, Hangzhou City, 311121 Zhejiang, China
Número Ref.: 14603054900

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Hangzhou AllTest Biotech Co., Ltd.